

Hero-NXT Car Seat

The Hero-NXT child car seat is a specifically modified car seat to accommodate the postural requirements of young children with additional needs. It is suitable for installation in most commercially available vehicles that have an ISOFIX connection. The car seat is designed for use on the back seat of the car. Due to its swivel base, the vehicle needs to have ISOFIX for installation.

Swivelling car seat for forward or rear-facing



ACCESSORIES

Security Tray
4181-5530-000



FUNCTIONAL BENEFITS

- The Hero-NXT can be used both forward and rear-facing direction.
- To support children with limited upper body control, the Hero-NXT has a tilt function to adjust the seat angle.
- The Hero-NXT is equipped with its own 5-point chest and shoulder strap that securely holds your child in place.

Package includes

- Black/Red upholstery
- Rear/Forward facing position
- Height adjustable back support
- 10° tilt with swivel base 360°
- Five-point harness
- Lateral supports as standard



Practical and functional



In direction of travel;
adjustable back height



Rotatable by 90°, 180° or
360° for ease of transfers



In rearward position



Lateral Supports
fitted as standard

Code	Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Crown Height	Distance Laterals	Total Height	Total Width	Total Depth	Capacity
4181-0010-020	28cm	21cm	25cm	21cm - 34cm	41cm - 53cm	12cm - 21cm	66cm - 73cm	43cm	63cm	18kg

*Measurements with standard cushion