


Hero-NXT



**SPECIAL
PURPOSE
CAR SEAT**

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements:	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Basic Packages

Includes:	Package Sizes
• Black/Red colour • 10° Tilt-in-space • 360° Swivel base • Rear/forward facing position • Five-point harness • Height-adjustable back support • Adjustable lateral side supports	One Size 4181-0010-020 

Specifications

Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Crown Height	Distance Laterals	Total Height	Total Width	Total Depth	Capacity
28cm	21cm	25cm	21cm - 34cm	41cm - 53cm	12cm - 21cm	66cm - 73cm	43cm	63cm	18kg

*measurements with standard cushion

Accessories

Lateral Supports (Recommended)		Abduction Block	
One Size, Pair 4181-3310-000		One Size 4181-1510-000	
5-Point Seat Belt Extension		Padded Tray Table	
One Size Available as special order only. 4081-6530-000		One Size 4181-5530-000	
Seat Protection Mat			
One Size 4181-0190-000			

Notes



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com

