

Ipai-LGT




CAR SEAT

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Footplate Height	
E Chest Width	
F Axilla Height	
User Weight	



Specifications

Seat Width (pelvis)	Seat Width (knee)	Seat Depth (min-max)	Back Height	Crown Height	Height Laterals	Distance Laterals	Weight	Capacity
310mm	340mm	320mm	380mm - 580mm	590mm - 790mm	200mm - 330mm	140mm - 280mm	7kg	36kg

Packages

Standard Package

Includes:

- Height Adjustable Headrest
- Lateral Supports
- 5-Point Harness
- ISOFIX Attachment
- Black Upholstery

One Size
4081-0021-060



Pen Key Release Upgrade

Upgrade Ipai-LGT with this reduced access buckle, supplied with emergency webbing slicer and release tool.
4081-6620-025



Accessories

Padded Tray Table

4080-5411-000



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY: