

Ipai-LGT Car Seat

The Ipai-LGT lives up to its name, characterised by a light weight easy-to-use design and providing mild body support and care for children with additional needs. Comes standard equipped with ISOFIX (for additional stability during transfers), thoracic supports and a 5-point harness that is padded at the shoulders. The large, stable tee pad within the lower portion of the harness assists in providing firm support across the child's pelvis.

Light weight, comfortable and safe!



Side Impact Protection System

Simple height adjustment providing comfort and safety

Lightweight option making transfer between vehicles easy

5 point comfort harness comes as standard

Package includes:

- Supportive 5 point harness with dual adjustment
- Backrest height adjustable
- Adjustable lateral supports

PEN KEY RELEASE UPGRADE

Upgrade IPAI LGT with this reduced access buckle, supplied with emergency webbing slicer and release tool.
4081-6620-025



FUNCTIONAL BENEFITS

- Simple height adjustment function allows the backrest 20cm of height range providing greater comfort and safety for the child.
- IPAI-LGT comes with ISOFIX connection points which can be used to provide additional stability to the seat during transfers. ISOFIX is not required for installation of the IPAI-LGT.
- Have the peace of mind that your child is kept safe with the Side Impact Protection System that can be unscrewed to secure the car seat against the side of the vehicle door, improving safety in the event of side impact to the vehicle.



Backrest Height



Adjustable depth ISOFIX, optional use for additional stability



Side Impact Protection System



Padded Tray Table
Code: 4080-5411-000

Code	Seat Width (rear)	Seat Width (front)	Seat Depth	Back Height	Head Support Height	Lateral Support Height	Distance between Laterals	Back Support Recline Angle	Total Height	Overall Width	Total Depth	Total Weight	Capacity
4081-0011-060	31cm	34cm	32cm	38-58cm	59-79cm	20-33cm	14-28cm	85°-162°	67-85cm	48cm	39cm	7kg	36kg