

Jenx Standz Size 1

Standz offers abducted hip positioning for children between the ages of 1-9. It can be used with or without abduction, depending on individual needs and goals.

Supine/Prone/Upright/Abduction Standing Frame

Built in angle gauges for frame angle and individual leg abduction



Tool free adjustments - cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort

Standz covers have been designed to wipe clean, they also contain an antimicrobial agent to aid infection control

Adjustable supports to accommodate knee flexion contractures and leg length discrepancies



Base Includes:

- Base with castors
- Hydraulic assisted tilt
- Pelvic and trunk support pads with strap
- Abducting legs
- Foot positioning plates
- Multiposition tray with extension wings
- Tray mount

The Standz is able to be used in the Prone or Supine configuration, depending on your child's needs



ACCESSORIES

Flat Pad Head Support

5240-4110-044



Multigrip Head Support

1 Link: 5356-4601-000
2 Link: 5356-4602-000
3 Link: 5356-4603-000



Flat Shoulder Support

5240-2750-044



Multigrip Body Support

5240-2721-044



Pelvic Support Board

5240-1901-044



Back Pad

5240-2402-044



Standard Knee Support

5240-1750-044



Multi-Axis Knee Support

5240-1770-044



Tray Holder

5240-7112-044



Knee Block

Size 1 5356-1751-000
Size 2 5242-1767-054




Sandals

Size 1 5356-6081-000
Size 2 5356-6082-000



Knee Block Mounting Bracket

5240-1720-000



Code	Angle Range	Age Range	Max User Height	Frame Width	Frame Length	Chest & Hip Pad Width
5240-0005-044	0-90°	1-9yrs	138cm	63.5cm	90cm	18-30cm
Knee Block Width	Height of Stander	Leg Abduction	Max User Weight	Foot Support Angle Dorsiflexion	Foot Support Angle Plantarflexion	
12cm	108cm	30°	45kg	10°	-25°	