

Jenx Standz Size 1

Standz easily accommodates a range of abducted hip positioning for children between the ages of 1-9.

Supine/Prone/Upright/Abduction Standing Frame

Built in angle gauges for frame angle and individual leg abduction



Tool free adjustments - cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort

Standz covers have been designed to wipe clean, they also contain an antimicrobial agent to aid infection control

Adjustable supports to accommodate knee flexion contractures and leg length discrepancies

Supine Package 1/2 Includes: Prone Package Includes:


- Gas Tilt
 - Angle Adjustable Footplates
 - Size 1/2 Sandals
 - Knee Supports with multi-axis technology
 - Incorporated Flexed Knee Pads
 - Pelvic Support Pads
 - Trunk Support Pads
 - Flat Shoulder Support (Package Opt 2)
 - 1/2 Link Multigrip Head Support
 - Wipe-clean Yellow Vinyl Covers
 - Tray with Extension Wings and Tray Holder
- Gas Tilt
 - Angle Adjustable Footplates
 - Size 1 Sandals
 - Standard Knee Supports
 - Pelvic Support Pads
 - Trunk Support Pads
 - Wipe-clean Yellow vinyl covers
 - Tray with Extension Wings

The Standz is able to be used in the Prone or Supine configuration, depending on your child's needs



ACCESSORIES

- Flat Pad Head Support**
5240-4110-044
- Multigrip Head Support**
1 Link 5356-4601-000
2 Link 5356-4602-000
3 Link 5356-4603-000



- Flat Shoulder Support**
5240-2750-044
- Multigrip Body Support**
5240-2721-044



- Pelvic Support Board**
5240-1901-044
- Back Pad**
5240-2402-044




- Standard Knee Support**
5240-1750-044
- Multi-Axis Knee Support**
5240-1770-044



- Tray Holder**
5240-7112-044
- Knee Block**
Size 1 5356-1751-000
Size 2 5242-1767-054



- Sandals**
Size 1 5356-6081-000
Size 2 5356-6082-000
- Knee Block Mounting Bracket**
5240-1720-000



Code	Angle Range	Age Range	Max User Height	Frame Footprint	Chest & Hip Pad Width	Knee Block Width	Height of Stander	Leg Abduction	Max User Weight	Foot Support Angle DorsiFlexion	Foot Support Angle PlantarFlexion
5240-0011-044 (prone) 5240-0020-044 (supine) 5240-0020-044 (supine opt 2)	0-90°	1-9 yrs	138cm	63.5 W x 90 Hcm	18-30cm	12cm	108cm	30°	45kg	10°	-25°