

Jenx Standz



ABDUCTION STANDER

Distributed by
medifab


NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

May 2025

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

User Measurements	
A Hip width	
B Axilla height	
C Chest width	
D Overall user height	
E Heel to knee	
F Hip height	
User weight	



Basic Packages

Standz Size 1, Yellow (Prone)	Standz Size 1, Yellow (Supine)
5240-0011-044	5240-0020-044
Includes: • Gas Tilt, • Angle Adjustable Footplates • Size 2 Sandals • Standard Knee Supports • Pelvic Support Pads • Trunk Support Pads • Wipe-clean Yellow vinyl covers • Tray and Tray extension wings	Includes: • Gas Tilt • Angle Adjustable Footplates • Size 1 Sandals • Knee Supports with Multi-axis technology and incorporated flexed knee pads • Pelvic Support Pads • Trunk Support Pads • 1 Link Multigrip Head Support • Wipe-clean Yellow vinyl covers • Tray with extension wings and a tray holder
Standz Size 1, Yellow (Opt 2, Supine)	
5240-0021-044	
Includes: • Gas Tilt • Angle Adjustable Footplates • Size 2 Sandals • Knee Supports with Multi-axis technology and incorporated flexed knee pads • Pelvic Support Pads • Trunk Support Pads • Flat shoulder support • 2 Link Multigrip Head Support • Wipe-clean Yellow vinyl covers • Tray with extension wings and a tray holder	

Product Details

Size 1 - Manual Tilt

Average Age	1 - 9 Years
Weight Capacity	45 kg
Hip Width	17.5cm - 30.5cm
Hip Height	45cm - 70cm
Chest Width	17.5cm - 30.5cm
Chest Height	58cm - 100cm
Max Height	138cm
Tilt Range (90° = Upright)	5° to 90°
Frame Width	64.5cm
Frame Length	89cm
Height of surface when horizontal	81cm
Abduction Range	30° Either Leg

Key

Blue	Suitable for Prone and Supine
Green	Suitable for Supine only
Yellow	Suitable for Prone only

Items outlined in blue are required and an option must be selected

Base unit



Standz Basic Package

Tool free adjustment standing frame includes base with castors, hydraulic assisted tilt, pelvic and trunk support pads with strap, abducting legs, foot positioning plates, multiposition tray with extension wings.

5240-0005-044




Knee support options

Standard Knee Support		Multi-Axis Knee Support	
Dual wall knee supports with accessories available for extra function (see below). A: 11.5cm, B: 12.5cm, C: 9.5cm & Min: 6.5cm 5240-1750-044		Width, height, depth and angle adjustable, with soft arched foam pads behind the knee and support straps to enable knee flexion. A: 13.5cm & B: 11.5cm 5240-1770-044	

Knee Blocks		Knee Block Mounting Bracket	
Flexible PU knee cups, complete with padded knee straps. Adjustable for height, width, and depth.			
Knee Blocks Size 1 Internal Support Size: 7.5cm W x 10cm H 5356-1751-011		Mounting bracket required when using knee block option for knee support. 5240-1720-000	
Knee Blocks Size 2 Internal Support Size: 9.5cm W x 12cm H 5356-1752-011			

Accessories for Standard Knee Support

Standard Strap	De-rotational Strap	Arched Foam Pads	Qty
(Required for Supine use). Holds knee in position and prevents knee flexion. Can also be used in Prone. 5240-1766-044	Adjustable strap for reducing internal or external rotation. 5240-1786-044	(Recommended for Supine use). Ensures a comfortable standing position when accommodating knee flexion. Supplied individually. 13.5cm W x 13.5 cm H x 5.5cm D 5353-1775-000	
			

Foot support options

Size 1 sandals		Size 2 sandals	
Ensure solid foot placement 8cm x 17cm long 5356-6081-000		Ensure solid foot placement 9.5cm x 22cm long 5356-6082-000	

Head support option

Flat Pad	
For children with reasonable head control. Can also be used as chest pad in prone. 19cm W x 10cm H 5240-4110-044	

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm W x 5cm D x 12cm H 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Shoulder support option

Flat Shoulder Support	Multigrip body support (for prone and supine)		
Shoulder support for supine use only (too wide for prone) 31.5cm W x 9.5cm H 5240-2750-044		Adjustable contour shoulder support 41cm max width 5240-2721-044	

Pelvic support options

Pelvic Support Board	Back Pad (Requires Pelvic support board)		
Lower back support 19.5cm W x 15cm H 5240-1901-044		Central back support 19.5cm W x 15cm H 5240-2402-044	

Accessories

Anti Tamper Conversion Kit

Allows for all wing knob adjustments to be converted to a tool fixing – ideal for environments where easy adjustment can be problematic.

5240-0170-000




Tray Elbow Blocks

Aids in supporting the arms and promoting a functional arm and hand position.

Only usable in conjunction with Tray Extension Wings.

10cm W x 14cm H

5356-7351-000



Tray holder (useful for Supine)

When in Supine, the tray needs to be removed to remove the child. Place the tray on the convenient tray holder to avoid tripping.

5240-7365-000



Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other markings on the paper.



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Standz_20250529_NZ