

Kidsflex Car Seat

Excellent posture support options and suitable for installation with most car types with ISOFIX. Available with a range of options and adjustments, allowing complete comfort and postural management for the child. Due to its swivel base, the vehicle needs to have ISOFIX for installation.

Fully adjustable specialist seating

Secure seat using the existing **3 point safety belt** fitted in the rear seat.

5 point harness comes standard with abdominal padding



Adjustable tilt angle and back support recline can be tailored to suit the users postural needs

Slimline seat cushions allow more space and headroom inside the car

Package includes:

- 5 Point harness
- Recline function
- Tilt angle adjustment (0 to 20°)
- ISOFIX attachment
- Standard head support
- Slim swivel base - left turn only
- Seat belt guide on the left




Size	Code	Max Weight	Seat Width	Seat Depth	Shoulder Height	Lateral Height	Lateral Width	Weight
2	4170-2022-060	36kg	31-35cm	34-39cm	48-55cm	25-29cm	23-32cm	8.9kg
3	4170-2023-090	75kg	36-40cm	40-46cm	55-64cm	30-36cm	29-40cm	12.3kg

ACCESSORIES

Support Tray Table

2: 4170-5532-000
3: 4170-5533-000



Foot Support

2/3: 4170-5242-000



Small Lateral Supports

Fixed: 4170-3301-000
Swing-away: 4170-3401-000



Large Lateral Supports

Fixed: 4170-3302-000
Swing-away: 4170-3402-000



FUNCTIONAL BENEFITS



Easily adjustable in seat depth.



Angle adjustable back support



Ample growth opportunities.