

Klaxon

Spare Parts & Accessories



Distributed by
medifab

KLAXON[®]

ORDERING GUIDE

October 2023

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
Chair Model	
Chair Width	
Client Weight	
Tube Diameter	

Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Connectors

Connector Front Side

	Small 1595-CF-10-39-0002	Standard 1595-CF-10-39-0000	Extra Long 1595-CF-10-39-0001
---	-----------------------------	--------------------------------	----------------------------------

Connector Rear Side

	Small 1595-CF-12-16-0001	Standard and Extra Long 1595-CF-12-16-0000	
---	-----------------------------	---	--

Complete Connector (Front Connector + Rear Connector + Front Hook)

	Small 1595-LSTD-0002	Standard 1595-LSTD-0000	Extra Long 1595-LSTD-0001
---	-------------------------	----------------------------	------------------------------

Complete Connector Tetra (Complete Connector + Quick-Release Lever Kit + Velcro Grip Kit)

	Standard 1595-LTET-0000	Extra Long 1595-LTET-0001
---	----------------------------	------------------------------

Clamping Rings (for Lateral Linking System only)

(requires 2 per chair)

		Left	Right
	25mm Round Shape	1595-AGLSX25	1595-AGLDX25
	32mm Round Shape	1595-AGLSX32	1595-AGLDX32
	25x35mm Oval Shape (Progeo)	1595-AGLSXOV	1595-AGLDXOV
	27x33mm Oval Shape (Sunrise)	1595-AGLSXOVS	1595-AGLDXOVS

Round-shaped, Convergent Tube Frames (1 piece - for chairs with tapered seat only)

	Diameter 25,4mm 1595-AGCP125	Diameter 32,4mm 1595-AGCP132
---	---------------------------------	---------------------------------

Round-Shaped, Parallel Tube Frames (1 piece - may require a set of RID spacers depending on tube size)

	Diameter 25,4mm 1595-AGP25	Diameter 32,4mm 1595-AGP32	Diameter 34,925mm 1595-AGP35
---	-------------------------------	-------------------------------	---------------------------------

Clamping Rings Continued

Clamping Ring for Oval-shaped, Parallel Tube Frames (1 piece)

	For Progeo 25mm x 35mm (WxH) 1595-AGPOV	For Sunrise 27mm x 33mm (WxH) 1595-AGPOVSX
---	---	--

Clamping Ring for Progeo Carbon Noir 2.0 and Duke


	(1 piece) 1595-AGPOVPC
---	---------------------------

Position-Adjustable Metal Lever for Clamping Ring	Pawl Kit for Clamping Ring (5 pieces)
---	---------------------------------------

(for Lateral Linking System only) 1595-CL-80-0000		(for Lateral Linking System only) 1595-KIT-AGL-0000	
--	---	--	---

Reducers


Reducer (To be used with Clamping Rings depending on frame tube size)

	Converts 25, 4 into 18/18, 4	1595-RID25-18-0000
	Converts 25, 4 into 20/20, 4	1595-RID25-20-0000
	Converts 25,4 into 23	1595-RID25-23-0000
	Converts 25,4 into 22/22,4	1595-RID25-22-0000
	Converts 32,4 into 28/28,4	1595-RID32-28-0000
	Converts 32,4 into 30/30,4	1595-RID32-30-0000

Hook System

Front Hook 2.0		Hook Lock-in System	
(for Lateral Linking System only)		(mounted at the forefront of Klick Connector) (for Lateral Linking System only)	
1595-CL-1C-0000		1595-CF-1C-0000	

Lateral Linking

Lateral Linking System			
Includes: Side Forks, Protective Sponge, Crossbeam, Brackets, Front Hook, Lower Pin (for Lateral Linking System only)			
1595-CL-KIT-0000			
Side Forks (Pair)		Protective Sponges (pair)	
(for Lateral Linking System only)		For side forks (for Lateral Linking System only)	
1595-CL-12-0001		1595-CL-12-SPONGE	
Lower Pin – Stainless Steel		Tetra Support Velcro Grip Kit with Logo	
(for Lateral Linking System only)		Replacement Strap for Tetra Connector (included in complete Tetra Connector)	
1595-CC-18-0010		1595-TEV02	

Crossbeam Options

Complete Crossbeam Range (Includes Crossbeam, Fixing Clamp and Aluminium Axle)

	230-270mm (Small Crossbeam) Seat width: 12-13"	1595-TRAV2327
	271-290mm (Medium Crossbeam) Seat width: 13-14"	1595-TRAV2729
	291-360mm (Medium Crossbeam) Seat width: 14-16"	1595-TRAV2936
	361-480mm (Medium Crossbeam) Seat width: 17-20"	1595-TRAV3648
	480-540mm (Medium Crossbeam) Seat width: 18-22"	1595-TRAV4854

Central Crossbeam

(for Lateral Linking System only)	
1595-CL-TRAV-0000	

Crossbeam

	Small 1595-CF-33-RID-0001	Medium 1595-CF-33-380-0001	Large 1595-CF-33-480-0001
---	------------------------------	-------------------------------	------------------------------

Fixing Clamp for Crossbeam (1 piece)

Brackets

(2 parts required per crossbeam) Qty:		For the connection of the Front hook 2.0 and the Central crossbeam (for Lateral Linking System only)	
1595-CF-32-0020		1595-CL-33-0000	

Aluminium Axle (1 piece)

(2 parts required per crossbeam)			
	Small 1595-CF-31-0001	Standard 1595-CF-31-0000	Extra Long 1595-CF-31-0002
Qty:	Qty:	Qty:	Qty:

Handlebar Options

Standard Handlebar		Handles			
	Standard handlebar for Klick Race, Power, Monster 1595-CC-25-0000	Standard grips for all models excluding Klick Hybrid Standard			
		1595-CC-26-0000			
Universal Clamping Ring for Tetra Handle		Tetra handlebar with Comfort Drive System			
(1 piece) 1595-T-HANDLE-00			(RT, P+T, R+T) 1595-TEMA2		
Tetra Handle (1 Piece)					
Horizontal		Vertical		Tripod	
1595-T-HANDLE-01		1595-T-HANDLE-02		1595-T-HANDLE-03	
Adjustable Handlebar Stem		Quick Release Adjustable Handlebar Stem			
1595-CC-23-22-24-0002		(ES, PS)			
		1595-CC-23-22-24-0001			
Tetra Support Quick Release Lever KIT					
(for Klick Race, Power, Monster)					
1595-TESG1					
Appendices for Tetra Handlebar		Handlebar Extension Standard			
1595-T-HORN-0000		1595-CC-25-0005			

Tyre Options

Spare Tyers

	20 x 210 (HS)	1595-400730
	FAT (FS)	1595-CM-3-2000
	Tread Tyre 14"	1595-CC-3-1400
	Tread Tyre 14" Smooth	1595-CC-3-1401

Fat Wheels (Pair) with Complete Mounting Kit

12,7mm 1595-FATWHEELS-01-12,7		12mm 1595-FATWHEELS-01-12	
----------------------------------	---	------------------------------	---

Spare Tyre Tubes

	14" 1595-CC-87-1400	20" x 2.00 (HS) 1595-400520	20" (FS) 1595-CM-87-2000
--	------------------------	--------------------------------	-----------------------------



Complete Splashguard Mounting Kit

Race & Power		Monster		Hybrid Standard	
14" Black 1595-CC-36-0000		(FS) 1595-FS-36-0000		(HS) 1595-490100	

Rim Brakes with Mounting Kit

1595-CC-8-9-10-0000	
---------------------	---

Disk Brakes with Mounting Kit

(MS, CS) 1595-MK-6-0000		(all models, excluding MS, CS) 1595-CC-6-7-0000	
----------------------------	---	--	---

Motor Options

Motor (set)

	250W brushless Motor (Klick Mini)	1595-MK-5-0000
	48V 250W (Brushless) Motor + 14" rim and spokes (P, P+L, PSL, PTL) (2021)	1595-CC-5-0003
	250W Motor (Brushless) + 14" rim and spokes (ES, ET, PS, PT) (2019)	1595-CC-4-5-1400
	500W Motor (Brushless) + rim and spokes (Klick Hybrid Standard)	1595-CH-4-5-2000
	1000W (Brushless) Motor (Klick Race)	1595-CC-5-0001
	1000W (Brushless) Motor + rim and spokes (Klick Monster)	1595-CM-4-5-2000

Battery Options

Spare Battery

	36V 5,8Ah (ES, ET) (2019)	1595-CC-1-5AH
	48V 11Ah (Klick Hybrid Standard)	1595-CH-1-11AH
	36V 5,8Ah (Klick Mini)	1595-MK-1-5AH
	Fly Battery 36V 2,5Ah (PS, PT) (2019)	1595-CC-1-2,5AH-PLUS
	Fly Battery 36V 2,5Ah (ES, ET) (2019)	1595-CC-1-2,5AH-SLIM
	Fly Battery 36V 2,5Ah (Mini Fly)	1595-MK-1-2,5AH
	Fly Race 48V 2,9Ah (P, R, M)	1595-CM-1-2,9AH
	Race 48V 11,6Ah (P, R, M)	1595-CM-1-11AH
	Race Lite 48V 5Ah (P, R, M)	1595-CM-1-5AH

Spare Battery Charger

	36V 5,8Ah (ES, ET, Fly version too)	1595-CC-89-0000
	48V 11,6Ah (RS, RT, FS, PSL, PTL, Fly version too)	1595-CC-88-0001
	36V 11,6Ah (PS, PT, Fly version too)	1595-CC-88-0000
	48V (P, R, M)	1595-CC-88-0002
	48V 11,6Ah (Klick Hybrid Standard)	1595-K500050
	36V 5,8Ah (Klick Mini)	1595-MK-88-0000

Control Options

Accelerator handle (ES, PS)

Left 1595-CC-28-0001		Right 1595-CC-28-0002	
-------------------------	---	--------------------------	---

Accelerator handle (RS, FS, PSL, P, R, M)

Left 1595-CC-28-0004		Right 1595-CC-28-0005	
-------------------------	---	--------------------------	---


Accelerator Lever

	(ES, PS, MS, CS) 1595-CC-28-0000	(RS, FS, PSL, P, R, M) 1595-CC-28-0003
--	-------------------------------------	---

Brake Levers with Electronic Sensor LS+RS

(ES, PS, MS, CS, HS) 1595-CC-30-0000		With parking brake LS+RS (RS, FS, PSL, P, R, M) 1595-CC-30-0002	
---	---	---	---

Electronic Control Unit

	(MS, CS)	1595-MK-13-ECU
	(RS, RT, FS)	1595-CC-2-ECU
	(P, R, M) 2021	1595-CC-2-2021
	(PSL, PTL)	1595-CC-3-ECU
	(HS)	1595-CH-13-0000
	(ES, ET, PS, PT)	1595-CC-13-ECU17

Control Options Continued

Left Thumb Throttle

Limited to 15 km/h (Klick Hybrid Standard only)

1595-CH-28-0001



Switch control	EBS Button		
with cruise, reverse, comfort and sport mapping (PSL, PTL)		(PSL, PTL)	
1595-CC-32-0000		1595-CC-32-0001	

Reverse Gear			
(ES, ET, PS, PT)		(RS, RT, FS)	
1595-CC-31-0000		1595-CC-31-0001	

Display Options

	(PSL, PTL)	1595-CC-27-0006
	(MS, CS)	1595-MK-27-0000
	(HS)	1595-CH-27-0000
	(RS, RT)	1595-CC-27-0002
	(FS)	1595-CC-27-0003
	(ES, ET, PS, PT)	1595-CC-27-0000
	(P, R, M)	1595-CC-27-0007

Accessories

Klaxon Starter Kit		Installation Kit	
(*Installation files NOT included: write to orders@klaxon-klick. com for training) 1595-TEMP/KLA004		Gold 1595-KIT-GOLD	
Shopping Bag Frame		Shopping Bag (Requires Shopping Bag Frame)	
1595-T-SHOPBAG		Left/Right Side (1 piece) 1595-SHOPBAG-0000	
			Front Side 1595-SHOPBAG-0001
			
Unit Control's Bag		Flybag	
(not HS) 1595-CC-90-0010		8" for Klick Mini 1595-FLYBAG01	
			14" for Klick Race and Power 1595-FLYBAG00
			
Foldable Stand (2019 version)			
8" and 14" wheel 1595-CC-CAV-0003		20" wheel 1595-CH-CAV-0001	
Phone Holder Handlebar Mount		Bell	
With Klaxon Logo 1595-PHONE-HOLDER		1595-209010	
Single Rear View Mirror		Quick Release Locking Pin with Ring Handle	
(Excluding HS) 1595-SPE-0000		(1 piece) 1595-DF-34-0001	

Accessories Continued

LED Light (Includes Battery and Mount)

	Front (White) Standard	1595-LED-WHITE-0000
	Front (White) (RS, RT, FS, PSL, PTL)	1595-CC-83-KIT
	Front (White) (P, R, M)	1595-CC-83-KIT-2021
	Back (Red)	1595-LED-RED-0000

Basket		Basket Linking Kit	
	10L 1595-BASKET-0000		(ES, PSL, RS, FS, P, R, M) 1595-KIT-TEL-0000
	16L 1595-BASKET-0001		(PS, PT) 1595-KIT-TEL-0001

Additional Traction Weight Kit

	(1 Piece) 1595-CC-93-0000	(2 pieces) 1595-CC-93-0000-2
---	------------------------------	---------------------------------



Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY: