



Quickie M6™
May 2025
EIF1

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓ Option # HCPCS Definition ✓ Option # HCPCS Definition

1. FRAME

Base Model

☐ 115WM21 K0007 Quickie M6

2. SIZE / SEATING

Frame Width

☐ 115W13 E2201 22"
☐ 115W14 E2202 24"
☐ 115W15 E2202 26"
☐ 115W16 K0108 28"
☐ 115W17 K0108 30"

Frame Length

☐ 115FL9 Regular (18"-19" depth)
☐ 115FL18 Long (18"-22" depth)

Seat Depth

☐ 115SD11 18"
☐ 115SD12 19"
☐ 115SD13 E2203 20" (Long frame only)
☐ 115SD14 E2203 21" (Long frame only)
☐ 115SD33 E2204 22" (Long frame only)

Seat Height

☐ 115HT65 17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster
☐ 115HT11 18.5" Only available w/ 5" casters
☐ 115HT14 20"

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth

☐ 115JC20 JAY Bundled Cushion
Cushion Name or Part # _____
Width: _____ Depth: _____

3. CASTERS

Casters

☐ 115CW6 5" Soft Roll
☐ 115CW36 E2214 8" x 2" Pneumatic
E2215
☐ 115CW37 E2214 8" x 2" Pneumatic w/Airless Insert
K0108

Seat Height Matrix

Caster	Front/Rear
5"	17" (Different Axle & Fork)
5"	18.5"
8"	20"

4. REAR WHEELS

Rear Wheel Type

☐ 115RW2 24" Spoke Wheel
☐ 115RW29 24" Spinergy Wheel

Rear Wheel Tire

☐ 115RT1 E2211 Pneumatic Tire
E2212
☐ 115RT2 E2211 Pneumatic Airless Insert Tire
E2213
☐ 115RT25 E2218 Primo Orion Poly

Handrim

☐ 115HR1 Aluminum Anodized
☐ 115HR4 Plastic Coated
☐ 115HR7 Omit Handrim

Adjustable Axle Plate

☐ 115AL2 Standard (for 18.5", 20" Seat Ht.)
☐ 115AL28 Low (only for 17" Seat Ht.)

Axle

☐ 115AX1 Quick Release

Rear Tire & Wheel Matrix

Tire	Rear Wheel	
	Spoke wheel	Spinergy wheel
Pneumatic	24	24
Pneumatic w/ Airless Insert	24	24
Primo Orion Poly	24	24

5. WHEEL LOCKS			
Wheel Lock		Extension Handle	
☐	115WL1		High Mount Push Wheel Lock
☐	115WL12		Do Not Mount Wheel Lock
☐	115WL10		Omit Wheel Lock
☐	115EH1	E0961	6" Pair

6. HANGERS AND FOOTREST			
Front Hangers		Extension Tubes	
☐	115H1		60° Swing Away
☐	115H146		70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>
☐	115H132		80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>
☐	115H41	E0990	ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>

Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
☐	115F1		Composite
☐	115F25	K0040	Adult Angle Adj
☐	115F5		Aluminum
Footplate Options			
☐	115HL1	E0951	Heel Loop

Lower Leg Length Range Matrix				
Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"
	Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"
	Long	16.00" - 16.50"	16.00" - 17.75"	NA
70°	No Ext Tube	NA	9.75" - 13.50"	NA
80°	No Ext Tube	NA	10.75" - 14.75"	NA
ELR	Short	NA	NA	13.00" - 17.00"
	Medium	NA	NA	14.25" - 17.00"
<i>All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle</i>				

7. BACKREST			
Backrest		Back Height Preset	
☐	115BT34		8° w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i>
☐	115B91		16-3/4"
☐	115B96		17-7/8"
☐	115B39		19"
☐	115B98		20-1/8"

8. ARMREST			
Armrest			
☐	115AA33	E0973	Dual Post Height Adj - Desk 10"
☐	115AA34	E0973	Dual Post Height Adj - Full 14"

9. KEY OPTIONS & ACCESSORIES			
Rear Anti-Tip		Positioning Belt	
☐	115AT1	E0971	Rear Anti-Tip
☐	115AT3		Omit Rear Anti-Tip
Miscellaneous			
☐	115A83		Tool Kit
☐	115A115	E0978	Belt 2" Aircraft Buckle Long
☐	115A116	E0978	Belt 2" Aircraft Buckle Padded Long
☐	115A142	E0978	Belt 1.5" Auto Buckle

10. COLOR

Standard Frame Colors			
☒	115C214		Black (Gloss)
☐	115C215		Black (Matte)
☐	115C211		Celestial Sapphire
☐	115C212		Desert Sage
☐	115C213		Desert Tan
☐	115C216		Metallic Shadow
☐	115C217		Midnight Blue
☐	115C218		Mystic Purple
☐	115C219		Petal Pink
☐	115C220		Pink Power
☐	115C221		Red Pepper
☐	115C222		Royal Blue
☐	115C223		Salute Green
☐	115C224		Scarlet Noir
☐	115C225		Silver Ore
☐	115C226		Solar Orange
☐	115C227		Starlight White
☐	115C228		Storm Gray
☐	115C229		Turbo Green
☐	115C230		Velocity Yellow

Personalized Frame Color			
☐	115C112		Desert Camouflage
☐	115C43		Mossy Oak Camouflage
☐	115C42		Stars & Stripes
☐	115C108		Zebra (<i>Choose above base color Std with White unless above base paint is chosen NA w/ Black, Matte Black</i>)