



Quickie Nitrum

July 2025

EIR12

NITRUM

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

ST: _____

ZIP: _____

Dealer Phone: () _____

Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: () _____

Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	✓	Option #	HCPCS	Definition
---	----------	-------	------------	---	----------	-------	------------

1. FRAME

Base Model

Weight capacity 275 lbs for Open frame. 300 lbs for Hybrid dual tube frame.

☐	112WM101	K0005	Quickie Nitrum - Open frame Use E1235 for sizes 14" or less
☐	112WM102	K0005	Quickie Nitrum Hybrid - Dual tube Use E1235 for sizes 14" or less.

Backrest Frame Type

Active Rigid - welded must select Active Axle plate and Active Fixed backrest.

☐	112FT68	Adjustable
☐	112FT69	Active Rigid - Welded

Frame Angle

See **FIGURE 1** Measure from front frame to floor.

☐	112FA1	75°
☐	112FA8	80°
☐	112FA2	85°

Frame Inset Per Side

See **FIGURE 2** Frame Inset is the distance from the outside of the seat tube to the outside of the footrest.

☐	112FF1	0"
☐	112FF2	1" (N/A w/ 12" wide)
☐	112FF3	2" (N/A with 12" & 13" wide)

Frame Length

See **FIGURE 3** measurement for set up in the APPENDIX page 7.

☐	112FS9	Standard
☐	112FS34	Classic

* **Standard Frame** - Seat Sling to front frame bend angle is 1.5" shorter than the classic position.
Classic Frame - Traditional seat sling to front frame bend angle dimensions. Front of back cane post to front edge of seat sling will remain the same for both Standard and Classic.

Transit

WC-19 and ISO 7176 standard compliance. N/A with Active Rigid.

☐	112TO45	K0108	Transit System
--------------------------	---------	-------	----------------

Nitrum - Open frame



Nitrum Hybrid - Dual tube frame



Power Assist

To order any Empulse Power Assist product please see the Order Form for Part #'s

Option	HCPCS	DESCRIPTION	RETAIL
☐	209MW91	E0986 Empulse R90 Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 2.5Ahr Removable Battery	
☐	209MW92	E0986 Empulse R90 Range Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 5.0Ahr Removable Battery	
☐	147MW95	N/A Empulse F35 Includes: F35 Stick, Subframe mount, Frame clamps kit, and Charger	



medifab

0800 543 343
solutions@medifab.com
medifab.com

✓ Option # HCPCS Definition

✓ Option # HCPCS Definition

2. SIZE / SEATING

Depth Adjustable Backrest (Optional)

See **Figure 4** Depth Adjust Backrest extends the frame rear of the seat sling, select a factory setting. N/A with 19"-20" seat depth. N/A with Active Rigid frame.

☐	112DP02	0 Setting Backrest at standard position, field adjustable to -1" or -2"
☐	112DP03	-1 Setting Backrest -1" offset, field adjustable to flush or -2"
☐	112DP04	-2 Setting Backrest -2" offset, field adjustable to flush or -1"




Figure 4

0 position

-1" position

-2" position

Seat Width

See **FIGURE 2** Measure from outside of seat tube to the outside of opposite seat tube.

☐	112S3	12"	☐	112S8	17"
☐	112S4	13"	☐	112S9	18"
☐	112S5	14"	☐	112S10	19"
☐	112S6	15"	☐	112S48	20"
☐	112S7	16"			(E2201)

Seat Depth

See **FIGURE 1** Measure from front of back post to front edge of seat sling. Seat sling is 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

☐	112SD5	12"	☐	112SD10	17"
☐	112SD6	13"	☐	112SD11	18"
☐	112SD7	14"	☐	112SD12	19"
☐	112SD8	15"	☐	112SD13	20"
☐	112SD9	16"			(E2203)

Figure 1





Figure 2



Seating Systems

☐	112OP01	Standard Seat Sling
☐	112OP14	Tension Adjustable seat sling
☐	112OP04	Omit Seat Sling

Seat Sling Options

Utility & catheter bags are sewn at the bottom of the upholstery. Catheter bag N/A in 12" depth. Utility bag and Catheter bag combination N/A at 15" width or less.

☐	1128A65	1 Utility Bag
☐	1128A66	1 Catheter bag

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth.

☐	112JC20	JAY Bundled Cushion Cushion Name or Part #
		Width: _____ Depth: _____

3. CASTERS / SEAT HEIGHT

Front Seat Height

See **FIGURE 5** Measure from the front of the seat rail to the floor.

☐	112SH15	16.5"	☐	112SH20	18.5"
☐	112SH16	17.0"	☐	112SH21	19.0"
☐	112SH17	17.5"	☐	112SH23	19.5"
☐	112SH18	18.0"	☐	112SH24	20.0"
			☐	112SH26	20.5"
			☐	112SH27	21.0"




Figure 5

Caster Wheels

☐	112CW1	3" Micro
☐	112CW3	4" Low Profile Polyurethane
☐	112CW31	4" Micro Lighted
☐	112CW72	4" x 1" Aluminum Poly
☐	112CW44	E2219 4" x 1.25 Semi Pneumatic
☐	112CW40	4" x 1.25 Spinergy CX6 cast carbon NEW!
☐	112CW5	5" Low Profile Polyurethane
☐	112CW51	5" Micro Lighted
☐	112CW6	E2219 5" Soft Roll (Grey) -N/A w/Frog Legs
☐	112CW65	5" x 1" Aluminum Poly
☐	112CW7	E2219 5" x 1.5" Semi Pneumatic
☐	112CW34	6" Poly
☐	112CW33	E2214 6" Pneumatic
☐	112CW35	E2215 6" x 1.5" Semi Pneumatic

Caster Color Wheels option

Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

☐	112CW64	4" x 1.5" Aluminum Soft Roll
☐	112CW66	5" x 1.5" Aluminum Soft Roll
☐	112CW68	6" x 1.5" Aluminum Soft Roll
	Colors:	☐ Match anodize color package
		☐ Match frame color
		☐ Black anodize caster only

Caster Housing Position

See **FIGURE 7** in page 7 for example. Options available on 0" inset only.

☐	112LH04	Position Setting 1 (Narrow)
☐	112LH05	Position Setting 2 (Wide)


✓	Option #	HCPCS	Definition	✓	Option #	HCPCS	Definition																																																																																																								
4. REAR WHEELS / SEAT HEIGHT																																																																																																															
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="6">Rear Seat Height</th> </tr> </thead> <tbody> <tr> <td colspan="6">See RSH matrix</td> </tr> <tr> <td>☐</td> <td>112RH11</td> <td>13.0"</td> <td></td> <td>☐</td> <td>112RH19</td> <td>17.0"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH12</td> <td>13.5"</td> <td></td> <td>☐</td> <td>112RH26</td> <td>17.5"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH13</td> <td>14.0"</td> <td></td> <td>☐</td> <td>112RH20</td> <td>18.0"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH14</td> <td>14.5"</td> <td></td> <td>☐</td> <td>112RH27</td> <td>18.5"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH15</td> <td>15.0"</td> <td></td> <td>☐</td> <td>112RH21</td> <td>19.0"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH16</td> <td>15.5"</td> <td></td> <td>☐</td> <td>112RH37</td> <td>19.5"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH17</td> <td>16.0"</td> <td></td> <td>☐</td> <td>112RH22</td> <td>20.0"</td> <td></td> </tr> <tr> <td>☐</td> <td>112RH18</td> <td>16.5"</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>								Rear Seat Height						See RSH matrix						☐	112RH11	13.0"		☐	112RH19	17.0"		☐	112RH12	13.5"		☐	112RH26	17.5"		☐	112RH13	14.0"		☐	112RH20	18.0"		☐	112RH14	14.5"		☐	112RH27	18.5"		☐	112RH15	15.0"		☐	112RH21	19.0"		☐	112RH16	15.5"		☐	112RH37	19.5"		☐	112RH17	16.0"		☐	112RH22	20.0"		☐	112RH18	16.5"																																	
Rear Seat Height																																																																																																															
See RSH matrix																																																																																																															
☐	112RH11	13.0"		☐	112RH19	17.0"																																																																																																									
☐	112RH12	13.5"		☐	112RH26	17.5"																																																																																																									
☐	112RH13	14.0"		☐	112RH20	18.0"																																																																																																									
☐	112RH14	14.5"		☐	112RH27	18.5"																																																																																																									
☐	112RH15	15.0"		☐	112RH21	19.0"																																																																																																									
☐	112RH16	15.5"		☐	112RH37	19.5"																																																																																																									
☐	112RH17	16.0"		☐	112RH22	20.0"																																																																																																									
☐	112RH18	16.5"																																																																																																													
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Axle</th> </tr> </thead> <tbody> <tr> <td>☐</td> <td>112AX1</td> <td></td> <td>Quick Release Steel</td> </tr> <tr> <td>☐</td> <td>112AX13</td> <td></td> <td>Quad Quick Release Steel Only available with Mountain Whls.</td> </tr> </tbody> </table>								Axle				☐	112AX1		Quick Release Steel	☐	112AX13		Quad Quick Release Steel Only available with Mountain Whls.																																																																																												
Axle																																																																																																															
☐	112AX1		Quick Release Steel																																																																																																												
☐	112AX13		Quad Quick Release Steel Only available with Mountain Whls.																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Axle Plate</th> </tr> </thead> <tbody> <tr> <td colspan="4">Extended COG N/A w/ Xtender Power Assist. 4-link suspension N/A with Hybrid frame type.</td> </tr> <tr> <td>☐</td> <td></td> <td>Adjustable COG (1"-4")</td> <td></td> </tr> <tr> <td>☐</td> <td></td> <td>Adjustable Extended COG (0"-3")</td> <td></td> </tr> <tr> <td>☐</td> <td></td> <td>4-Link Suspension axle plate</td> <td></td> </tr> <tr> <td>☐</td> <td></td> <td>Active Rigid (.5" to 4") Adjustment range 1/2" before and after in 1/4" increments except the .5" and 4" positions which are maxed.</td> <td></td> </tr> </tbody> </table>								Axle Plate				Extended COG N/A w/ Xtender Power Assist. 4-link suspension N/A with Hybrid frame type.				☐		Adjustable COG (1"-4")		☐		Adjustable Extended COG (0"-3")		☐		4-Link Suspension axle plate		☐		Active Rigid (.5" to 4") Adjustment range 1/2" before and after in 1/4" increments except the .5" and 4" positions which are maxed.																																																																																	
Axle Plate																																																																																																															
Extended COG N/A w/ Xtender Power Assist. 4-link suspension N/A with Hybrid frame type.																																																																																																															
☐		Adjustable COG (1"-4")																																																																																																													
☐		Adjustable Extended COG (0"-3")																																																																																																													
☐		4-Link Suspension axle plate																																																																																																													
☐		Active Rigid (.5" to 4") Adjustment range 1/2" before and after in 1/4" increments except the .5" and 4" positions which are maxed.																																																																																																													
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Rear Suspension (Optional)</th> </tr> </thead> <tbody> <tr> <td colspan="4">4-Link Rear Suspension COG (0"-3") Single Post only available up to 2" COG. Plastic Side Guards only available up to 2.5" COG. N/A w/Xtender. Max weight capacity of 265 lbs.</td> </tr> <tr> <td>☐</td> <td>112SV08</td> <td></td> <td>Weight from 0 lbs to 175 lbs</td> </tr> <tr> <td>☐</td> <td>112SV09</td> <td></td> <td>Weight from 176 lbs. to 265 lbs</td> </tr> </tbody> </table>								Rear Suspension (Optional)				4-Link Rear Suspension COG (0"-3") Single Post only available up to 2" COG. Plastic Side Guards only available up to 2.5" COG. N/A w/Xtender. Max weight capacity of 265 lbs.				☐	112SV08		Weight from 0 lbs to 175 lbs	☐	112SV09		Weight from 176 lbs. to 265 lbs																																																																																								
Rear Suspension (Optional)																																																																																																															
4-Link Rear Suspension COG (0"-3") Single Post only available up to 2" COG. Plastic Side Guards only available up to 2.5" COG. N/A w/Xtender. Max weight capacity of 265 lbs.																																																																																																															
☐	112SV08		Weight from 0 lbs to 175 lbs																																																																																																												
☐	112SV09		Weight from 176 lbs. to 265 lbs																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Rear Wheel</th> </tr> </thead> <tbody> <tr> <td colspan="4">Ultra Lightweight only available w/ Aluminum Handrim. Max weight capacity on Proton wheels is 220 lbs, Spinergy CLX is 265 lbs.</td> </tr> <tr><td>☐</td><td>112RW103</td><td></td><td>20" Mag Wheel 5 Spoke</td></tr> <tr><td>☐</td><td>112RW104</td><td></td><td>22" Mag Wheel 5 Spoke</td></tr> <tr><td>☐</td><td>112RW102</td><td></td><td>24" Mag Wheel 5 Spoke</td></tr> <tr><td>☐</td><td>112RW90</td><td></td><td>20" Lite Spoke</td></tr> <tr><td>☐</td><td>112RW91</td><td></td><td>22" Lite Spoke</td></tr> <tr><td>☐</td><td>112RW92</td><td></td><td>24" Lite Spoke</td></tr> <tr><td>☐</td><td>112RW93</td><td></td><td>25" Lite Spoke</td></tr> <tr><td>☐</td><td>112RW94</td><td></td><td>26" Lite Spoke</td></tr> <tr><td>☐</td><td>112RW06</td><td></td><td>22" Ultra Lightweight Spoke Wheel</td></tr> <tr><td>☐</td><td>112RW10</td><td></td><td>24" Ultra Lightweight Spoke Wheel</td></tr> <tr><td>☐</td><td>112RW16</td><td></td><td>25" Ultra Lightweight Spoke Wheel</td></tr> <tr><td>☐</td><td>112RW28</td><td></td><td>22" Spinergy SPOX Black Spokes</td></tr> <tr><td>☐</td><td>112RW95</td><td></td><td>24" Spinergy SPOX Black Spokes</td></tr> <tr><td>☐</td><td>112RW96</td><td></td><td>25" Spinergy SPOX Black Spokes</td></tr> <tr><td>☐</td><td>112RW97</td><td></td><td>26" Spinergy SPOX Black Spokes</td></tr> <tr><td>☐</td><td>112RW32</td><td></td><td>22" Spinergy LX Black Spokes</td></tr> <tr><td>☐</td><td>112RW98</td><td></td><td>24" Spinergy LX Black Spokes</td></tr> <tr><td>☐</td><td>112RW99</td><td></td><td>25" Spinergy LX Black Spokes</td></tr> <tr><td>☐</td><td>112RW40</td><td></td><td>26" Spinergy LX Black Spokes</td></tr> <tr><td>☐</td><td>112RW152</td><td></td><td>24" Spinergy CLX Carbon Wheel</td></tr> <tr><td>☐</td><td>112RW153</td><td></td><td>25" Spinergy CLX Carbon Wheel</td></tr> <tr><td>☐</td><td>112RW148</td><td></td><td>24" Proton Wheel</td></tr> <tr><td>☐</td><td>112RW149</td><td></td><td>25" Proton Wheel</td></tr> <tr><td>☐</td><td>112RW50</td><td></td><td>24" Mountain</td></tr> </tbody> </table>								Rear Wheel				Ultra Lightweight only available w/ Aluminum Handrim. Max weight capacity on Proton wheels is 220 lbs, Spinergy CLX is 265 lbs.				☐	112RW103		20" Mag Wheel 5 Spoke	☐	112RW104		22" Mag Wheel 5 Spoke	☐	112RW102		24" Mag Wheel 5 Spoke	☐	112RW90		20" Lite Spoke	☐	112RW91		22" Lite Spoke	☐	112RW92		24" Lite Spoke	☐	112RW93		25" Lite Spoke	☐	112RW94		26" Lite Spoke	☐	112RW06		22" Ultra Lightweight Spoke Wheel	☐	112RW10		24" Ultra Lightweight Spoke Wheel	☐	112RW16		25" Ultra Lightweight Spoke Wheel	☐	112RW28		22" Spinergy SPOX Black Spokes	☐	112RW95		24" Spinergy SPOX Black Spokes	☐	112RW96		25" Spinergy SPOX Black Spokes	☐	112RW97		26" Spinergy SPOX Black Spokes	☐	112RW32		22" Spinergy LX Black Spokes	☐	112RW98		24" Spinergy LX Black Spokes	☐	112RW99		25" Spinergy LX Black Spokes	☐	112RW40		26" Spinergy LX Black Spokes	☐	112RW152		24" Spinergy CLX Carbon Wheel	☐	112RW153		25" Spinergy CLX Carbon Wheel	☐	112RW148		24" Proton Wheel	☐	112RW149		25" Proton Wheel	☐	112RW50		24" Mountain
Rear Wheel																																																																																																															
Ultra Lightweight only available w/ Aluminum Handrim. Max weight capacity on Proton wheels is 220 lbs, Spinergy CLX is 265 lbs.																																																																																																															
☐	112RW103		20" Mag Wheel 5 Spoke																																																																																																												
☐	112RW104		22" Mag Wheel 5 Spoke																																																																																																												
☐	112RW102		24" Mag Wheel 5 Spoke																																																																																																												
☐	112RW90		20" Lite Spoke																																																																																																												
☐	112RW91		22" Lite Spoke																																																																																																												
☐	112RW92		24" Lite Spoke																																																																																																												
☐	112RW93		25" Lite Spoke																																																																																																												
☐	112RW94		26" Lite Spoke																																																																																																												
☐	112RW06		22" Ultra Lightweight Spoke Wheel																																																																																																												
☐	112RW10		24" Ultra Lightweight Spoke Wheel																																																																																																												
☐	112RW16		25" Ultra Lightweight Spoke Wheel																																																																																																												
☐	112RW28		22" Spinergy SPOX Black Spokes																																																																																																												
☐	112RW95		24" Spinergy SPOX Black Spokes																																																																																																												
☐	112RW96		25" Spinergy SPOX Black Spokes																																																																																																												
☐	112RW97		26" Spinergy SPOX Black Spokes																																																																																																												
☐	112RW32		22" Spinergy LX Black Spokes																																																																																																												
☐	112RW98		24" Spinergy LX Black Spokes																																																																																																												
☐	112RW99		25" Spinergy LX Black Spokes																																																																																																												
☐	112RW40		26" Spinergy LX Black Spokes																																																																																																												
☐	112RW152		24" Spinergy CLX Carbon Wheel																																																																																																												
☐	112RW153		25" Spinergy CLX Carbon Wheel																																																																																																												
☐	112RW148		24" Proton Wheel																																																																																																												
☐	112RW149		25" Proton Wheel																																																																																																												
☐	112RW50		24" Mountain																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Spinergy Spoke Colors (Optional)</th> </tr> </thead> <tbody> <tr> <td colspan="4">Please select one of the following Spinergy wheel options first. Color spokes N/A with CLX carbon wheels and Spinergy Flexrim .</td> </tr> <tr> <td colspan="4">Spinergy Spox spokes color:</td> </tr> <tr> <td>☐ Black</td> <td>☐ Blue</td> <td>☐ Red</td> <td>☐ White ☐ Yellow</td> </tr> <tr> <td colspan="4">Spinergy LX spokes color:</td> </tr> <tr> <td>☐ Black</td> <td>☐ Blue</td> <td>☐ Green</td> <td>☐ Orange ☐ Pink</td> </tr> <tr> <td>☐ Purple</td> <td>☐ Red</td> <td>☐ White</td> <td>☐ Yellow</td> </tr> </tbody> </table>								Spinergy Spoke Colors (Optional)				Please select one of the following Spinergy wheel options first. Color spokes N/A with CLX carbon wheels and Spinergy Flexrim .				Spinergy Spox spokes color:				☐ Black	☐ Blue	☐ Red	☐ White ☐ Yellow	Spinergy LX spokes color:				☐ Black	☐ Blue	☐ Green	☐ Orange ☐ Pink	☐ Purple	☐ Red	☐ White	☐ Yellow																																																																												
Spinergy Spoke Colors (Optional)																																																																																																															
Please select one of the following Spinergy wheel options first. Color spokes N/A with CLX carbon wheels and Spinergy Flexrim .																																																																																																															
Spinergy Spox spokes color:																																																																																																															
☐ Black	☐ Blue	☐ Red	☐ White ☐ Yellow																																																																																																												
Spinergy LX spokes color:																																																																																																															
☐ Black	☐ Blue	☐ Green	☐ Orange ☐ Pink																																																																																																												
☐ Purple	☐ Red	☐ White	☐ Yellow																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="5">Rear Wheel Tire</th> </tr> </thead> <tbody> <tr> <td>☐</td> <td>112RT1</td> <td></td> <td>Pneumatic</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT2</td> <td>E2211 E2213</td> <td>Pneumatic Airless Insert N/A with Xtender.</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT24</td> <td>E2211 E2213</td> <td>Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT3</td> <td></td> <td>Full Poly Tire</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT38</td> <td>E2218</td> <td>SHOX Performance Poly</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT25</td> <td>E2218</td> <td>Primo Orion Poly</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT39</td> <td>E2211 E2212</td> <td>High-Pressure Clincher Black</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT21</td> <td>E2211 E2212</td> <td>Schwalbe Marathon Plus Evolution</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT42</td> <td>E2211 E2212</td> <td>Schwalbe One</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT40</td> <td>E2211 E2212</td> <td>Schwalbe Down Town Two Grip Only available in 24".</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT41</td> <td>E2211 E2212</td> <td>Schwalbe Right Run Only available in 24".</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT27</td> <td>E2211 E2212</td> <td>Knobby V-TRAK by Kenda</td> <td></td> </tr> <tr> <td>☐</td> <td>112RT12</td> <td>E2211 E2212</td> <td>Mountain Only available in 24".</td> <td></td> </tr> </tbody> </table>								Rear Wheel Tire					☐	112RT1		Pneumatic		☐	112RT2	E2211 E2213	Pneumatic Airless Insert N/A with Xtender.		☐	112RT24	E2211 E2213	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material		☐	112RT3		Full Poly Tire		☐	112RT38	E2218	SHOX Performance Poly		☐	112RT25	E2218	Primo Orion Poly		☐	112RT39	E2211 E2212	High-Pressure Clincher Black		☐	112RT21	E2211 E2212	Schwalbe Marathon Plus Evolution		☐	112RT42	E2211 E2212	Schwalbe One		☐	112RT40	E2211 E2212	Schwalbe Down Town Two Grip Only available in 24".		☐	112RT41	E2211 E2212	Schwalbe Right Run Only available in 24".		☐	112RT27	E2211 E2212	Knobby V-TRAK by Kenda		☐	112RT12	E2211 E2212	Mountain Only available in 24".																																			
Rear Wheel Tire																																																																																																															
☐	112RT1		Pneumatic																																																																																																												
☐	112RT2	E2211 E2213	Pneumatic Airless Insert N/A with Xtender.																																																																																																												
☐	112RT24	E2211 E2213	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material																																																																																																												
☐	112RT3		Full Poly Tire																																																																																																												
☐	112RT38	E2218	SHOX Performance Poly																																																																																																												
☐	112RT25	E2218	Primo Orion Poly																																																																																																												
☐	112RT39	E2211 E2212	High-Pressure Clincher Black																																																																																																												
☐	112RT21	E2211 E2212	Schwalbe Marathon Plus Evolution																																																																																																												
☐	112RT42	E2211 E2212	Schwalbe One																																																																																																												
☐	112RT40	E2211 E2212	Schwalbe Down Town Two Grip Only available in 24".																																																																																																												
☐	112RT41	E2211 E2212	Schwalbe Right Run Only available in 24".																																																																																																												
☐	112RT27	E2211 E2212	Knobby V-TRAK by Kenda																																																																																																												
☐	112RT12	E2211 E2212	Mountain Only available in 24".																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Xtender Power Assist</th> </tr> </thead> <tbody> <tr> <td colspan="4">RSH 16.5"-19" with 22" wheels, RSH 17.5"-20" with 24" wheels 0 camber, pneumatic tire. 22" & 24" wheels only. Anti-Tips no charge. N/A w/ Transit and 4-link Suspension. Dual Drive mode available with B4ME. Max weight capacity of 265 lbs. N/A with Active Rigid frame.</td> </tr> <tr> <td>☐</td> <td>112XT6</td> <td>E0986</td> <td>22" Xtender Power Assist JWX-2 NICKEL</td> </tr> <tr> <td>☐</td> <td>112XT7</td> <td>E0986</td> <td>22" Xtender Power Assist JWX-2 LITHIUM</td> </tr> <tr> <td>☐</td> <td>112XT4</td> <td>E0986</td> <td>24" Xtender Power Assist JWX-2 NICKEL</td> </tr> <tr> <td>☐</td> <td>112XT5</td> <td>E0986</td> <td>24" Xtender Power Assist JWX-2 LITHIUM</td> </tr> <tr> <td colspan="4">Please select one of the following handrim options</td> </tr> <tr> <td>☐</td> <td>112XH5</td> <td></td> <td>Plastic Coated</td> </tr> <tr> <td>☐</td> <td>112XH6</td> <td></td> <td>Stainless Steel</td> </tr> </tbody> </table>								Xtender Power Assist				RSH 16.5"-19" with 22" wheels, RSH 17.5"-20" with 24" wheels 0 camber, pneumatic tire. 22" & 24" wheels only. Anti-Tips no charge. N/A w/ Transit and 4-link Suspension. Dual Drive mode available with B4ME. Max weight capacity of 265 lbs. N/A with Active Rigid frame.				☐	112XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	☐	112XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	☐	112XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	☐	112XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	Please select one of the following handrim options				☐	112XH5		Plastic Coated	☐	112XH6		Stainless Steel																																																																				
Xtender Power Assist																																																																																																															
RSH 16.5"-19" with 22" wheels, RSH 17.5"-20" with 24" wheels 0 camber, pneumatic tire. 22" & 24" wheels only. Anti-Tips no charge. N/A w/ Transit and 4-link Suspension. Dual Drive mode available with B4ME. Max weight capacity of 265 lbs. N/A with Active Rigid frame.																																																																																																															
☐	112XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL																																																																																																												
☐	112XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM																																																																																																												
☐	112XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL																																																																																																												
☐	112XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM																																																																																																												
Please select one of the following handrim options																																																																																																															
☐	112XH5		Plastic Coated																																																																																																												
☐	112XH6		Stainless Steel																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Xtender - Other Options</th> </tr> </thead> <tbody> <tr> <td>☐</td> <td>112XA10</td> <td></td> <td>Spare NiMH Power Pack - NICKEL</td> </tr> <tr> <td>☐</td> <td>112XA11</td> <td></td> <td>Spare Li-ion Power Pack - LITHIUM</td> </tr> </tbody> </table>								Xtender - Other Options				☐	112XA10		Spare NiMH Power Pack - NICKEL	☐	112XA11		Spare Li-ion Power Pack - LITHIUM																																																																																												
Xtender - Other Options																																																																																																															
☐	112XA10		Spare NiMH Power Pack - NICKEL																																																																																																												
☐	112XA11		Spare Li-ion Power Pack - LITHIUM																																																																																																												
<table border="1" style="width: 100%;"> <thead> <tr> <th colspan="4">Xtender Wheel Options</th> </tr> </thead> <tbody> <tr> <td>☐</td> <td>112XA8</td> <td></td> <td>SPARE 22" Lite Spoke Pneu Alum</td> </tr> <tr> <td>☐</td> <td>112XA9</td> <td>E2211 E2213</td> <td>SPARE 22" Lite Spoke Pneu Plastic</td> </tr> <tr> <td>☐</td> <td>112XA2</td> <td></td> <td>SPARE 24" Lite Spoke Pneu Alum</td> </tr> <tr> <td>☐</td> <td>112XA4</td> <td>E2211 E2213</td> <td>SPARE 24" Lite Spoke Pneu Plastic</td> </tr> </tbody> </table>								Xtender Wheel Options				☐	112XA8		SPARE 22" Lite Spoke Pneu Alum	☐	112XA9	E2211 E2213	SPARE 22" Lite Spoke Pneu Plastic	☐	112XA2		SPARE 24" Lite Spoke Pneu Alum	☐	112XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic																																																																																				
Xtender Wheel Options																																																																																																															
☐	112XA8		SPARE 22" Lite Spoke Pneu Alum																																																																																																												
☐	112XA9	E2211 E2213	SPARE 22" Lite Spoke Pneu Plastic																																																																																																												
☐	112XA2		SPARE 24" Lite Spoke Pneu Alum																																																																																																												
☐	112XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic																																																																																																												

✓	Option #	HPCPS	Definition	✓	Option #	HPCPS	Definition
7. BACKREST							
Adjustable Back Type Options				Active Rigid Back Type Options			
☐ 112BT96 Folding w/ Double Lock Down Angle Adjustable				☐ 112BT97 Active Fixed Backrest			
☐ 112BT122 Freestyle Backrest Folding w/Double Lock Down				Setting ☐ -4 ☐ -2 ☐ 0 ☐ +2 ☐ +4			
Choose Option 1 or 2							
1. Aluminum Height Adjustable Backrest				OR	2. Ultra Lightweight Aluminum Fixed Backrest		
☐ 112B159 Setting 10" - 12" Back Height ☐ 10" ☐ 11" ☐ 12"					Select back height. 8", 9", or 10" only available with omit push handles.		
☐ 112B139 Setting 12" - 14" Back Height ☐ 12" ☐ 13" ☐ 14"					☐ 112BB21 8" ☐ 112BB17 14"		
☐ 112B138 Setting 14" - 16" Back Height ☐ 14" ☐ 15" ☐ 16"					☐ 112BB20 9" ☐ 112BB18 15"		
☐ 112B137 Setting 16" - 20" Back Height ☐ 16 ☐ 17 ☐ 18 ☐ 19 ☐ 20					☐ 112BB13 10" ☐ 112BB19 16"		
					☐ 112BB14 11" ☐ 112BB10 17"		
					☐ 112BB15 12" ☐ 112BB11 18"		
					☐ 112BB16 13" ☐ 112BB12 19"		
Back release options							
Twist-lock bar N/A with Freestyle back or Height Adjustable backrest. N/A with 12"-13" width.							
Back release options N/A with Active Rigid backrest							
☐ 112OD24 Push-bar release							
☐ 112OD25 Twist-lock release							
Towel Bar							
☐ 112QT01 Standard Style - 4" Offset							
☐ 112QT11 Narrow Profile - 3" Offset <i>Limited w/Solid Backs, Not suggested with 19" & 20" wide.</i>							
Push Handles							
Integral push handles N/A with Active Rigid backrest.							
☐ 112PH6 Integral Push Handles							
☐ 112PH5 Fold Down Push Handles N/A w/ 10"-12" adj height & 8"-12" fixed height. N/A w/ standard EXO upholstery.							
☐ 112PH2 Omit Push Handles							
Backrest Options							
Telescoping Push Handles have "Double D" tubing to eliminate rotation. Available with Omit push handle and EXO / EXO Pro Tension Adjustable upholstery only.							
☐ 112A192 Telescoping Push Handles							
☐ 112A195 Telescoping Push Handles - Tall							
JAY J3 Back							
☐ 112JB01 J3 Back To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s. N/A w/Freestyleback.							

Freestyle Backrest Option 3			
Accessories such as head supports & lateral supports are not available with Freestyle backrest. The Freestyle back comes standard with a small lumbar positioning component Available with Plastic side guards and limited ranges with Carbon fiber sideguard with fender. Chair Seat depth will be reduced by 1" with backrest.			
Backrest Type			
☐	112JB06	E2613	JAY J3 Back LT Aluminum model (J3PAXXXYYZZ)
☐	112JB07	E2613	JAY J3 Back LT Carbon Fiber model (J3CFPAXXXYYZZ)
Backrest Height			
Medium size backrest not available in Carbon Fiber.			
☐	112QF04		Small 6.5" shell hgt, 9" back height
☐	112QF05		Medium 8" shell hgt, 10.5" back height
☐	112QF06		Tall 9.5" shell hgt, 12" back height
Backrest Color			
Not applicable with carbon fiber back or Zebra color.			
☐	112GZ01		Black
☐	112GZ02		Matching Frame- Standard Colors
☐	112GZ03		Matching Frame- Personalized Color

Backrest Width			
☐	112JQ01		14" Backrest Width
☐	112JQ02		16" Backrest Width
☐	112JQ03		18" Backrest Width
☐	112JQ04		20" Backrest Width

Height Ranges	Dim "A"	Dim "B"
Small	12.0" - 14.5"	3.0" - 5.5"
Medium	12.0" - 16.0"	1.5" - 5.5"
Tall	12.0" - 17.5"	0.5" - 5.5"



✓	Option #	HCPCS	Definition	✓	Option #	HCPCS	Definition
---	----------	-------	------------	---	----------	-------	------------

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away

Some limitations when used with the Angle Adjustable Back in the rear most setting
Please select a setting, will default to 11" if no selection is made.

☐	112AA162		Padded Swing-Away Fixed height (13")	
☐	112AA18	E0973	Padded Swing-Away STD (9" to 14")	
	Setting	☐ 9"	☐ 10"	☐ 11" ☐ 12" ☐ 13" ☐ 14"

Armrest - Height Adjustable

☐	112AA11	E0973	Single Post Height Adjustable - Std Range 8" - 13". Adjust in 1" increments	
☐	112AA29	E0973	Single Post Height Adjustable - Low Range 6.5" - 10.5". Adjust in 1" increments	

Armpad - For Height Adjustable

☐	112AM23		Desk Length (10") Classic	
☐	112AM21		Desk Length (10") Contour	
☐	112AM24		Full Length (14") Classic	
☐	112AM22		Full Length (14") Contour	
☐	112AM27		Ergo Pad (12")	

Side Guards

Increases overall width of chair. Some interference may exist w/ padded SA arms at lowest settings
Not available w/Height Adjustable Arms. **Alum Flip sideguard option N/A with Twist-lock fold release, Freestyle back and Active Rigid.**

☐	112SG33	K0108	Plastic Standard	
☐	112SG29	K0108	Plastic Low	
☐	112SG35	K0108	Aluminum	
☐	112SG36	K0108	Aluminum w/ Fender	
☐	112SG12	K0108	Carbon Fiber	
☐	112SG41	K0108	Carbon Fiber Long w/LG Fender	
☐	1129A05	K0108	Flip Back Option Aluminum side guard only.	

Side Guard - Color Option

Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber options.

☐	112GE01		Black	
☐	112GE02		Matching Frame- Standard Colors	
☐	112GE03		Matching Frame- Personalized Color	

9. KEY OPTIONS AND ACCESSORIES

Anti Tips

Anti Tip Single Side N/A 12W-14W range, N/A w/ Rear suspension. Included with Xtender.

☐	112AT1	E0971	Anti Tips (Pair)	
☐	112AT14	E0971	Anti Tip Single Side - Flip Up ☐ Left ☐ Right (Facing from rear)	
☐	112AT83	E0971	Tip assist ☐ Right Mnt ☐ Left Mnt	
☐	112AT20	E0971	Sport Anti Tip / 5th wheel	

Active lifestyle options

LED lights mounted on caster arms, battery & controller on leg strap.
Leg strap included with LED lights.

☐	112A1		Backpack - Adult	
☐	112A2		Backpack - Kids	
☐	1121A01		Seat Pouch Black	
☐	1128A47		Mobile phone pocket	
☐	1129A06		Twin LED lights	
☐	112A999		Quickie Water Bottle & Holder	

Spoke Guards

Only available with Lite spoke. N/A w/ Xtender Power Assist.

☐	1124A01	K0065	Quickie Black Spoke Guard	
☐	1124A02	K0065	Quickie Clear Spoke Guard	

Positioning Belt

☐	112A117	E0978	Belt 2" Aircraft Buckle w/Clamps	
☐	112A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	
☐	112A141	E0978	Belt 1.5" Auto Buckle w/Clamps	
☐	112A804		Positioning Belt Clamps Only	

Accessories

☐	112A2100	E2207	Crutch Holder	
☐	1128A21		Tool Kit	
☐	112NP01		Gel Swing Away Pad	
☐	1129A07		Frame Protections pads	

10. COLOR

NEW!

Standard Frame Color

Handrim color option only available if Aluminum Coated colored option was selected.

			Handrim	
☐	112C214		Black (Gloss)	☐
☐	112C215		Black (Matte)	☐
☐	112C211		Celestial Sapphire	☐
☐	112C212		Desert Sage	☐
☐	112C213		Desert Tan	☐
☐	112C216		Metallic Shadow	☐
☐	112C217		Midnight Blue	☐
☐	112C218		Mystic Purple	☐
☐	112C219		Petal Pink	☐
☐	112C220		Pink Power	☐
☐	112C221		Red Pepper	☐
☐	112C222		Royal Blue	☐
☐	112C223		Salute Green	☐
☐	112C224		Scarlet Noir	☐
☐	112C225		Silver Ore	☐
☐	112C226		Solar Orange	☐
☐	112C227		Starlight White	☐
☐	112C228		Storm Gray	☐
☐	112C229		Turbo Green	☐
☐	112C230		Velocity Yellow	☐

Personalized Frame Color

☐	112C112		Desert Camouflage	
☐	112C43		Mossy Oak Camouflage	
☐	112C42		Stars & Stripes	
☐	112C108		Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Blackk, Matte Black.	

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

☐	112D16		Black	
☐	112D17		Blue	
☐	112D18		Orange	
☐	112D19		Red	
☐	112D20		Silver	

Anodized Color Package

Forks, Rear Axle Clamp, and backrest bracket. N/A with Frog Legs.

☐	112EB07		Black	
☐	112EB08		Blue	
☐	112EB13		Gold	
☐	112EB09		Lime Green	
☐	112EB10		Purple	
☐	112EB11		Red	
☐	112EB12		Ti Grey	

Frame Graphics

Graphic designs on front frame tubes. N/A with transit & Matte Black Cherry.

☐	1128A24		Bloom	
☐	1128A25		Dragon	
☐	1128A26		Ghost Flames	
☐	1128A27		Retro Squares	
☐	1128A28		Tribal	
☐	1128A33		Floating Dots Black	
☐	1128A34		Floating Dots Matte Black	
☐	1128A78		Floating Dots Blue	
☐	1128A76		Floating Dots Orange	
☐	1128A35		Floating Dots Silver (grey)	
☐	1128A77		Floating Dots Red	

10. COLOR (CONT)

Built For Me

To inquire about any unique personalized for Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011



Frame graphics



11. APPENDIX

Handrims	Wheels					
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain	Proton
Aluminum Anodized / Painted	20,22,24	20,22,24,26	22,24,25	22,24,25,26	24	24,25
Plastic Coat	20,22,24	22,24	-	24,26	-	-
Natural Fit - Standard Grip	24	22,24,25,26	-	22,24,25,26	-	-
Natural Fit - Super Grip	24	22,24,25,26	-	22,24,25,26	-	-
Q Grip Neoprene Coated	-	24,25	-	24,25	-	-
Ellipse 3R	24	24,25	-	24,25	-	-
Carbolife - Gekko	-	24,25	-	24,25	-	24,25
Carbolife - Curve L Tetra grip	-	24,25	-	24,25	-	24,25
Omit Handrim	20,22,24	-	-	22,24,25,26	24	24,25

Tire	Wheel				
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Proton
Pneumatic	20,22,24	20,22,24,26		22,24,26	24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26		22,24,26	
Pneum w/ Airless Insert - LTW	20,22,24	20,22,24		20,22,24	
Full Poly	20,22,24	20,22,24		22,24	
SHOX Performance Poly	20,22,24	20,22,24,25		22,24,25	
Primo Orion Poly	20,22,24	20,22,24		20,22,24	
High-Pressure Clinch Blk	22,24	22,24,26	24	22,24,26	24
Schwalbe Marathon Plus	24	24,25,26	22,24,25	24,25,26	24,25
Schwalbe One	24	24,25	24,25	24,25	24,25
Schwalbe Down Town	24	24	24	24	24
Schwalbe Right Run	24	24	24	24	24
Knobby V Trak by Kenda	24	24	24	24	24
Mountain		24		24	

20"= 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590




Figure 6




Figure 7

Lower Leg Length



Figure 3

