



Zippie X'Cape
July 2025
EIFZ1

Date: _____	SUBMITTING FOR: ☐ Quote ☐ Order
Dealer Acct #: _____	
Dealer: _____	PO#: _____
Dealer Contact: _____	Ship To: _____
Dealer Address: _____	Attention: _____
Dealer City: _____ ST: _____ ZIP: _____	Address: _____
Dealer Phone: () _____ Fax: () _____	Address: _____
Confirm Via: ☐ Fax ☐ Email	Ship To City: _____ ST: _____ ZIP: _____
Mark For: _____	Ship To Phone: () _____ Fax: () _____

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	✓	Option #	HCPCS	Definition
---	----------	-------	------------	---	----------	-------	------------

1. FRAME

A. Zippie X'CAPE - TRANSIT INCLUDED 165 lbs. Weight Capacity	B. Zippie X'CAPE - TRANSIT OPTIONAL 165 lbs. Weight Capacity
A. FRAME STYLE (Includes True Fit) (N/A w/ Xtender)	B. FRAME STYLE (Includes True Fit)
☐ 121WM87 E1238 X'CAPE - Transit Incl. without seating (JAY seating opts available)	☐ 121WM83 E1238 X'CAPE - Transit Opt. without seating (JAY seating opts available)
☐ 121WM84 E1236 X'CAPE - Transit Incl. w/ Seating (Sling and Back Uph.)	☐ 121WM85 E1236 X'CAPE - Transit Opt. w/ Seating (Sling and Back Uph.)
Note: HCPCS Codes shown as submitted. Some sizes fit other HCPC codes i.e. 15"-18" widths = K0005	Note: HCPCS Codes shown as submitted. Some sizes fit other HCPC codes i.e. 15"-18" widths = K0005

TRANSIT SYSTEM WC-19

☐ 121TO48 K0108 Transit Opt. (N/A w/ Xtender)	☐ 121TO40 Omit Transit
--	---

FRAME TYPE

☐ 121FT64 Zippie-Swing Away Frame	☐ 121FT65 Zippie-80° Fixed Front Frame
--	---

2. SIZE/SEATING

2A. FRAME WIDTH

8", 9" and 10" widths N/A with Xtender

☐ 121W48 8" (Frame Width Growth 9"-10")	☐ 121W53 13" (Frame Width Growth 14"-15")
☐ 121W49 9" (Frame Width Growth 10"-11")	☐ 121W54 14" (Frame Width Growth 15"-16")
☐ 121W50 10" (Frame Width Growth 11"-12")	☐ 121W55 15" (Frame Width Growth 16"-17")
☐ 121W51 11" (Frame Width Growth 12"-13")	☐ 121W56 16" (Frame Width Growth 17"-18")
☐ 121W52 12" (Frame Width Growth 13"-14")	

2B/C. SLING DEPTH, FRAME DEPTH AND BACK POST POSITION. See matrix on page 2.

8"-11" Frame Depth with 20", 22" or 24" wheels and Personalized Seating must have hublock or Fixed Front frame. When selecting Frame Depth with JAY seating, add 1" to accommodate seating. Example: need 12" depth, order 13" depth.

2B. Sling Depth when choosing with Seating Sling must be up to 3" less than or equal to the Frame Depth (C.)

☐ 121GJ1 8"	☐ 121GJ7 14"
☐ 121GJ2 9"	☐ 121GJ8 15"
☐ 121GJ3 10"	☐ 121GJ9 16"
☐ 121GJ4 11"	☐ 121GJ10 17"
☐ 121GJ5 12"	☐ 121GJ11 18"
☐ 121GJ6 13"	

2C. Frame Depth and Back Post Position with Growth Note: Must be equal or greater than Sling Depth if chosen with seat Sling

☐ 121FD40 8" (Frame depth growth 9"-11", N/A w/ 24" wheels)	☐ 121FD46 14" (Frame depth growth 15"-17")
☐ 121FD41 9" (Frame depth growth 10"-12", N/A w/ 24" wheels)	☐ 121FD47 15" (Frame depth growth 16"-18")
☐ 121FD42 10" (Frame depth growth 11"-13")	☐ 121FD48 16" (Frame depth growth 17"-19")
☐ 121FD43 11" (Frame depth growth 12"-14")	☐ 121FD49 17" (Frame depth growth 18"-20")
☐ 121FD44 12" (Frame depth growth 13"-15")	☐ 121FD50 18" (Frame depth growth 19"-20")
☐ 121FD45 13" (Frame depth growth 14"-16")	

SEATING

Seat Sling

☐	121OP01	Standard Seat Sling	
☐	121OP04	Omit Seat Sling (available when base is ordered with seating)	

JAY Cushions (Optional)

To order any other JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out w/ desired width and depth

☐	121JC20	JAY Bundled Cushion	
		Cushion Name or Part#	
		Width (in.):	Depth (in.):

Seat Sling to Frame Depth Matrix

		Difference from front of frame to front of sling										
	Frame Length	Sling Length										
		8"	9"	10"	11"	12"	13"	14"	15"	16"	17"	18"
	8"	0"										
	9"	0"	0"									
	10"	0"	0"	0"								
	11"	1"	1"	1"	0"							
	12"		2"	2"	1"	0"						
	13"			3"	2"	1"	0"					
	14"				3"	2"	1"	0"				
	15"					3"	2"	1"	0"			
	16"						2"	1"	0"	0"		
	17"							1"	0"	0"	0"	
	18"								0"	0"	0"	0"

3.CASTERS/SEAT HEIGHT

Caster Wheel

☐	121CW30	3" Micro Light Up	
☐	121CW3	4" Low-Profile Polyurethane	
☐	121CW31	4" Micro Light Up	
☐	121CW64	4" x 1-1/2" Aluminum Soft Roll	
☐	121CW5	5" Low-Profile Polyurethane	
☐	121CW51	5" Micro Light Up	
☐	121CW6	E2219 5" Soft Roll (Grey)	
☐	121CW7	E2219 5"x1-1/2" Semi-Pneumatic	
☐	121CW66	5"x1-1/2" Aluminum Soft Roll	
☐	121CW34	6" Poly	
☐	121CW35	E2219 6"x1-1/2" Semi-Pneumatic	

Caster Forks

☐	121CF65	Multiple Position	
--------------------------	---------	-------------------	--

Front Seat Heights Matrix

Rear Wheel Size	Front Caster			
	☐ 3"	☐ 4"	☐ 5"	☐ 6"
☐ 12" ☐ 16" ☐ 20" ☐ 22"	☐ 13.0"	☐ 13.0"	☐ 14.0"	☐ 15.0"
	☐ 13.5"	☐ 13.5"	☐ 14.5"	☐ 15.5"
	☐ 14.0"	☐ 14.0"	☐ 15.0"	☐ 16.0"
	☐ 14.5"	☐ 14.5"	☐ 15.5"	☐ 16.5"
	☐ 15.0"	☐ 15.0"	☐ 16.0"	☐ 17.0"
	☐ 15.5"	☐ 15.5"	☐ 16.5"	☐ 17.5"
	☐ 16.0"	☐ 16.0"	☐ 17.0"	☐ 18.0"
	☐ 16.5"	☐ 16.5"	☐ 17.5"	
	☐ 17.0"	☐ 17.0"		
	N/A	☐ 17.5"		
☐ 24"	N/A	N/A	☐ 14.0"	☐ 15.0"
			☐ 14.5"	☐ 15.5"
	☐ 14.0"	☐ 14.0"	☐ 15.0"	☐ 16.0"
	☐ 14.5"	☐ 14.5"	☐ 15.5"	☐ 16.5"
	☐ 15.0"	☐ 15.0"	☐ 16.0"	☐ 17.0"
	☐ 15.5"	☐ 15.5"	☐ 16.5"	☐ 17.5"
	☐ 16.5"	☐ 16.5"	☐ 17.0"	☐ 18.0"
	☐ 17.0"	☐ 17.0"	☐ 17.5"	
	N/A	☐ 17.5"		
			N/A	N/A

	Standard Frame
	Hemi Fame (N/A w/ 90° hanger, 80° hanger footplate in lowest position may cause interference w/caster)
	Depends on Caster setting

FRONT SEAT HEIGHTS ARE SHOWN IN THE BOXES IN THE FRONT SEAT HEIGHTS MATRIX. DEPENDING ON CG SETTING, SEAT DEPTH, AND FRAME SIZE SOME OF THE SIZES MAY NOT BE AVAILABLE

4.REAR WHEELS/SEAT HEIGHT

XTENDER POWER ASSIST OPTION

22" & 24" wheels only. RSH 14" to 16" for 22", RSH 15.0" to 17" for 24". 0° camber, pneumatic tire. Anti-Tips no charge. Dual drive mode available through B4Me.

Xtender Power Assist

☐	121XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	
☐	121XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	
☐	121XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	
☐	121XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	

Please select one of the following handrim options

☐	121XH5		Plastic Coated	
☐	121XH6		Stainless Steel	

Xtender - Other Options

☐	121XA10		Spare NiMH Power Pack - NICKEL	
☐	121XA11		Spare Li-ion Power Pack - LITHIUM	

Xtender Wheel Options

☐	121XA8		SPARE 22" Lite Spoke Pneu Alum	
☐	121XA9	E2211 E2213	SPARE 22" Lite Spoke Pneu Plastic	
☐	121XA2		SPARE 24" Lite Spoke Pneu Alum	
☐	121XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic	

Rear Wheels

N/A w/ Xtender & 5 Spoke Mag Wheels. N/A w/ Hub Lock. Std Mag Wheels only available with One Arm Drive. 24" wheels. N/A with 8" and 9" frame depths. 8"-12" Frame Depth with 20-24" Rear Wheels and ELR's is ONLY available with Hub Lock Wheels.

☐	121RW48		12" Mag (Not available with handrims or one arm drive)	
☐	121RW14		16" Mag (Not available with handrims or one arm drive)	
☐	121RW3		20" Mag Wheel	
☐	121RW8		22" Mag Wheel	
☐	121RW1		24" Mag Wheel	
☐	121RW103		20" Mag Wheel 5 Spoke	
☐	121RW104		22" Mag Wheel 5 Spoke	
☐	121RW102		24" Mag Wheel 5 Spoke	
☐	121RW74		20" Lite Spoke Wheel	
☐	121RW75		22" Lite Spoke Wheel	
☐	121RW76		24" Lite Spoke Wheel	
☐	121RW06		22" Ultra light wt. Spoke Wheel 2.0	
☐	121RW10		24" Ultra light wt. Spoke Wheel 2.0	
☐	121RW28		22" Spinergy SPOX Black Spokes	
☐	121RW95		24" Spinergy SPOX Black Spokes	
☐	121RW32		22" Spinergy LX Black Spokes	
☐	121RW98		24" Spinergy LX Black Spokes	

Rear Wheel Tire

☐	121RT3		Full Poly Tire	
☐	121RT6	E2211 E2212	High Pressure Clincher	
☐	121RT27	E2211 E2212	Knobby V-TRAK by Kenda	
☐	121RT38	E2218	SHOX Performance Poly	
☐	121RT1		Pneumatic	
☐	121RT2	E2211 E2213	Pneumatic Airless Insert	
☐	121RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	
☐	121RT28	E2211 E2212	Schwalbe Speed	
☐	121RT4		Low Poly	

Tire	Wheel	5 Spoke Mag	Std. Mag	Lite Spoke	Ultra Light-weight	Spinergy
Full Poly Tire		20", 22", 24"	20", 22", 24"	20", 22", 24"		
High Pressure Clincher				22", 24"	24"	22", 24"
SHOX Performance Poly		20", 22", 24"		20", 22", 24"		22", 24"
Knobby V-TRAK				24"		24"
Pneumatic		20", 22", 24"	12", 16", 20", 22", 24"	20", 22", 24"		22", 24"
Pneumatic Airless Insert		20", 22", 24"	12", 16", 20", 22", 24"	20", 22", 24"		22", 24"
Schwalbe Marathon Plus				24"	22", 24"	22", 24"
Schwalbe Speed Run				24"	24"	24"
Low Poly			12"	24"		

22" = 489 mm 24" = 540 mm

Handrim

N/A w/ Xtender. Natural Fit N/A w/ 20" and 22". Carbolife N/A w/ 5 Spoke Mag and Mag Wheel. Carbolife and Ellipse NA w/ 20" and 22".

☐	121HR1		Aluminum Anodized	
☐	121HR4		Plastic Coated (N/A on Ultra and 22" Spinergy LX, SPOX)	
☐	121HR22		Natural Fit - Standard Grip (N/A on Ultra Lite or 20" Lite Spoke)	
☐	121HR23		Natural Fit - Super Grip (N/A on Ultra Lite or 20" Lite Spoke)	
☐	121HR39		Ellipse 3R	
☐	121HR40		Carbolife-Gekko (N/A Mag Wheel Spoke)	
☐	121 HR41		Carbolife-Curve L-Tetra Grip	
☐	121HR7		Omit Handrim (Option may not be available on all rear wheels)	

Axle Plate

☐	121AL33		Standard	
--------------------------	---------	--	----------	--

4. REAR WHEELS/SEAT HEIGHT Cont.

One Arm Drive

NA w/ Xtender. NA w/ 12"-16" wheels. Only avail. On 12"-16" widths and 20", 22" and 24" mag wheels. NA with Natural Fit handrims. NA with SPHA arms

☐	121AL16	E0958	One Arm Drive	
☐	121OA2		Standard Handrim	
☐	121HM1		☐ Right Hand Mount	
☐	121HM2		☐ Left Hand Mount	
☐	121OA1		Small Handrim (N/A w/ 20" rear wheels)	
☐	121HM1		☐ Right Hand Mount	
☐	121HM2		☐ Left Hand Mount	

Axle

☐	121AX1		Quick Release (N/A w/ One Arm Drive)	
☐	121AX11		Quad Quick Release (Only Available with 12" and 16" Mag Wheels)	

Center of Gravity Preset (dependent on wheels, casters & depth, From standard frame depth +2" forward and - depth growth up to -3")

☐	121CC28		-0.5"	
☐	121CC26		0 (Flush w/ backrest)	
☐	121CC16		0.5" (Only avail 9"-18" deep)	
☐	121CC20		1" (Only avail 10"-18" deep)	
☐	121CC15		1.5" (Only avail 10"-18" deep)	
☐	121CC17		2" (Only avail 16"-18" deep)	
☐	121CC18		2.5" (Only avail 16"-18" deep)	

Camber

☐	121CI1		0°	
☐	121CI2		3° (N/A w/ Single Post Ht. Adj. Arm)	

Rear Seat Height (see Allowable Squeeze chart)

☐	121RH32		11"	
☐	121RH33		11.5"	
☐	121RH9		12"	
☐	121RH10		12.5"	
☐	121RH11		13"	
☐	121RH12		13.5"	
☐	121RH13		14"	
☐	121RH14		14.5"	
☐	121RH15		15"	
☐	121RH16		15.5"	
☐	121RH17		16"	
☐	121RH18		16.5"	
☐	121RH19		17"	
☐	121RH26		17.5"	
☐	121RH20		18"	

Rear Wheel	Rear Seat Height															
	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	
12" low poly																
12" pneu																
16"																
20"																
22"																
24"																
24" pneu																

When tires are not listed default to the wheel size. All Heights are subject to vary slightly

Allowable Squeeze (may vary +/- 1/2")

Frame Depth	Max Squeeze
8"-9"	1.0"
10"-12"	1.5"
13"-15"	2.0"
16"-18"	2.5"

5. WHEEL LOCK

Hub Lock

Only available with 12" or 16" mag wheels or Lite Spoke wheels.
Not available with One Arm Drive.

☐	121LK01		12" Hub Lock	
☐	121LK02		16" Hub Lock	
☐	121LK04		20"-24" Hub Lock Lite Spoke	

Hub Lock Mount Position

Foot lock may limit folding width.

☐	121HH01		Foot lock Mount	
☐	121HH02		Seat Mount	
☐	121HH03		Back Mount (back ht. may be limited)	

Wheel Lock

High Mount will be mounted in low position on 12 and 16 rear wheels.
Some may not be available depending on wheel size and position.

☐	121WL1		High Mount Push	
☐	121WL2		High Mount Pull	
☐	121WL30		Compact	
☐	121WL10		Omit	

Extension Handle (Only avail w/ High Mount Wheel Locks)

☐	121EH1	E0961	6" Pair	
--------------------------	--------	-------	---------	--

Swing Away 8"-12" Frame Depth Wheel Lock go/no-go matrix

Frame Depth	C.G.	20" wheel	22" wheel	24" wheel	Frame Depth	C.G.	20" wheel	22" wheel	24" wheel
8"	-0.5"				11	-0.5"			
	0"					0"			
	0.5"					0.5"			
	1"					1"			
9"	-0.5"				12	1.5"			
	0"					-0.5"			
	0.5"					0"			
	1"					0.5"			
10"	-0.5"					1"			
	0"					1.5"			
	0.5"					No-go. User must use hublock.			
	1"					Go			
	1.5"								

6. HANGER AND FOOTREST

Z-finity Hangers/Footrest System

N/A on Fixed Frame.

Some interference may exist in certain configurations.

☐	121H162	Z-finity Hanger	
☐	121H160	Elevating Hanger - Ext Mount <i>Includes ELR pad. N/A w/ Z-finity. Some limitations may exist.</i>	
☐	121H144	Omit Hanger	

Z-finity Footrest Mount (Must Choose see Illustration)

☐	121DI22	Direct Hanger Mt/Fixed Frame - High	
☐	121DI23	Direct Hanger Mt/Fixed Frame - Low	
☐	121DI24	Extension Tube Mount	

Z-finity Footrest Pivot Arm (Must Choose, use matrix for leg length)

☐	121DI20	Short 2"	
☐	121DI21	Long 3-1/4"	

Footplate

Some interference may exist in certain configurations.

☐	121F320	K0040	Multi-Angle Adjustable - Pair (N/A w/ 8"-9" widths)	
☐	121F321	K0040	Angle Adjustable Platform Footboard	
☐	121F6		Omit Footplate	

Extension Tube Length (When Choosing Ext. Tube Mount)

☐	121E62		Hemi	
☐	121E1		Short	
☐	121E2		Medium	
☐	121E5		Long	

Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)

N/A with 8" and 9" width. Footplate mounted in lowest position may cause interference w/ caster.

☐	121H157		80° Swing In/Out Front Mount	
☐	121H158		90° Swing In/Out Front Mount	
Extensions (may interfere with some seat heights & casters)				
☐	121E23		2" (not mounted)	
☐	121E24		4" (not mounted)	


Footplate

☐	121F9	K0040	Kids Angle Adjustable (only available w/ 10"-16" frame widths)	
☐	121F8	K0040	Adult Angle Adj. (only available with 10"-16" widths)	

Swing In/Out Hanger Leg Lengths

(Seat sling to footrest measured)


80° Hanger	Min	Max
No Ext. Tube	5.7"	9.7"
2" Ext. Tube	11.2"	12.2"
4" Ext. Tube	11.2"	14.1"
90° Hanger	Min	Max
No Ext. Tube	6.0"	10.0"
2" Ext. Tube	11.5"	12.5"
4" Ext. Tube	11.5"	14.5"



Z-finity Matrix

Direct Hanger Mount

Extension Tube Mount



Direct Hanger Mounted			
Pivot Arm	Hanger Mounting Position	Min	Max
Short	High Mount	0.5"	4.5"
	Low Mount	3.5"	8.0"
Long	High Mount	-1.0"	5.5"
	Low Mount	2.5"	9.0"

Extension Tube Mounted			
Pivot Arm	Extension Tube Length	Min	Max
Short	Hemi Short	5.8"	14.1"
	Short	6.5"	15.1"
	Medium	8.4"	17.0"
	Long	10.3"	19.0"
Long	Hemi Short	4.6"	15.3"
	Short	5.2"	16.3"
	Medium	7.2"	18.3"
	Long	9.2"	20.2"


ELR Ext Mt.(MED Ext Tube)

Pivot Arm	Min	Max
Short	7.0"	18.0"
Long	6.0"	19.0"

Footrest Options

☐	121HL8	K0038	Leg Strap	
☐	121HL11	E0951	Heel Loops (N/A w/ Kids Angle Adj. Z-finity. N/A on widths below 14")	
☐	121IG2		Neoprene Impact Guards (N/A with footrest in highest mounted position)	

7. BACKRESTS									
BACKREST									
☐	121BT10	Standard Backrest (Angle adj./folding)				☐	121BT101	E2398	Dynamic Backrest (N/A w/ Fixed Back)
BACK TYPE (Choose Height Adjustable or Fixed Height)									
☐	121BT120	Height Adjustable (see matrix below)				OR	☐	121BT53	Fixed Height
Height Adjustable (See Backrest Height Range matrix) <i>Must choose Lower Base Tube, Backrest Style, and Upper Back Tube.</i>					Fixed Height				
Lower Base Tube (See Backrest Height Matrix Below)					<i>Fixed Height back N/A with upholstery</i>				
☐	121BP58	Short			☐	121BP9	Standard (16" Height)		
☐	121BP59	Medium			☐	121BP10	Medium (20" Height)		
☐	121BP60	Tall (N/A with Dynamic Backrest)			☐	121BP11	Tall (23" Height)		
Backrest Styles (See Images Below)					Fixed Height				
☐	121PH16	Straight Non Push Handle			☐	121BP61	Standard w/ Folding PH (16" Height)		
☐	121PH6	Integral Push handle			☐	121BP62	Medium w/ Folding PH (20" Height)		
☐	121PH17	Folding Push Handle with Ht. Adj. Back (8"-10" widths, stroller handle will face outward)			☐	121BP63	Tall w/ Folding PH (23" Height)		
☐	121PH15	Omit Back Tube (See omit matrix)							
Upper Back Tube (See Backrest Height Matrix Below)									
☐	121B245	Short							
☐	121B5	Medium							
Height Adjustable Backrest									
Preset Back height (Choose for Set Up) <i>No Upholstery w/ 19"-24"</i>									
☐	9"	☐	15"	☐	20"				
☐	11"	☐	16"	☐	21"				
☐	12"	☐	17"	☐	22"				
☐	13"	☐	18"	☐	23"				
☐	14"	☐	19"	☐	24"				
Backrest Height Ranges									
Upper Back Tube (Choose One)									
			Straight non-push handle Short	Straight non-push handle Medium	Integral push handle Short	Integral push handle Medium	Folding Push Handle Short	Folding Push Handle Medium	Omit
Lower Base Tube (Choose one)	Standard	Short	11"-15"	15"-19"	11"-15"	15"-19"	12"-16"	16"-20"	9"
		Medium	15"-19"	18"-23"	15"-19"	18"-23"	16"-20"	19"-24"	13"
		Tall	18"-22"	18"-24"	18"-22"	18"-24"	19"-23"	19"-24"	16"
	Dynamic	Short	14"-18"	18"-22"	14"-18"	18"-22"	15"-19"	19"-23"	12"
		Medium	18"-22"	21"-24"	18"-22"	21"-24"	19"-23"	22"-24"	16"



7. BACKRESTS Cont.

Back Upholstery

Uph N/A with Fixed Height back. Std Uph N/A w/Flip-Up, Adj. Locking Flip-Up and Length Adj Locking Flip-up. 3DX N/A on 8"-11" widths. Ballistic N/A on 8"-9" widths.

☐	121BU1		Standard Upholstery	
☐	121BU71	E2611	Backrest Tension Adj. 3DX Vented	
☐	121BU70	E2611	Backrest Tension Adj. Ballistic	
☐	121BU66		Omit Back Upholstery	

JAY Back (Optional)

To order any other JAY back on-chair, please See the JAY Price Lst. Requires Part #

☐	121JB02	JAY Bundled Backrest	
		Backrest Name or Part #	
		Width (in.):	Height (in.):

Stabilizer Bar

N/A w/ Dynamic back and Fixed Back. Only available w/ Integral push handle.

☐	121SZ1	Stabilizer Bar	
--------------------------	--------	----------------	--

8. ARMRESTS/SIDEGUARDS

Armrest and Armpad - Height Adjustable

Flip-Up, Adj. Locking Flip-Up and Length Adj. Locking Flip-up N/A with Std Back Upholstery. Low RSH in combination with large wheels may result in arm interference. SPHA arms N/A with One Arm Drive or Xtender.

☐	121AA58	E0973	T-Post Xtra Tall(with sideguard)	
☐	121AA59	E0973	T-Post 2X Tall(with sideguard)	
☐	121AA3		Flip-Up	
☐	121AA1	E0973	Adjustable Locking Flip-Up	
☐	121AA7	E0973	Length Adjust. Locking Flip-Up	
☐	121AA11	E0973	Single-Post Ht Adj. Std. (N/A w/ 3° camber or 8"-9" frame depths)	
☐	121AA29	E0973	Single-Post Height Adj. Low (N/A w/ 3° camber or 8"-9" frame depths)	
☐	121AA14		Omit Armrest	

Side Guard

N/A with height adjustable armrests and T-Posts

☐	121SG33	K0108	Plastic Std	
☐	121SG29	K0108	Plastic Low	


Armpad

N/A with T-Post and Adj. Locking Flip-Up and Foam only available with Length Adj. Locking Flip-Up.

☐	121AM23		Desk Classic	
☐	121AM24		Full Classic (N/A w/ 10"-13" Depth)	
☐	121AM9		Foam Short (N/A w/ T Post)	
☐	121AM10		Foam Long (N/A w/ T Post)	


Single Post Height Adjustable Armrest

Pad	Armrest Size	Transfer Tube Length (D)	Range (H)	Pad Length (L)	
				Desk	Full
Classic	STD	8"	9-1/2" to 14-1/2"	10-1/4"	14"
	LOW	5-1/4"	8" to 12"		



T-POST

Armrest Size	Range (H)	Pad Length (L)
XTALL	6" to 7"	7"
2X TALL	7-1/2" to 8-1/2"	



9. ACCESSORIES

Back Pack (one Backpack Black/Charcoal already included)

☐	121A1		Black Utility Backpack	
☐	121A2		Backpack Black/Charcoal Kids	

Miscellaneous

☐	1218A03		Caddy (Only avail w/ Ext. Mt. hangers)	
☐	1213A01		Tool Kit	
☐	121A2100	E2207	Crutch Holder	

Positioning Belt

☐	121A823	E0978	Belt 1" Dual Adj. Sqz Blk Pad	
☐	121A824	E0978	Belt 2" Aircraft Buckle	
☐	121A825	E0978	Belt 2" Aircraft Buckle Padded	
☐	121A826	E0978	Belt 1.5" Auto Buckle	
☐	121A814	E0978	1" Comfort Fit APB5116	

Upper Extremity Support Tray

N/A with 8"-9" widths

☐	121A29	E0950	X-Small (Suggested 10"-12" chair width)	
☐	121A30	E0950	Small (Suggested 13"-14" chair width)	
☐	121A31	E0950	Medium (Suggested 15"-17" chair width)	
☐	121A32	E0950	Large (Suggested 18"-20" chair width)	

Canopy (SunShade)

Need 2" back post clearance to mount canopy hardware

☐	1218A06	K0108	Canopy (Sunshade) 10"-14" Width	
☐	1218A07	K0108	Canopy (Sunshade) 15"-20" Width	

Seat Pouch

☐	1211A01		Black	
--------------------------	---------	--	-------	--

Anti-Tip

☐	121AT1	E0971	Rear Anti-Tip	
☐	121AO1	E0971	Front Anti-Tip (N/A with transit opt.)	

Spoke Guard

N/A with 12" and 16" rear wheels or 5 spoke mag or Xtender

☐	1214A01	K0065	Black	
☐	1214A02	K0065	Clear Zippie	

Designer Spoke Guard

Only available with 22" or 24" Spoke rear wheels

☐	1214A10	K0065	Basketball	
☐	1214A11	K0065	Chrome Wheels	
☐	1214A12	K0065	Graffiti	
☐	1214A13	K0065	Happy Lions	
☐	1214A14	K0065	Hot Wheels	
☐	1214A15	K0065	Mice n Cheese	
☐	1214A16	K0065	Princess	
☐	1214A17	K0065	Space Galactica	
☐	1214A18	K0065	Super Hero	
☐	1214A19	K0065	Zebra	

10. COLOR (Standard with frame and cross brace color choices)

Frame Color (Must choose)

Standard Colors

☐	190C214		Black (Gloss)	
☐	190C215		Black (Matte)	
☐	190C211		Celestial Sapphire	
☐	190C212		Desert Sage	
☐	190C213		Desert Tan	
☐	190C216		Metallic Shadow	
☐	190C217		Midnight Blue	
☐	190C218		Mystic Purple	
☐	190C219		Petal Pink	
☐	190C220		Pink Power	
☐	190C221		Red Pepper	
☐	190C222		Royal Blue	
☐	190C223		Salute Green	
☐	190C224		Scarlet Noir	
☐	190C225		Silver Ore	
☐	190C226		Solar Orange	
☐	190C227		Starlight White	
☐	190C228		Storm Gray	
☐	190C229		Turbo Green	
☐	190C230		Velocity Yellow	

Cross Brace Color

Standard Colors

☐	190C627		Celestial Sapphire - Cross Brace	
☐	190C628		Desert Sage - Cross Brace	
☐	190C629		Desert Tan - Cross Brace	
☐	190C630		Gloss Black - Cross Brace	
☐	190C631		Matte Black - Cross Brace	
☐	190C632		Metallic Shadow - Cross Brace	
☐	190C633		Midnight Blue - Cross Brace	
☐	190C634		Mystic Purple - Cross Brace	
☐	190C635		Petal Pink - Cross Brace	
☐	190C636		Pink Power - Cross Brace	
☐	190C637		Red Pepper - Cross Brace	
☐	190C638		Royal Blue - Cross Brace	
☐	190C639		Salute Green - Cross Brace	
☐	190C640		Scarlet Noir - Cross Brace	
☐	190C641		Silver Ore - Cross Brace	
☐	190C642		Solar Orange - Cross Brace	
☐	190C643		Starlight White - Cross Brace	
☐	190C644		Storm Gray - Cross Brace	
☐	190C645		Turbo Green - Cross Brace	
☐	190C646		Velocity Yellow - Cross Brace	

Personalized Colors (Must choose cross brace std. color)

☐	121C43		Mossy Oak Camouflage	
☐	121C42		Stars and Stripes	
☐	121C108		Zebra (Choose standare base color) N/A Black Matte/Gloss	

Additional Options Through Built-4-Me

To inquire about any unique personalization for Zippie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 - 9011

Commonly Requested Options Include:

- ✓ Solid Seat
- ✓ Grade Aid
- ✓ Drum Brake

