

Multiroller



**SPECIAL
PURPOSE
STROLLER**

SCRIPTING AND
ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Standard Package

Multiroller Stroller Package

Stroller frame with easy to fold function, with upholstered seat, height adjustable flip up footplate, interchangeable front end stroller frame wheelset, drum brakes and canopy.

2610-1000-000



Specifications

Push Handle Height (cm)	86 – 120
Overall Width (cm)	68
Wheel Diameter, Front/Buggy (cm)	18
Wheel Diameter, Rear (cm)	41
Folded size (with rear wheels on) (cm)	50 (H) x 60 (W) x 84 (L)
Folded size (with rear wheels removed) (cm)	34 (H) x 56 (W) x 84 (L)
Load Capacity Stroller (kg)	75
Weight of Stroller (kg)	9
Seat Width (cm)	30 – 38
Seat Depth (cm)	33
Lower Leg Length (cm)	20 – 30
Foot Support (cm)	38 (W) x 23 (L)
Back Support Width (cm)	30 – 40
Back Support Height (cm)	29 – 48

Accessory Options

Running Wheel Add On		5 Point Harness	
2610-8416-000		Standard Buckle 2610-6500-000	
		Pen Key Release 2610-6500-025	
Seat Depth Extension		Infill for Smaller Users	
2610-0281-000		2610-0393-000	
Foot Support Extension			
Extends footplate adjustment by 10cm (40cm maximum).			
2610-5370-000			
2 Point Hip Belt – Centre Release – Small		2 Point Hip Belt – Pen Key Release – Small	
1405-6211-028		1405-6211-025	
Anti-Tip		Storage Basket	
2610-8420-000		Maximum weight 3kg. 2881-8051-200	

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Multiroller_20250904