



	Folding Rear Wheel Drive (RWD)	Rear Wheel Drive (RWD) Non Expandable		Mid Wheel Drive (MWD) Non Expandable	MWD Paediatric Expandable	MWD Expandable						RWD Expandable		Front Wheel Drive (FWD) Expandable	Standing	4WD	
	Q50R Carbon	Q100R	Q200R	Q400M	Z500M Mini	Q500M Mini	Q500M	Q700M	Magic 360	Frontier V6	Centro Glide	Q700R	Frontier V4	Q700F	Q700 UP	Extreme X8	XT4
Contracts	N/A	Accessible	Accessible	Enable	Accessible & Enable	Accessible & Enable	Accessible & Enable	N/A	N/A	N/A	N/A	Enable	N/A	Enable	N/A	N/A	N/A
Seat Width	Fixed 17.5"	17"	17"	15" - 20"	10" - 16"	16 - 20"	15" - 20"	16 - 22"	14" - 22"	14" - 24"	14.5" - 21.5"	15" - 20"	14" - 24"	15" - 20"	16" - 20"	16" - 24"	12" - 24"
Seat Depth	Fixed 16.5"	16" - 20"	16" - 20"	16" - 22"	10" - 19"	16 - 22"	16" - 22"	14 - 20"	14" - 22"	14" - 24"	12" - 25"	16" - 22"	14" - 24"	16" - 22"	16" - 22"	16" - 24"	13" - 24"
Floor to Seat Pan Height Range	19" (approx.)	17" (432mm)	17" (435mm)	16" - 19"	17" - 19"	17" - 19"	16" - 19"	16.5" - 19"	17" - 19"	17" - 19"	18"	16.5" - 19"	17" - 19"	16.5" - 19"	17.7" (450mm)	18.5" - 20.5"	19" (480mm)
Power Seat Elevation (maximum height)	N/A	N/A	N/A	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)	12" (300mm)
Power Recline (maximum angle)	N/A	N/A	N/A	N/A	N/A	130° Anti-Shear	130° or 150° Anti-Shear	130° or 150° Anti-Shear	160°	160°	138°	130° or 150° Anti-Shear	160°	130° or 150° Anti-Shear	90° - 175° Anti-Shear	150°	135°
Speed Range	6kph 12km range per battery (optional 2nd Battery)	6kph	10kph	10kph 40km range	10kph 30km range	10kph 30km range	10kph 40km range	10kph 40km Range	10kph 35km Range	10kph 35km Range	10kph	10kph 40km Range	10kph 35km Range	10kph 40km Range	10kph 40km range	10kph	10 kph (5 kph when lifted)
Overall Width	570mm / 22.5"	540mm / 21"	580mm / 22.8"	620mm / 24"	520mm / 20.5"	520mm / 20.5"	620mm / 24"	622mm / 24.5"	610mm / 24" (with Urban Wheels) 635mm / 25" (with Cross Over Wheels) 660mm / 26" (with Off Road Wheels)	655mm / 25.8" (Compact option w/Urban Wheels) 660mm / 26" (Compact option w/Off Road Wheels) 695mm / 27" (with Hybrid Wheels) 710mm / 28" (with All Terrain Wheels)	620mm / 24"	622mm / 24.5"	640mm / 27" (with Hybrid Wheels) 710mm / 28" (with All Terrain Wheels)	622mm / 24.5"	630mm / 24.8"	700mm / 27.5"	700mm / 27.5"
Overall Length	923mm / 36.3" (including footrest as non-removable)	1010mm / 39.5"	1070mm / 42.1"	800mm / 31.5"	800mm / 31.5"	800mm / 31.5"	800mm / 31.5"	840mm / 33"	980mm / 38.5"	980mm / 38.5" (Compact option) 1050mm / 41.5"	902mm / 35.5" (with no foot support)	940mm / 37"	940mm / 37"	940mm / 37"	860mm / 33.8"	1030mm / 40.5"	1000mm / 39"
Chair Weight	14.5kg (without batteries) (Single battery = 1.4kg)	95kg (209.4lbs)	98kg (216lbs)	165kg	145kg	145kg	165kg	175kg	149kg	144kg	155kg	175kg	145kg	175kg	140/175kg	165kg	185kg/220kg
Battery Size	24V, 10Ah Lithium-Ion (optional 2nd battery)	55Ah	55Ah	60Ah	50Ah	50Ah	60Ah/72Ah	72Ah	73Ah	73Ah	12V 70 - 90Ah	72Ah	73Ah	72Ah	60Ah/72Ah	70Ah or 90Ah	115Ah
Turning Radius	800mm / 31.5"	1600mm / 63'	1700mm / 67'	530mm / 21"	480mm / 19"	480mm / 19"	530mm / 21"	560mm / 22"	555mm / 21.5"	580mm / 23"	450mm / 18"	735mm / 29"	520mm / 20.5"	735mm / 29"	560mm / 22"	1420mm / 56"	1410mm / 55.5"
Drive Wheel Size	8"	14"	14"	13"	13"	13"	13"	14"	14" (3 sizes)	14" (2 sizes)	14"	14"	14" (2 sizes)	14"	14"	14"	14"
Arm Support Options	• Flip-back armrests	• Lift-out Removable	• Lift-out Removable	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Single Post Height Adjustable • Dual Post Flip Back Height Adjustable • Cantilever Flip Back Arm	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Standard Single Stern Armrests • Flexi Armrests with Dual Locking Post	• Standard Single Stern Armrests • Flexi Armrests with Dual Locking Post	• Adjustable Height and Flip Back	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Standard Single Stern Armrests • Flexi Armrests with Dual Locking Post	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Flip Back Height Adjustable Reclining Armrests	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm	• Cantilever Flip Back • Single Post flip-back or removable • Reclining flip-back with 2pt post locking arm
Lower Leg Support Options	• Carbon-fibre flip-up footplate, non-removable	• Swing Away • Angle Adjustable Footplate	• Swing Away • Angle Adjustable Footplate	• Fixed or Power Centremount - Single or dual flip up • Swing Away 70° • Swing Away Manual Elevating • Stump Support	• Centremount - Fixed or Angle Adjustable • Swing Away - 70, 80 or 90° • Swing Away Manual Elevating	• Fixed or Power Centremount - Single or dual flip up • Swing Away; Power Elevating or Manual Elevating with Composite Angle Adjustable Footplates • Swing Away 70° - Aluminium Angle Adjustable flip-up footplate • Stump Support	• Fixed or Power Centremount - Single or dual flip up • Swing Away; Power Elevating or Manual Elevating with Composite Angle Adjustable Footplates • Swing Away 70° - Aluminium Angle Adjustable flip-up footplate • Stump Support	• Fixed or Power Centremount - Single or dual flip up • Swing Away; Power Elevating or Manual Elevating with Composite Angle Adjustable Footplates • Swing Away 70° - Aluminium Angle Adjustable flip-up footplate • Stump Support	• Fixed and Power Elevating Centre Mount Legrests - Single Flip Up Footplate • Swing Away Manual or Power Elevating Legrests - Aluminium Angle Adjustable Footplate • Stump Support	• Fixed and Power Elevating Centre Mount Legrests - Single Flip Up Footplate • Swing Away Manual or Power Elevating Legrests - Aluminium Angle Adjustable Footplate • Stump Support	• Swing Away 70° Leg Hangers • Fixed and Power Elevating Centre Mount Legrests • Swing Away Individual Manual or Power Elevating legrests • Stump Support • Angle Adjustable Footplates, also Extra Long available	• Fixed or Power Centremount - Single or dual flip up • Swing Away; Power Elevating or Manual Elevating with Composite Angle Adjustable Footplates • Swing Away 70° - Aluminium Angle Adjustable flip-up footplate • Stump Support	• Fixed and Power Elevating Centre Mount Legrests - Single Flip Up Footplate • Swing Away Manual or Power Elevating Legrests - Aluminium Angle Adjustable Footplate	• Fixed or Power Centremount - Single or dual flip up • Swing Away; Power Elevating or Manual Elevating with Composite Angle Adjustable Footplates • Swing Away 70° - Aluminium Angle Adjustable flip-up footplate • Stump Support	• Power Elevating Centre Mount Legrest • Split • Height adjustable • Angle adjustable	• Fixed and Power Elevating Centre Mount Legrests - Single Flip Up Footplate • Swing Away Manual or Power Elevating Legrests - Aluminium Angle Adjustable Footplate	• Fixed and Power Elevating Centre Mount Legrests - Single Flip Up Footplate • Swing Away Manual or Power Elevating Legrests - Aluminium Angle Adjustable Footplate
Max User Weight	136kg	125kg	136kg	136kg	136kg	136kg	136kg	136kg	160kg	155kg	175kg (AS 3695 Full compliance - Class B Indoor/Outdoor)	136kg	155kg	160kg	120kg	155kg - 182kg	160kg with power lift/tilt 182kg with fixed frame

*Lift seat with back canes - Easy for therapist to fit any back support. Need an option to upgrade to pro seat for higher weight capacity and shear recline.