

Piper Stroller



**FULLY FOLDABLE
STROLLER**

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

September 2025

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size Options & Specifications

	Piper Standard Package
Code	2510-0003-034
Approx. Age of User	4-8 years
Max. Load (kg)	50
Seat Width (cm)	23/27/34*
Seat Depth (cm)	24/29*
Seat to Floor Height (cm)	45
Back Support Height (cm)	68
Footplate Height (cm)	10-30
Approx. Stroller Weight (kg)	14.5
Folded Dimensions (cm)	100 (L) x 31 (W) x 30 (H)
Backrest Angle (°)	90-135

*full range requires use of Seat Inlay accessory

Standard Package Includes:

- Removable upholstery
- Comfortable 5-point support harness
- Back support recline adjustment
- Height-adjustable and detachable 1 piece foot support
- Solid tyres
- Removable wheels
- Swivel locks on front wheels
- Basket under stroller for storage



Accessories

Seat Inlay

2510-6080-000



Upper Limb Support

2510-6140-000



Footplate Height Reducing Spacer

(Reduces footplate height by 5 or 10cm).
2510-6170-000



Sun Canopy

2510-7519-000



Rain Cover – must be used with Sun Canopy

2510-7522-000

Pen Key Release Hip Belt

(30cm Pad).
1405-6210-025



Pen Key Release 5 Point Harness

2880-6511-025

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Piper_20250904_AU