

Q50R/carbon

Travel lighter



RIDICULOUSLY LIGHT: 14.5 KG

One of the world's lightest carbon fiber folding power wheelchairs.

MODERN DESIGN

Sleek, ergonomic design with premium materials – the epitome of comfort and style. Available in 5 accent colors: Blue, Orange, Purple, Red, and Yellow.

GO UP TO 24 KM/H

Bring an (optional) second battery along for the ride to double your range to up to 24 km/h!

EFFORTLESS FOLDING

Flip the clip & push the backrest down. You don't even need to remove the battery.

COMPACT FOLDED SIZE

Folds to a compact size of 77 cm length x 57 cm width x 37 cm height, to fit small car trunks, buses, and closets.

SMOOTH RIDE

Its carbon fiber frame can naturally flex over uneven surfaces, providing unique shock-absorbing characteristics to deliver one of the smoothest, most efficient rides around.

LIGHT BUT INCREDIBLY STRONG

Can support a user weight of up to 136 kg.

NARROW 57 CM WIDTH

Its slim width easily fits through narrow doorways and navigates busy rooms.

RELIABLE MOBILITY


250W high-efficiency brushless motors provide the power and reliability you need.



TECHNICAL DATA

DIMENSIONS	93 cm Long x 57 cm Wide
SEAT DIMENSIONS	43 cm Wide X 42 cm Deep
BACK HEIGHT	46.5 cm
MAX. RANGE	12 km/h (one battery), 24 km/h (with optional additional battery)
MAX. SPEED	3.7 MPH
BATTERY SIZE	10Ah lithium-ion
DRIVE WHEEL POSITION	Rear
WHEEL SIZE	20 cm Drive Wheel; 16.5 cm Castors
USER WEIGHT CAPACITY	136 kg
PRODUCT WEIGHT	Chair (without batteries): from 14.5 kg, Battery (single): 1.4 kg

Refer to order form as specifications may vary based on configuration.



Learn more



This document can be viewed in PDF
format at www.SunriseMedical.com.



MK-130131 Rev A ©02.2024



800.333.4000

www.SunriseMedical.com

