



Quickie 2 Lite

January 2022

EIQ2X

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

ST: _____

ZIP: _____

Dealer Phone: () _____

Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: () _____

Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base model

E1236 coding still applies for EIQ2X with sizes 14" or less.

☐	002WM58	K0005	Quickie 2 Lite QLT	\$2,343.00
--------------------------	---------	-------	--------------------	------------

Transit Option

☐	002TO49	K0108	Transit System WC19	\$345.00
--------------------------	---------	-------	---------------------	----------

Frame Style (265 lbs. Weight Capacity)

See figure 1. Measure from front frame to floor

☐	002FS31	Swing Away Frame	NC
☐	002FS32	Fixed Front 70"	NC
☐	002FS33	Fixed Front 80"	NC

Frame Inset (Fixed Front Only)

See **FIGURE 2** Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2.

☐	002FF1	0"	NC
☐	002FF2	1" N/A for 14" and 15" Wide using composite angle adj footrest	NC

2. SIZE / SEATING

Frame Width

See **FIGURE 2** Measured from outside of seat tube to the outside of opposite seat tube.

☐	002W5	14"	NC
☐	002W6	15"	NC
☐	002W7	16"	NC
☐	002W8	17"	NC
☐	002W9	18"	NC
☐	002W10	19"	NC
☐	002W11	E2201 20"	\$537.00

Seat Depth

See **FIGURE 1** and **Seat Squeeze Matrix** Measured from back canes to front edge of upholstery.

☐	002SD7	14" Front of seat sling 2" behind front edge of frame.	NC
☐	002SD8	15"	NC
☐	002SD9	16"	NC
☐	002SD10	17"	NC
☐	002SD11	18"	NC
☐	002SD72	19"	\$81.00
☐	002SD13	E2203 20"	\$590.00

Frame Depth

If growth is desired Seat Depth (*) may be ordered up to 2" smaller on swing away frames and up to 3" on fixed front frames. Please refer to matrix in page 8 of Appendix regarding frame Depth Vs seat Depth.

☐	002FD5	14"	NC
☐	002FD6	15"	NC
☐	002FD7	16" (*)	NC
☐	002FD8	17" (*)	NC
☐	002FD9	18" (*)	NC
☐	002FD10	19" (*)	NC
☐	002FD36	20" (*)	NC

Seat Sling

☐	002SL1	Standard Seat Sling	NC
☐	002SL6	Omit Seat Sling	-\$15.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.



Cushion name or part number must be filled out with desired width and depth

☐	002JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____	Depth: _____

Cushion

☐	002CU1	K0669 2" Cushion	\$161.00
--------------------------	--------	------------------	----------

2. SIZE / SEATING (CONT)



Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
14"-15"	2.0"
16"-19"	2.5"
20"	3.0"

3. CASTERS / SEAT HEIGHT

Front Seat Height

20" Max height for Fixed front. 20.5"-21" Only available in Swing Away

☐ 002SH42	14.5"	NC	☐ 002SH18	18.0"	NC
☐ 002SH12	15.0"	NC	☐ 002SH20	18.5"	NC
☐ 002SH13	15.5"	NC	☐ 002SH21	19.0"	NC
☐ 002SH14	16.0"	NC	☐ 002SH23	19.5"	NC
☐ 002SH15	16.5"	NC	☐ 002SH24	20.0"	NC
☐ 002SH16	17.0"	NC	☐ 002SH26	20.5"	NC
☐ 002SH17	17.5"	NC	☐ 002SH27	21.0"	NC

Caster Fork

☐ 002CF65		Multiple Position	NC
☐ 002CF69		Multiple Position Single Sided Fork	\$182.00
☐ 002CF53	E1015	Frog Legs Suspension (0-100 lbs)	\$467.00
☐ 002CF54	E1015	Frog Legs Suspension (101-190 lbs)	\$467.00
☐ 002CF55	E1015	Frog Legs Suspension (191-265 lbs)	\$467.00


Caster Wheels

☐ 002CW1		3" Micro	NC
☐ 002CW3		4" Low Profile Polyurethane	NC
☐ 002CW31		4" Micro Lighted	\$75.00
☐ 002CW72		4" x 1" Aluminum Polyurethane	\$118.00
☐ 002CW44	E2219	4" x 1.25" Semi Pneumatic	\$59.00
☐ 002CW64		4" x 1.5" Aluminum Soft Roll	\$199.00
☐ 002CW5		5" Low Profile Polyurethane	NC
☐ 002CW51		5" Micro Lighted	\$92.00
☐ 002CW6	E2219	5" Soft Roll (Grey) N/A w/ Frog Legs	\$97.00
☐ 002CW65		5" x 1" Aluminum Polyurethane	\$118.00
☐ 002CW7	E2219	5" x 1.5" Semi-Pneumatic	\$92.00
☐ 002CW66		5" x 1.5" Aluminum Soft Roll	\$209.00
☐ 002CW34		6" Polyurethane	NC
☐ 002CW33	E2214 E2215	6" Pneumatic	\$81.00
☐ 002CW35	E2219	6" x 1.5" Semi Pneumatic	\$92.00
☐ 002CW68		6" x 1.5" Aluminum Soft Roll	\$220.00

Caster Housing

☐ 002CP15		Integral Caster Housing	NC
☐ 002CP12		Hemi Integral Cstr Housing Subtract 2" from Standard Front STFH Matrix	NC

Multi Position Single side Frog Leg



Front Seat Height Matrix

Integral Caster Housing (subtract 2" from FSH for Hemi Housing)											
Caster Size	*FSH - SA Front Frame	*FSH - Fixed Front Frame	Fork Hole (1 = Bottom Hole)								
			Multi Pos Low			Multi Pos Std			Frog Legs		
			Spacer Size:			Spacer Size:			Spacer Size:		
			0"	0.5"	1"	0"	0.5"	1"	0"	0.5"	1"
3"	16.5	15.5	2								
	17.0	16.0	1	2							
	17.5	16.5		1	2	2					
	18.0	17.0			1	1	2		2		
	18.5	17.5					1	2	1	2	
	19.0	18.0						1		1	2
4"	19.5	18.5									1
	20.0	19.0									
	17.0	16.0	2								
	17.5	16.5	1	2		3					
	18.0	17.0		1	2	2	3		3		
	18.5	17.5			1	1	2	3	2	3	
5"	19.0	18.0					1	2	1	2	3
	19.5	18.5						1			
	20.0	19.0								1	2
	20.5	19.5									1
	19.0	18.0				2					
	19.5	18.5				1	2				
6"	20.0	19.0					1	2	1		
	20.5	19.5						1		1	
	21.0	20.0									1

*FSH values given are estimates. Actual seat to floor height may vary based on seat depth, squeeze & caster selection.

Rear Seat Height						
13" & 13.5" N/A w/ transit						
☐	002RH11	13.0"	NC	☐	002RH18	16.5"
☐	002RH12	13.5"	NC	☐	002RH19	17.0"
☐	002RH13	14.0"	NC	☐	002RH26	17.5"
☐	002RH14	14.5"	NC	☐	002RH20	18.0"
☐	002RH15	15.0"	NC	☐	002RH27	18.5"
☐	002RH16	15.5"	NC	☐	002RH21	19.0"
☐	002RH17	16.0"	NC	☐	002RH37	19.5"

Axle			
☐	002AX1	Quick Release Steel	NC
☐	002AX11	Quad Quick Release Mountain Wheel only	\$135.00

Xtender Power Assist (265 weight capacity)			
☐	002XT6	E0986 22" Xtender Power Assist JWX-2 NICKEL	\$8,637.00
☐	002XT7	E0986 22" Xtender Power Assist JWX-2 LITHIUM	\$9,496.00
☐	002XT4	E0986 24" Xtender Power Assist JWX-2 NICKEL	\$8,637.00
☐	002XT5	E0986 24" Xtender Power Assist JWX-2 LITHIUM	\$9,496.00

Please select one of the following handrim options			
☐	002XH5	Plastic Coated	NC
☐	002XH6	Stainless Steel	NC

Xtender - Other Options			
☐	002XA10	Spare NiMH Power Pack - NICKEL	\$1,498.00
☐	002XA11	Spare Li-ion Power Pack - LITHIUM	\$1,927.00

Xtender Wheel Options			
☐	002XA8	SPARE 22" Lite Spoke Pneu Alum	NC
☐	002XA9	E2211 E2213 SPARE 22" Lite Spoke Pneu Plastic	\$215.00
☐	002XA2	SPARE 24" Lite Spoke Pneu Alum	NC
☐	002XA4	E2211 E2213 SPARE 24" Lite Spoke Pneu Plastic	\$215.00

Rear Wheel Type			
See availability matrix. Not applicable w/ Xtender Power Assist Option.			
Ultra Lightweight Wheels only available w/ Alum Anodized Handrim			
☐	002RW103	20" Mag Wheel 5 Spoke	NC
☐	002RW104	22" Mag Wheel 5 Spoke	NC
☐	002RW102	24" Mag Wheel 5 Spoke	NC
☐	002RW90	20" Lite Spoke Wheel	\$215.00
☐	002RW91	22" Lite Spoke Wheel	\$215.00
☐	002RW92	24" Lite Spoke Wheel	NC
☐	002RW93	25" Lite Spoke Wheel	NC
☐	002RW94	26" Lite Spoke Wheel	\$215.00
☐	002RW77	22" Ultra Lightweight Spoke Wheel	\$268.00
☐	002RW78	24" Ultra Lightweight Spoke Wheel	\$268.00
☐	002RW28	22" Spinergy SPOX Black	\$907.00
☐	002RW95	24" Spinergy SPOX Black	\$907.00
☐	002RW96	25" Spinergy SPOX Black	\$907.00
☐	002RW97	26" Spinergy SPOX Black	\$907.00
☐	002RW32	22" Spinergy LX Black	\$1,138.00
☐	002RW98	24" Spinergy LX Black	\$1,138.00
☐	002RW99	25" Spinergy LX Black	\$1,138.00
☐	002RW40	26" Spinergy LX Black	\$1,138.00
☐	002RW50	24" Mountain	\$268.00

Rear Wheel Type (Cont)			
* Omit Rear wheel will default chair configuration set up to a 24" wheel.			
☐	002RW25	Omit Rear Wheel and Axle *	-\$55.00
☐	002RW49	Omit Rear Wheel (Use w/24" Xtender)	NC
☐	002RW111	Omit Rear Wheel (Use w/22" Xtender)	NC

NEW! Rear Wheel Tire			
☐	002RT1	Pneumatic	NC
☐	002RT2	E2211 E2213 Pneumatic Airless Insert N/A with Xtender.	\$92.00
☐	002RT24	E2211 E2213 Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$105.00
☐	002RT3	Full Poly	NC
☐	002RT38	E2218 SHOX Performance Poly	\$102.00
☐	002RT25	E2218 Primo Orion Poly	\$116.00
☐	002RT39	E2211 E2212 High-Pressure Clincher Black	\$102.00
☐	002RT27	E2211 E2212 Knobby V-TRAK	\$102.00
☐	002RT21	E2211 E2212 Schwalbe Marathon Plus Evolution	\$151.00
☐	002RT12	E2211 E2212 Mountain Tire Only available on 24" Mountain Rear Wheel	\$102.00

Handrim			
Natural Fits available in 22", 24", and 25". Ultra lightweight wheels ONLY available w/ alum ano handrim. Natural Fit w/ 5 spoke MAG only available in 24". Not all handrims available on all wheels			
☐	002HR1	Aluminum Anodized	NC
☐	002HR4	Plastic Coated -N/A on every rim	\$118.00
☐	002HR22	Natural Fit - Standard Grip	\$478.00
☐	002HR25	Natural Fit - No Thumb Grip	\$414.00
☐	002HR23	Natural Fit - Super Grip	\$532.00
☐	002HR30	Natural Fit LT- Standard Grip	\$478.00
☐	002HR31	Natural Fit LT- Super Grip	\$532.00
☐	002HR32	Natural Fit LT - No Thumb Grip	\$414.00
☐	002HR26	The Surge - Oval w/Gription Strip	\$429.00
☐	002HR27	The Surge LT - Oval w/Gription Strip	\$462.00
☐	002HR28	Q-Grip Neoprene Coated	\$376.00
☐	002HR24	Spinergy Flex Rim Available on 24" & 25" Spinergy LX wheel	\$751.00
☐	002HR7	Omit Handrim Option may not be available on all wheels	NC
☐	002HR8	8 Oblique N/A w/ Xtender, Only available with 24" or 26" or Lite Spoke rims.	\$349.00

Rear Seat to Floor Height Matrix																
Note: These values are representative of a 16" seat depth with 1" of seat squeeze and standard axle plate. Due to variances caused by seat angle, rear seat height ranges may shift up or down by as much as 1/2". Call customer service or online configurator for more information.																
Wheel/Tire	13.0	13.5	14.0	14.5	15.0	15.5	16.0	16.5	17.0	17.5	18.0	18.5	19.0	19.5		
20" All	X	X	X	X	X	X	X	X	N/A	N/A	N/A	N/A	N/A	N/A		
22" Poly	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A	N/A	N/A		
22" Pneu	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A	N/A		
24" Pneu	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A		
24" All but Pneu	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A		
25" SHOX	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A	N/A		
25" All but SHO	N/A	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X	N/A		
26" All	N/A	N/A	N/A	N/A	N/A	N/A	X	X	X	X	X	X	X	X		

4. REAR WHEELS / SEAT HEIGHT (cont)

Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Tire					
Pneumatic	20,22,24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		24,26
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Pressure Clincher black			22,24,26	22,24	22,24,26
Knobby V-TRAK			24		24
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Perform. Poly	20,22,24		24,25		24,25
20"= 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590					

Camber			
☐	002C11	0°	NC
☐	00217CI	2° (default)	NC
☐	002CI3	4°	NC

Center of Gravity Preset			
☐	002CC2	1.0" (default)	NC
☐	002CC3	1.5"	NC
☐	002CC4	2.0"	NC
☐	002CC5	2.5"	NC
☐	002CC6	3.0"	NC
☐	002CC7	3.5"	NC

5. WHEEL LOCKS

Wheel Lock			
No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth			
☐	002WL46	High Mount Push	NC
☐	002WL50	High Mount Pull	NC
☐	002WL75	Ergo Scissor	\$118.00
☐	002WL64	Compact	\$151.00
☐	002WL24	Unilateral Wheel Lock Right	\$322.00
☐	002WL23	Unilateral Wheel Lock Left	\$322.00
☐	002WL10	Omit	NC

Wheel Lock Option			
☐	002WL12	Do Not Mount	NC

Extension Handle			
N/A w/ ergo scissor or compact wheel locks			
☐	002EH7	E0961 Extension Handle (Pair)	\$81.00

Grade Aid			
Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options			
☐	002GA1	E0974 Grade Aid	\$161.00

6. HANGER & FOOTREST

Swing - Away Hangers			
Footplate mounts with Extension Tube using clamping hardware			
☐	002H145	60° Swing In/Out Hanger	NC
☐	002H139	70° Swing In/Out Hanger	NC
☐	002H141	80° Swing In/Out Hanger	NC
☐	002H32	Omit Hanger and Footplate	-\$15.00

Footrest Extension Tubes			
See Extension matrix for leg length range			
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC

Extension Matrix

Frame Type	Footrest Angle	Inset	Footrest Type and Extension							
			Composite/ Aluminum			Platform		Tubular/ Autofold		
			Short	Medium	Long	Short	Long	Short	Long	
			Min	Max						
Fixed Front	70°	0°	13.5	17.5	13.5	16.5	15.5	12	12.5	
			13.5	19.5	15	16.5	20.5	16.5	18	
		1°	13.5	17.5	13.5	16.5	17	12	13.5	
			13.5	19.5	16	16.5	20.5	16.5	18	
	80°	0°	14.5	17.5	14.5	18	16.5	13	13.5	
			14.5	19.5	15.5	20.5	18	18	19.5	
		1°	14.5	17.5	14.5	18	16.5	13	14.5	
			14.5	19.5	17	20.5	18	18	19.5	
		60°	12	13.5	12	13	12.5			
			12	15.5	12	17				
		70°	12	16.5	12	14				
			12	18.5	12	18.5				
	80°	Min	12	16.5	12	13.5				
		Max	15.5	17.5	14	17.5				

Z-finity Hangers/Footrest System			
(Maximum weight capacity for Z-Finity footrests is 165 lbs)			
Some interference may exist in certain configurations.			
☐	002H162	Z-finity Hanger	\$376.00

Z-finity Footrest Mount			
☐	002DI22	Direct Hanger Mt/Fixed Frame - High	NC
☐	002DI23	Direct Hanger Mt/Fixed Frame - Low	NC
☐	002DI24	Extension Tube Mount	NC

Z-finity Footrest Pivot Arm			
☐	002DI20	Short 2"	NC
☐	002DI21	Long 3 1/4"	NC

Z-finity Footplate			
Some interference may exist in certain configurations.			
☐	002F320	K0040 Multi Angle Adjustable (Pair) (not avail. w/8"-9" widths)	\$199.00
☐	002F321	K0040 Angle Adjustable Platform Footboard	\$230.00
☐	002F6	Omit Footplate	NC

Z-finity ExtensionTube Length (When Choosing Ext. Tube Mount)			
☐	002E62	Hemi Short	NC
☐	002E1	Short	NC
☐	002E2	Medium	NC
☐	002E5	Long	NC

Option #	HCPCS	Definition	Price	Option #	HCPCS	Definition	Price																																																																								
6. HANGER & FOOTREST (cont)																																																																															
Footplates				Lower Leg Length Setting																																																																											
Heel loop or leg strap not included, See Footplate Options Section																																																																															
☐	002F1		Composite (14"-20"w)	NC	☐	002LO16	16.5"	NC																																																																							
☐	002F271	K0040	Composite Angle Adj-(16"-20"w)	\$194.00	☐	002LO17	17.0"	NC																																																																							
☐	002F3	K0040	Aluminum Angle Adj (14"-20"w)	\$204.00	☐	002LO18	17.5"	NC																																																																							
☐	002F48	K0040	Platform Flip Up	\$204.00	☐	002LO19	18.0"	NC																																																																							
☐	002F245	K0040	Tubular Flip Up Fixed front frame only	\$107.00	☐	002LO20	18.5"	NC																																																																							
☐	002F250	K0040	Auto Fold Angle Adj Fixed front frame only	\$215.00	☐	002LO21	19.0"	NC																																																																							
Footplate Options				☐	002LO13	15.0"	NC	☐	002LO22	19.5"	NC																																																																				
☐	002HL1	E0951	Heel Loop	\$43.00	☐	002LO14	15.5"	NC	☐	002LO23	20.0"	NC																																																																			
☐	002HL5	K0038	Leg Strap	\$43.00	☐	002LO15	16.0"	NC	☐	002LO24	20.5"	NC																																																																			
7. BACKREST / SEATING																																																																															
Back Type				JAY Back (Optional)																																																																											
☐	002BT10		Standard (9"-20")	NC	☐	002JB01	J3 Back (See J3 Order Form)	Price List																																																																							
☐	002BT84		Angle Adjustable Back (3" & 4")	\$199.00																																																																											
Back Post Options				Military Service Patches (Optional)																																																																											
Standard Back type options:				☐	002VU2		Army Patch	\$107.00																																																																							
☐	002PH6		Integral Push Handles	NC	☐	002VU4	Air Force Patch	\$107.00																																																																							
☐	002PH11		Integral Push Handles w/8°Medium & Tall only, N/A w/ Angle Adj Back	\$38.00	☐	002VU12	Coast Guard Patch	\$107.00																																																																							
☐	002PH2		Omit Push Handles	NC	☐	002VU1	Marines Patch	\$107.00																																																																							
☐	002PH12		Omit Push Handles w/8° Medium & Tall only, N/A w/ Angle Adj Back	NC	☐	002VU3	Navy Patch	\$107.00																																																																							
☐	002PH17		Folding Push Handle with Ht. Adj. Back N/A w/Angle Adj Back	\$242.00	☐	002VU5	Purple Heart Patch	\$107.00																																																																							
Angle Adj Back type options:				☐	002VU9		Purple Heart with Army Patch	\$166.00																																																																							
☐	002PH5		Fold-Down Push Handles N/A w/ standard Upholstery and J3 Back	\$107.00	☐	002VU11	Purple Heart with Air Force Patch	\$166.00																																																																							
				☐	002VU13		Purple Heart with Coast Guard Patch	\$166.00																																																																							
				☐	002VU8		Purple Heart with Marines Patch	\$166.00																																																																							
				☐	002VU10		Purple Heart with Navy Patch	\$166.00																																																																							
Back Height (Std or Angle Adjustable)				Angle Adjustable Back Setting																																																																											
See Matrix				☐	002AH1		3" Pivot Point Setting	NC																																																																							
☐	002B1		Low	NC	☐	002AH2		4" Pivot Point Setting	NC																																																																						
☐	002B5		Medium	NC	Select desired setting																																																																										
☐	002B9		Tall	NC	☐	-14°	☐	-12°	☐	-9°	☐	-7°	☐	-1°	☐	0°	☐	+4°	☐	+5°																																																											
Back Height Setting (Std or Angle Adjustable)				Backrest Height Matrix																																																																											
17" max height for shipping. 21" only available w/ Ang Adj Back																																																																															
☐	002UX4	9"	NC	☐	002UX18	16"	NC	<table><tr><th>Back Type</th><th>Push Handles</th><th>Backrest Height</th><th>Avail. Ht. Range</th></tr><tr><td rowspan="6">Standard - (Fixed)</td><td rowspan="3">Integral</td><td>Low</td><td>9"-12"</td></tr><tr><td>Medium</td><td>13"-16"</td></tr><tr><td>Tall</td><td>17"-20"</td></tr><tr><td rowspan="3">Fold Down</td><td>Medium</td><td>12"-15"</td></tr><tr><td>Tall</td><td>16"-19"</td></tr><tr><td rowspan="3">Omit</td><td>Low</td><td>9"-12"</td></tr><tr><td>Medium</td><td>13"-16"</td></tr><tr><td>Tall</td><td>17"-20"</td></tr><tr><td rowspan="6">Angle Adjustable (3" Pivot)</td><td rowspan="3">Integral</td><td>Low</td><td>13"-16"</td></tr><tr><td>Medium</td><td>17"-20"</td></tr><tr><td>Tall</td><td>N/A</td></tr><tr><td rowspan="3">Fold</td><td>Medium</td><td>16"-19"</td></tr><tr><td>Low</td><td>13"-16"</td></tr><tr><td>Medium</td><td>17"-20"</td></tr><tr><td rowspan="3">Omit</td><td>Tall</td><td>N/A</td></tr><tr><td rowspan="6">Angle Adjustable (4" Pivot)</td><td rowspan="3">Integral</td><td>Low</td><td>14"-17"</td></tr><tr><td>Medium</td><td>18"-21"</td></tr><tr><td>Tall</td><td>N/A</td></tr><tr><td rowspan="3">Fold</td><td>Medium</td><td>17"-20"</td></tr><tr><td>Low</td><td>14"-17"</td></tr><tr><td>Medium</td><td>18"-21"</td></tr><tr><td rowspan="3">Omit</td><td>Tall</td><td>N/A</td></tr><tr><td rowspan="3">Folding Push Handle - (Stroller Handle w/ Knuckle)</td><td rowspan="3">N/A</td><td>Low</td><td>9"-12"</td></tr><tr><td>Medium</td><td>13"-16"</td></tr><tr><td>Tall</td><td>17"-20"</td></tr></table>				Back Type	Push Handles	Backrest Height	Avail. Ht. Range	Standard - (Fixed)	Integral	Low	9"-12"	Medium	13"-16"	Tall	17"-20"	Fold Down	Medium	12"-15"	Tall	16"-19"	Omit	Low	9"-12"	Medium	13"-16"	Tall	17"-20"	Angle Adjustable (3" Pivot)	Integral	Low	13"-16"	Medium	17"-20"	Tall	N/A	Fold	Medium	16"-19"	Low	13"-16"	Medium	17"-20"	Omit	Tall	N/A	Angle Adjustable (4" Pivot)	Integral	Low	14"-17"	Medium	18"-21"	Tall	N/A	Fold	Medium	17"-20"	Low	14"-17"	Medium	18"-21"	Omit	Tall	N/A	Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	Low	9"-12"	Medium	13"-16"	Tall	17"-20"
Back Type	Push Handles	Backrest Height	Avail. Ht. Range																																																																												
Standard - (Fixed)	Integral	Low	9"-12"																																																																												
		Medium	13"-16"																																																																												
		Tall	17"-20"																																																																												
	Fold Down	Medium	12"-15"																																																																												
		Tall	16"-19"																																																																												
		Omit	Low	9"-12"																																																																											
Medium	13"-16"																																																																														
Tall	17"-20"																																																																														
Angle Adjustable (3" Pivot)	Integral	Low	13"-16"																																																																												
		Medium	17"-20"																																																																												
		Tall	N/A																																																																												
	Fold	Medium	16"-19"																																																																												
		Low	13"-16"																																																																												
		Medium	17"-20"																																																																												
Omit	Tall	N/A																																																																													
	Angle Adjustable (4" Pivot)	Integral	Low	14"-17"																																																																											
			Medium	18"-21"																																																																											
Tall			N/A																																																																												
Fold		Medium	17"-20"																																																																												
		Low	14"-17"																																																																												
		Medium	18"-21"																																																																												
Omit	Tall	N/A																																																																													
	Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	Low	9"-12"																																																																											
			Medium	13"-16"																																																																											
Tall			17"-20"																																																																												
☐	002UX6	10"	NC	☐	002UX20	17"	NC																																																																								
☐	002UX8	11"	NC	☐	002UX22	18"	NC																																																																								
☐	002UX10	12"	NC	☐	002UX24	19"	NC																																																																								
☐	002UX12	13"	NC	☐	002UX26	20"	NC																																																																								
☐	002UX14	14"	NC	☐	002UX28	21" (*)	NC																																																																								
☐	002UX16	15"	NC																																																																												
Backrest Upholstery																																																																															
Std EXO upholstery and Tension adjustable binding is always black.																																																																															
☐	002BU1		Standard Upholstery	NC																																																																											
☐	002BU95		Standard Upholstery EXO	NC																																																																											
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue				\$15.00																																																																											
☐	002BU70	E2611	Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$338.00																																																																											
☐	002BU89	E2611	Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$355.00																																																																											
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue				\$15.00																																																																											
☐	002BU71	E2611	Backrest Tension Adjust 3DX Vented *Reference 186BU71 for billing	\$355.00																																																																											
☐	002BU82	E2611	Jay GO Backrest	\$437.00																																																																											
☐	002BU14		Omit Back Upholstery	-\$15.00																																																																											

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
8. ARMREST / SIDE GUARDS									
Padded Swing Away					Omit Armrest				
Some limitations when used w/ angle adj back in the rear most setting									
☐	002AA18	E0973	Padded Swing-Away STD (9" to 13")	\$242.00	☐	002AA14		Omit Armrest	NC
		Setting	☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"						
Height Adjustable					Armpad for Height Adjustable				
4 Point (FS) J3 mounting bracket will cause interference with Flip Back Height Adj armrest.					Ergo armpad only available with Single post Height adj & Flip back armrest				
☐	002AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$306.00	☐	002AM23		Desk Length - Classic (10")	NC
☐	002AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$306.00	☐	002AM24		Full Length - Classic (14")	NC
☐	002AA65	E0973	Flip Back Height Adj LITE (8.5"-13.5")	\$338.00	☐	002AM27		Ergo Pad (9.5"-13.5")	NC
					Side Guard				
☐	002SG33	K0108	Plastic Standard (8.5" High)	\$156.00	☐	002SG29	K0108	Plastic Low (6.5" High)	\$156.00

9. KEY OPTIONS & ACCESSORIES									
Anti-Tip					Additional Cover & Pack Accessories				
☐	002AT1	E0971	Rear Anti-Tip (NC w/ Xtender)	\$161.00	☐	0028A02		Gel Cover Hanger/Frame	\$118.00
Backpack					☐	0028A01		Gel Cover Swing Away Arm	\$161.00
☐	002A1		Quickie Backpack Black	\$81.00	Positioning Belts				
☐	002A2		Quickie Kids Backpack Black	\$64.00	☐	002A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$75.00
Seat Pouch					☐	002A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$118.00
☐	0021A01		Seat Pouch Black /Charcoal	\$54.00	☐	002A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$128.00
JAY Positioning					☐	002A120		Positioning Belt Slide Clamps	\$54.00
☐	002K1	E1020	JAY Lower Extremity S/A support	\$424.00	Spoke Guards				
			☐ Left ☐ Right Position		<i>Only available w/ Mountain, and Lite Spoke wheels</i>				
☐	002K2		Omit Hanger	-\$10.00	<i>Not available with 5 spoke MAG wheels</i>				
			☐ Left ☐ Right		☐	0024A01	K0065	Quickie Black Spoke Guards	\$123.00
☐	002K7		Pad Size	NC	☐	0024A02	K0065	Quickie Clear Spoke Guards	\$123.00
			☐ 6" x 8" ☐ 8" x 10"		Miscellaneous				
					☐	0028A03		Caddy	\$54.00
					☐	002IG2		Impact Guard Neoprene	\$97.00
					☐	002LS14	K0038	Extra Leg Strap Black	\$38.00
					☐	0028A21		Tool Kit	\$59.00
					☐	002A999		Quickie Water Bottle & Holder	\$59.00

10. COLOR									
NEW! Standard Frame Colors					Matte Frame Colors				
☐	002C140		Aztec Gold	NC	☐	002C116		Matte Black	NC
☐	002C2		Black (Not avail as base w/Zebra)	NC	☐	002C151		Matte Black Cherry	NC
☐	002C141		Black Cherry	NC	☐	002C150		Matte Electric Blue	NC
☐	002C10		Black Opal	NC	☐	002C149		Matte Evergreen	NC
☐	002C38		Blue Opal	NC	☐	002C152		Matte Purple	NC
☐	002C36		Candy Blue	NC	Personalized Frame Color				
☐	002C14		Candy Purple	NC	☐	002C112		Desert Camouflage	\$424.00
☐	002C35		Candy Red	NC	☐	002C43		Mossy Oak Camouflage	\$424.00
☐	002C142		Electric Blue	NC	☐	002C42		Stars & Stripes	\$424.00
☐	002C143		Evergreen	NC	☐	002C108		Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Black	\$424.00
☐	002C44		Glossy White	NC	Anodized Color Package				
☐	002C126		Glow	NC	<i>Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.</i>				
☐	002C144		Green Apple	NC	☐	002EB07		Black - Caster only	\$54.00
☐	002C147		Hot Sparkle Pink	NC	☐	002EB08		Blue	\$107.00
☐	002C148		Mauve Pink	NC	☐	002EB13		Gold	\$107.00
☐	002C103		Rootbeer	NC	☐	002EB09		Lime Green	\$107.00
☐	002C114		Sparkle Silver	NC	☐	002EB10		Purple	\$107.00
☐	002C145		Sunrise Orange	NC	☐	002EB11		Red	\$107.00
☐	002C113		Titanium Color Paint	NC	☐	002EB12		Ti Grey	\$107.00
☐	002C8		Yellow	NC					



Sunrise Medical (US) LLC
2842 Business Park Ave. • Fresno, CA 93727 • USA