

Quickie Iris



ULTIMATE TILT-IN-SPACE MANUAL WHEELCHAIR

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details

Assessor Name		Client Name	
Phone		Funding	
Email		Trial Approval #	
Delivery Address			

Package & Size Options

Medium Package (Depth 16"-20")	Long Package (Depth 18"-22")
Wheelchair Width 14" 1567-0014-000	Wheelchair Width 14" 1567-0114-000
Wheelchair Width 16" 1567-0016-000	Wheelchair Width 16" 1567-0116-000
Wheelchair Width 18" 1567-0018-000	Wheelchair Width 18" 1567-0118-000
Wheelchair Width 20" 1567-0020-000	Wheelchair Width 20" 1567-0120-000

Medium Package Includes:	Long Package Includes:
• True-fit growth system • Tie downs • Width-adjusting system • Medium seat frame • 6 x 1.5" semi-pneumatic front castors • 16" rear mag wheels, pneumatic • 80° legrest hangers with angle adjustable footplates • Single post height-adjustable arm supports with full length pads • Height-adjustable push handle • Growable seat pan • Anti-tips • Black opal frame colour • Foot operated tilt • Foot operated brakes • Medium seat depth range: 16-20" • Spex 4-point belt	• True-fit growth system • Tie downs • Width-adjusting system • Long seat frame • 6" x 1.5" semi-pneumatic front castors • 16" rear mag wheels, pneumatic • 80° legrest hangers with angle-adjustable footplates • Single post height-adjustable arm supports with full length pads • Height-adjustable push handle • Growable seat pan • Anti-tips • Black opal frame colour • Foot operated tilt • Foot operated brakes • Long seat depth range: 18-22" • Spex 4-point belt



Specifications

Seat Width (R)		
Seat Depth (S)		
Back Height (T)		
Seat to Footplate		
Seat to Floor		

Castor Upgrades

Short Forks for Low Seat to Floor	8 x 2" Pneumatic	8 x 2" Pneumatic with Airless Insert	8 x 1.5" Semi Pneumatic
1568-8415-100	1568-8401-010	1568-8402-020	1568-8403-010


Rear Wheel Upgrades

16" Mag with Pneumatic Airless Type	20" 5 Spoke Mag	22" 5 Spoke Mag	24" 5 Spoke Mag
1568-01-9568-116	Pneumatic 1568-9520-110	Pneumatic 1568-9522-110	Pneumatic 1568-9524-110
	Airless Insert 1568-9520-120	Airless Insert 1568-9522-120	Airless Insert 1568-9424-120
24" Spoke Wheel with Pneumatic Tyre	24" Spoke Wheel with Pneumatic Airless Tyre	BMX Tyre Upgrade	
1568-01-9524-010	1568-01-9524-020	1460-PART-055	

Wheel Lock


Change Wheel Lock to Push to Lock	Push to Lock as additional to footlock	Push to Lock 6" Extensions	Drum Brakes 16" Pneumatic	Drum Brakes 24" Pneumatic
1568-8451-000	1568-8455-000	1579-161930	1568-8452-16P/16S	1568-8453-24P/24S

Footrest Hangers

Offset/Flared Hanger Receiver (Pair) Adds 2" width to legrests.		1568-5151-000					
Change to 90° Front Mount Hanger (Pair)		Change to 80° Front Mount Hanger (Pair)		Change to 70° Front Mount Hanger (Pair)		Change to 70° Extension Hanger (Pair)	
(When purchased with wheelchair). 1568-5090-000		(When purchased with wheelchair). 1568-5080-000		(When purchased with wheelchair). 1568-5070-000		(When purchased with wheelchair). 1568-5071-000	
Front Mount Footrest Extension Tube		Elevating Legrests		Articulating Legrests			
2" (Pair) 1568-5122-000		14-15" Wide 1568-5034-000		14-15" Wide 1568-5044-000			
4" (Pair) 1568-5124-000		16-17" Wide 1568-5036-000		16-17" Wide 1568-5046-000			
6" (Pair) 1568-5126-000		18-20" Wide 1568-5038-000		18-20" Wide 1568-5048-000			

Footplates

Heel Loops for Footplates		Footplate Padding		Spex Impact Guards	
12-14" Wide 1568-5312-14		10-11" Wide 1415-5310-300		70° Hangers 1415-5360-300	
15-16" Wide 1568-5315-16		12-13" Wide 1415-5312-300		80 & 90° Hangers 1415-5361-300	
17-18" Wide 1568-5317-18		14-15" Wide 1415-5314-300		Articulating Hangers 1415-5363-300	14-15" Wide
19-20" Wide 1568-5319-20		16-17" Wide 1415-5316-300			16-17" Wide
		18-22" Wide 1415-5318-300			18-20" Wide

Contracture Footrest			
Short (10.5" - 13" Leg Lengths)	Medium (12.5" - 15" Leg Lengths)	Long (14.5" - 17" Leg Lengths)	
14" 1568-5170-14S	14" 1568-5170-14M	14" 1568-5170-14L	
16" 1568-5170-16S	16" 1568-5170-16M	16" 1568-5170-16L	
18" 1568-5170-18S	18" 1568-5170-18M	18" 1568-5170-18L	
20" 1568-5170-20S	20" 1568-5170-20M	20" 1568-5170-20L	

Footplates (Continued)

Jay Foot Box - Standard Height		Spex Calf Strap	
		Available as Padded (1) or Neoprene (0) - replace the X when ordering	
14" 1453-5014-STD		14" 1402-6X72-017	
16" 1453-5016-STD		16" 1402-6X73-017	
18" 1453-5018-STD		17" 1402-6X74-017	
20" 1453-5020-STD		18" 1402-6X75-017	
		20" 1402-6X76-017	

Spex Calf Panels	
14-15" 1402-6314-017	
16-17" 1402-6316-017	
18-19" 1402-6318-017	
20" 1402-6320-017	
Custom 1402-6329-017	

Wheelchair Canes

Dynamic Canes		Reclining Back Canes	
Preferred Set-Up 1568-2170-000		Suits sizes 14" - 20" 1568-243356-000	
Soft			
Medium			
Hard			

Arm Supports

Length Adjustable Locking Flip-Up Arm Supports	Desk Length Arm Support Pad (Each)	Quickie Iris Dual Post Flip-Back Arm Support Upgrade
Full length pad (back cane mounted). 1568-5651-000	(When purchased with wheelchair). 1568-5532-000	1568-5622-099

Armrest (Continued)

Spex Comfi Padded Style Armrests		Spex Armrest Attachment Upgrade
Size 1 Narrow 1156-5811-410	Size 1 Wide 1156-5821-410	For mounting to existing wheelchair armrest mounts. 1156-5880-100
Size 2 Narrow 1156-5812-410	Size 2 Wide 1156-5822-410	
Size 3 Narrow 1156-5813-410	Size 3 Wide 1156-5823-410	

Tray

Spex Clear Tray – 5ml	Spex Black Tray – 5ml	Zippie Style Tray	Travel Tray
Small 1311-7111-000 (suits 10-13" wheelchairs)	Small 1311-7101-000	X-Small 1560-7110-100	Small 1311-7441-209
Medium 1311-7112-000 (suits 14-17" wheelchairs)	Medium 1311-7102-000	Small 1560-7113-100	Medium 1311-7442-209
Large 1311-7113-000 (suits 18-21" wheelchairs)	Large 1311-7103-000	Medium 1560-7115-100	Large 1311-7443-209
		Large 1560-7118-100	
		X-Large 1560-7121-100 	

Tray Accessories

Elbow Protraction Guards		Neoprene Elbow Pads, 8mm
Guards for Small Tray 1311-7351-000		Small 1311-7311-400
Guards for Medium Tray 1311-7352-000		Medium 1311-7312-400
Guards for Large Tray 1311-7353-000		Large 1311-7313-400

Tray Padding			
	8mm Neoprene Padding	10mm CCF Foam Padding	
Small Tray	1311-7431-291	1311-7431-360	
Medium Tray	1311-7432-291	1311-7432-360	
Large Tray	1311-7333-291	1311-7343-360	

Canopy Options

Iris Canopy	Folding Sunhood	UV Shade	Rain Cover
Suits 10" - 14" Chairs 1568-125320	450mm Wide 1446-7545-000	1446-7710-000	Bingo stroller style. 2880-7651-000
Suits 15" - 20" Chairs 1568-125321	550mm Wide 1446-7555-000		
			

Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other markings on the paper.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Quickie Iris_20250127