

SAFE2GO

by spex



VEHICLE HARNESS

spex


SCRIPTING AND ORDERING GUIDE

20251120


Step 1 – Select a Vest Package

Note: If additional vests are required, please select these in the ‘Optional Extras’ section at the end of this script form.

Chest circumference (cm):



Sitting shoulder height (cm):



A) Standard Size Vest

Please select size from the sizing table below:

Select Vest Size:	1	2	3	4	5	6	7
Crotch Strap Size	S	S	M	M	L	L	L
Chest circumference (cm)	65-69	70-74	75-79	80-85	86-94	95-104	105-115
Sitting shoulder height (cm)	36-46	40-50	40-50	49-59	53-63	56-66	58-68



- **Note:** If one of your measurements is outside the size range, order a user-matched size.
- If you require a user-matched size vest, please continue to step 1) B.

Package Includes: Vest, crotch strap, webbing slicer kit, and buckle release tool.

Product Code: 4610-100X-025
(Replace X with vest size number – e.g. 1 - 7)

B) User-Matched Size Vest

If you don't require a user-matched size vest, skip this step.



Chest circumference (cm):

Waist circumference (cm):

Sitting shoulder height (cm):

- This includes a user-matched crotch strap, sized based on the vest size.

Package Includes: Vest, crotch strap, webbing slicer kit, and buckle release tool.


Product Code: 4610-1099-025

Step 2) Select Vehicle Mount


There are two styles of vehicle mount available. Check which type is best suited to your vehicle. Select extra mounts and crotch straps to use the harness across multiple vehicles.

A – TOP TETHER VEHICLE MOUNT:

Check that your vehicle has a top tether anchor point and ISOFIX connection points.



Measure the distance between the top tether ISOFIX anchor point and the back face of the seat:



- If this measurement is more than 19cm, the Top Tether Vehicle Mount will fit in the vehicle.
- The maximum distance is 128cm.
- There is no user-matched option for the Top Tether Vehicle Mount.

Top Tether Vehicle Mount

Product Code:
4610-2010-000

Qty:



Extra Crotch Straps

Small
4610-6101-000

Qty:

Medium
4610-6102-000

Qty:

Large
4610-6103-000


Qty:



B – VERTICAL SEAT VEHICLE MOUNT



- Check that you can pass your hand through the gap in between the seat and the backrest.
- Some vehicles have a flap over the gap which can be disconnected.



- Measure the circumference of the backrest part of the seat.
- If this measurement is between 100cm and 170cm, then the mount will fit in the vehicle.

Vertical Seat Vehicle Mount

Product Code:
4610-2020-000

Qty:

Extra Crotch Straps

Small
4610-6101-000

Qty:

Medium
4610-6102-000


Qty:

Large
4610-6103-000

Qty:


User-Matched Size:

If your seat circumference doesn't fit inside these measurements and you would like to order a user-matched size, please provide the following measurements:



Circumference of the backrest (cm):

Product Code: 4610-2029-000



Measure from the gap to the back face of the seat, provide this as one measurement:

Height and Thickness of the backrest (cm):

Optional Extras

Vests Only

Standard,
Product Code:
4610-300X-025
(Replace X with vest size
numbers 1 - 7)

Size:

User-Matched,
Product Code:
4610-3099-025



Crotch Straps Only

Small
4610-6101-000

Qty:

Medium
4610-6102-000

Qty:

Large
4610-6103-000

Qty:



Buckle Cap

Product Code:
4610-9000-025



Pen Popper Release Tool with Retractor

Product Code:
1400-6010-025



Seat Belt Cutter Kit

Product Code:
4604-0010-017



User-Matched Crotch Strap

Please provide the following measurements:

Waist circumference (cm):

Measure around the waist, just above the top of the pelvis.
The person should sit/stand tall with arms a little away
from the body.

4610-6199-000

Production Options

**Replace Spex Safe2Go
Vehicle Harness Pen
Popper Buckles with
Centre Release Buckles**
1409-4600-028



spex

solutions@spexseating.com
spexseating.com

DISTRIBUTED BY:

Script: Safe2Go_20251120