

Sharky



SPECIALISED TILT-IN-SPACE COMMODE

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

January 2025

Client Details

Client/User Name:	
Date:	
DOB:	

Get it Faster

9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

Product Details

Options Available	Size 1	Size 2
35cm Back Support Height	6081-0010-100	N / A
45cm Back Support Height	6081-0010-000	6082-0015-100
55cm Back Support Height	N / A	6082-0015-000
Seat Width (cm)	24, 27, 30	29, 32, 35
Seat Depth (cm)	27 – 36	27 – 40
Maximum Traversable Toilet Height (cm)	52	
Tilt Adjustment Angle (°)	-10 to +45	
Back Support Recline Angle (°)	72 – 120	
Minimum Distance Between Lateral Supports (cm)	19	24
Leg Elevation Angle (°)	78 – 138	
Foot Support Angle (°)	80 – 100	
Arm Support Height (cm)	14.5 – 20.5	
Wheel Diameter (mm)	75	
Load Capacity (kg)	60	
Frame Length (cm)	70	
Frame Depth (cm)	56	



Sharky Size 1



Sharky Size 2

Upgrade Kit to Sharky Size 2

Upgrade Sharky mobility base from size 1 to size 2 using this kit.
6085-8081-086



Mobility Base (required)

(Sizes are specific for each size Sharky)

Sharky mobility base with height adjustable compact folding base, -10 to 45° tilt in space and goes over toilet heights from 42 - 53cm.	
For Sharky Size 1 6081-0015-000	For Sharky Size 2 6082-0010-000



Head Support

Suitable for use on either size Sharky

Sharky headrest, multi adjustable. 6085-4750-000	
---	--

Lateral Supports

Both sizes suitable for use on either size Sharky

Sharky Lateral Pads. 11 x 6cm 6085-3310-000	
Sharky Lateral Pads. 14 x 8cm 6085-3320-000	

Hip Support (required)

Options for Size 1 Sharky

24cm		27cm		30cm	
Sharky seat width reducing hip supports for seat width of 24cm. 6085-1110-000		Sharky seat width reducing hip supports for seat width of 27cm. 6085-1120-000		Sharky seat width reducing hip supports for seat width of 30cm. 6085-1130-000	

Options for Size 2 Sharky

29cm		32cm		35cm	
Sharky seat width reducing hip supports for seat width of 29cm. 6085-1210-000		Sharky seat width reducing hip supports for seat width of 32cm. 6085-1220-000		Sharky seat width reducing hip supports for seat width of 35cm. 6085-1230-000	

Seat Options

Seat Infill Pad		Armrests	
Sharky removable seat infill pad for hygiene opening 6085-0401-000		Sharky removable and height adjustable armrests 6085-5500-000	
Medial Thigh Support (Both sizes suitable for use on either size Sharky)			
Size 1 8cm wide x 12cm deep 6085-1500-000	Size 2 10cm wide x 14cm deep 6085-1520-000		
Splash Guard		Commode Pan	
Sharky removable splash guard 6085-7675-000		Sharky removable commode pan with retainer 6085-8210-000	

Belt Options (recommended)

Pelvi-Loc (for use on either size Sharky)		Standard Hip Belt (for use on either size Sharky)	
Pelvi-Loc 3D for maximum pelvic security		Sharky hip belt with neoprene pads	
X-Small 1430-6340-024		Size 1 Pad Width 12.5cm 6085-1910-000	
Small 1430-6341-024		Size 2 Pad Width 16.5cm 6085-1920-000	

Foot Positioning

Foot Support System			
Sharky removable foot support system (These options are specific for each size Sharky).	Size 1 For Size 1 Sharky 6085-1810-000	Size 2 For Size 2 Sharky 6085-1820-000	
Calf Strap		Ankle Huggers (for use on either size Sharky)	
For Size 1 Sharky 6085-6160-000		Size 1 Ankle circumference 20-37cm approx 6085-6101-000	
For Size 2 Sharky 6085-6260-000		Size 2 Ankle circumference 28-47cm 6085-6102-000	

Harnesses


Shoulder Harness (for use on either size Sharky)		
Sharky neoprene shoulder harness.		
Size 1 24.5cm high x 19cm wide 6085-6610-000	Size 2 31.5cm high x 21.5cm wide 6085-6620-000	

Handles

Sharky Push Handles	
6085-8600-000	

Accessories

Toilet Adapter		
Size 1 (for size 1 Sharky) 6085-0600-000	Size 2 (for size 2 Sharky) 6085-0620-000	

Bath Subframe		
Sharky bath subframe with seat angle adjustment. Minimum width of bath: 36cm		
Size 1 (for size 1 Sharky) 6085-0500-000	Size 2 (for size 2 Sharky) 6085-0520-000	

Tray Table			
Sharky depth adjustable, removable tray table. Note: Arm Supports required for use with Tray Table	Size 1 (for size 1 Sharky) 6085-7100-000	Size 2 (for size 2 Sharky) 6085-7200-000	



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com

A large, abstract image of a Sharky product, likely a baby walker or stroller, in blue and white. The word "SHARKY" is printed in blue, uppercase letters on a white diagonal band across the product.

SHARKY

DISTRIBUTED BY:

Script: Sharky_20250120_AU