

Simple Stuff Works



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

July 2025

Client details

Client/User name:			
DOB:		Date:	
Therapist:			

Notes

Further Information

CHOOSE A **PACKAGE** TO USE ON **FOAM MATTRESS**

Packages

Infant Package	Child Package	Adult Package
For infants from 0 to 5 years of age • 2x cot Bed Mesh • 2x cot Toppers • 2x cot cotton sheets • 2x small Lateral Supports • 2x small Wedges • Small Horseshoe pillow • Simple Stabiliser Size 1 • Travel bag 7150-7001-C00	For children from 5 to 12 years of age • 2x single Bed Mesh • 2x single Toppers • 2x single cotton sheets • 2x medium Lateral Supports • 2x small Wedges • 2x small Lateral Supports • Simple Stabiliser Size 2 • Small Terry Towelling Pillow • Travel bag 7150-7002-S00	For adults from 12+ years of age • 2x single Bed Mesh • 2x single Toppers • 2x single cotton sheets • 2x medium Lateral Supports • 2x medium Wedges • 2x large Lateral Supports • Large Terry Towelling Pillow • Travel bag 7150-7003-S00



OR REFER TO **PAGE 3** TO **CONFIGURE** YOUR OWN

Sheet Sets

Cot Bed	Single Bed	
Set Includes: • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x cot temperate toppers • 2 x cot fitted sheets 7150-7021-C00	Set Includes: • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x single temperate toppers • 2 x single fitted sheets 7150-7022-S00	

Sheets


Stabilising Mesh (Set of 2)		
Cot Bed Base layer: 1500 mm (L) x 700 mm (W) Top layer: 1800 mm (L) x 1300 mm (W) 7150-7011-C00	Single Bed Base layer: 1800 mm (L) x 1000 mm (W) Top layer: 1800 mm (L) x 1500 mm (W) 7150-7012-S00	

Temperate Topper (Set of 2)		
Cot Bed 1000 mm (L) x 590 mm (W) 7150-7051-C00	Single Bed 1000 mm (L) x 1000 mm (W) 7150-7052-S00	

Fitted Sheets			
Cot Bed 1400 mm (L) x 700 mm (W) 7150-7121-C00	Single Bed/King Single 1900 mm (L) x 900 mm (W) 7150-7122-S00	Double/Queen Bed 2030 mm (L) x 1530 mm (W) 7150-7123-S00	


Flat Sheet
Single 1620 mm (L) x 2600 mm (W) 7150-7126-S00

Wedges

Wedges (pair)		
With padded covers.		
Small A: 220 mm B: 110 mm C: 110 mm 7150-7351-000	Medium A: 360 mm B: 80 mm C: 90 mm 7150-7352-000	

Lateral Supports

Padded Lateral Supports		
Also available in an option for air mattresses, designed to grip directly to the mattress surface.		
Small (Foam) 7150-7201-000	Medium (Foam) 7150-7202-000	Large (Foam) 7150-7203-000
Small (Air) 7150-7211-000	Medium (Air) 7150-7212-000	Large (Air) 7150-7213-000
A: 180 mm B: 120 mm C: 180 mm	A: 220 mm B: 140 mm C: 220 mm	A: 280 mm B: 220 mm C: 280 mm




Simple Stabiliser

Supports hips and thighs in neutral supine alignment			
Size 1	Size 2	Size 3	Size 4
H: 110mm D: 110mm W: 420mm To fit approx leg circumference of 250mm. Suggested age range 2 – 4yrs.	H: 190mm D: 350mm W: 460mm To fit approx leg circumference of 350mm. Suggested age range 2 – 10yrs.	H: 250mm D: 400mm W: 650mm To fit approx leg circumference of 550mm. Suggested age range 5 – 13yrs.	H: 280mm D: 400mm W: 650mm To fit approx leg circumference of 600mm. Suggested age range 13 – adult.
With breathable inner & terry cotton outer 7150-7341-000	With breathable inner & terry cotton outer 7150-7342-000	With breathable inner & terry cotton outer 7150-7343-000	With breathable inner & terry cotton outer 7150-7344-000
With wipeable incontinence inner & terry cotton outer 7150-7345-000	With wipeable incontinence inner & terry cotton outer 7150-7346-000	With wipeable incontinence inner & terry cotton outer 7150-7347-000	With wipeable incontinence inner & terry cotton outer 7150-7348-000
Spare terry cotton cover 7150-7221-000	Spare terry cotton cover 7150-7222-000	Spare terry cotton cover 7150-7223-000	Spare terry cotton cover 7150-7224-000


Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3
A: 270mm B: 380mm C: 280mm D: 80mm E: 130mm F: 80mm	A: 400mm B: 560mm C: 340mm D: 100mm E: 200mm F: 90mm	A: 400mm B: 610mm C: 400mm D: 110mm E: 150mm F: 70mm
With breathable inner & terry cotton outer (standard) 7150-7971-000	With breathable inner & terry cotton outer (standard) 7150-7972-000	With breathable inner & terry cotton outer (standard) 7150-7973-000
With wipeable incontinence inner & terry cotton outer 7150-7981-000	With wipeable incontinence inner & terry cotton outer 7150-7982-000	With wipeable incontinence inner & terry cotton outer 7150-7983-000
		
Spare terry cotton cover 7150-7991-000	Spare terry cotton cover 7150-7992-000	Spare terry cotton cover 7150-7993-000

Foot Support

Provides positioning & skin protection for the feet (pair)

Size 1 A: 210mm (max) B: 200mm 7150-5351-000	Size 2 A: 250mm (max) B: 220mm 7150-5352-000	Size 3 A: 300mm (max) B: 300mm 7150-5353-000	
--	--	--	---

Pillows

Roll Pillow

Small A: 150mm B: 400mm	Large A: 180mm B: 500mm	
With breathable inner & terry cotton outer 7150-7851-000	With breathable inner & terry cotton outer 7150-7852-000	
With wipeable incontinence inner & terry cotton outer 7150-7855-000	With wipeable incontinence inner & terry cotton outer 7150-7856-000	
Spare terry cotton cover 7150-7857-000	Spare terry cotton cover 7150-7858-000	

Neck Support Pillow

(Suitable for air & foam mattresses). 300mm x 350mm x 100mm 7150-4050-000	
Spare terry cotton cover 7150-4055-000	

Horseshoe Pillow

Small A: 460mm B: 190mm	Large A: 500mm B: 200mm	
With breathable inner & terry cotton outer 7150-7711-000	With breathable inner & terry cotton outer 7150-7712-000	
With wipeable incontinence inner & terry cotton outer 7150-7721-000	With wipeable incontinence inner & terry cotton outer 7150-7722-000	
Spare terry cotton cover 7150-7731-000	Spare terry cotton cover 7150-7732-000	

Terry Towling Pillow

Small L: 480mm W: 350mm 7150-7501-000	Large L: 680mm W: 450mm 7150-7502-000	
Spare terry cotton cover 7150-7511-000	Spare terry cotton cover 7150-7512-000	



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: SimpleStuffWorks_20250707