

Sirocco Power Add-Ons



The Sirocco Power Add-On is a revolutionary solution for wheelchair users looking to enhance their mobility with power, control, and confidence. Unlike conventional single-wheel power add-ons, the Sirocco's unique dual-wheel design provides superior stability, traction and handling, making it ideal for both urban and off-road environments. Designed and manufactured in France, the Sirocco range features cutting-edge technology, including electromagnetic braking, active cruise control, and an intuitive HD colour display. Whether navigating city streets

or tackling challenging terrain, the Sirocco Power Add-On ensures a smooth and safe ride. With its easy-to-use universal attachment system, it seamlessly connects to rigid and folding manual wheelchairs without the need for modifications or adding any extra bracketry to the wheelchair. Available in multiple models, the Sirocco Extreme offers high-performance power with two 1500W motors, making it perfect for off-road adventurers, while the Sirocco Voyageur delivers a versatile, everyday mobility solution.

FEATURES AND BENEFITS:

- **Dual-Wheel Design:** Offers enhanced stability and control, providing an alternative for users who want extra confidence in varied terrains.
- **Powerful Performance:** Equipped with dual motors (1500W/720W each) for exceptional off-road capability.
- **Minimalist Attachment:** Ultra-fast, clip-on attachment system leaves wheelchair with no extra bracketry installed.
- **Universal Compatibility:** Easily attaches to most manual wheelchairs, including folding models, without modifications.
- **Advanced Braking System:** Programmable electromagnetic braking with ABS for increased safety and responsiveness.
- **User-Friendly Controls:** HD colour LCD screen with active cruise control and customisable throttle configuration.
- **Long-Lasting Battery Life:** Up to 100km range, depending on terrain and usage conditions.
- **Designed for Adventure:** Built to handle rugged environments, perfect for outdoor enthusiasts and active users.

Sirocco Extreme Power Add-On

UNSTOPPABLE POWER. ULTIMATE ADVENTURE



Material	Aluminium and Steel
Motor	2 x 1500W motors
Wheel	2x 16" with pneumatic all terrain tyres
Display	HD Colour LCD Screen
Battery	2x Li-ion 48V 16AH (Charger included)
Handlebar	Aluminium - adjustable in angle and height
Brake	Electromagnetic brake + Hydraulic disc brake
Size	W 75cm x H 92cm X D 44cm
Connection	Lateral frame clamp
Weight	30kg without attachments
Range	Up to 100km (based on use on flat terrain, 80kg user travelling at a speed of 10km/h or less).
Max. Weight	200kg

PACKAGE INCLUDES:

- Wheelchair attachment
- Aluminium chassis and 2x shock absorbers
- HD colour control screen
- Active cruise control
- Automatic LED headlight with horn
- Reverse gear
- Programmable electromagnetic brake with ABS
- 2x 16" light alloy wheels with pneumatic all terrain tyres
- Mudguard
- 2x 1500 Watt motors
- Twist throttle and hydraulic brake with parking lock
- 2x Lithium Ion 696Wh batteries
- Battery charger

Sirocco Voyageur Power Add-On

EFFORTLESS MOBILITY. BOUNDLESS FREEDOM



Material	Aluminium and Steel
Motor	2 x 720W motors
Wheel	2x 12" with pneumatic all terrain tyres
Display	HD Colour LCD Screen
Battery	1x Li-ion 48V 16AH (2x batteries for long range model)
Handlebar	Aluminium - adjustable in angle and height
Brake	Electromagnetic brake + Hydraulic disc brake
Size	W 69cm x H 78cm X D 44cm
Connection	Lateral frame clamp
Weight	22kg without attachments
Range	Up to 100km (based on use on flat terrain, 80kg user travelling at a speed of 10km/h or less).
Max. Weight	120kg

- PACKAGE INCLUDES:
- Wheelchair attachment
 - HD colour screen
 - Active cruise control
 - Automatic LED headlight
 - Reverse gear
 - Programmable electromagnetic brake with ABS
 - 2x 12" light alloy wheels with pneumatic all terrain tyres
 - Fenders
 - 2x 720 Watt motors
 - Twist throttle and hydraulic brake with parking lock
 - 1x Li-ion 48V 16AH (2x batteries for long range model)
 - Battery charger

