

Smirthwaite Combi



BATHROOM CHAIR

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size Options

	Size 1	Size 2	Size 3
Code	6115-1101-000	6115-1102-000	6115-1103-000
Seat Unit Measurements			
Seat height (cm)	21.5-29	27.5-37.5	30-40
Seat width (cm)	29	32.5	36.5
Seat depth (cm)	20.5-29	24.5-35.5	31-38.5
Back support height (cm)	32.5	37	40
Load capacity (kg)	35	45	55
Product weight (kg)	11	13	14

Base package includes:

- Upholstered Seat Unit with Adjustable Armrests
- Velcro Pelvic Strap
- Plastic Commode Seat with Lid
- Moulded Potty Seat



Accessories

Hand Rail Toggle

Size 1 6115-7201-000

Size 2 6115-7202-000

Size 3 6115-7203-000




Splash Guard

Used in conjunction with Hand Rail Toggle.

Size 1 6115-7601-000

Size 2 6115-7602-000



Reversible Footboard

Size 1 6115-5021-000

Size 2 6115-5022-000

Size 3 6115-5023-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com

DISTRIBUTED BY:

Script: Combi_20250123