

spex

Shaping Better Lives



CLINICAL WHEELCHAIR SEATING

SCRIPTING AND ORDERING GUIDE

Spex Seat Base

BASE PACKAGE:

- Single piece seat base plate
- Either Quick-Release or Fixed-Mount hardware set
- Pair of retainer brackets

Weight: 2.4kg for Seat Base with Quick-Release, 2.3kg for Seat Base with Fixed Mount (16x18" 40x45cm)



2. SELECT HARDWARE

Fixed-Mount



1008-WWLL-FXB



Quick-Release



1008-WWLL-QRB

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
10" 25cm	12" 30cm	1008-1012-XXX	☐
	14" 36cm	1008-1014-XXX	☐
12" 30cm	12" 30cm	1008-1212-XXX	☐
	14" 35cm	1008-1214-XXX	☐
	16" 40cm	1008-1216-XXX	☐
14" 35cm	18" 45cm	1008-1218-XXX	☐
	14" 35cm	1008-1414-XXX	☐
	16" 40cm	1008-1416-XXX	☐
15" 38cm	18" 45cm	1008-1418-XXX	☐
	17" 43cm	1008-1517-XXX	☐
	20" 50cm	1008-1520-XXX	☐
16" 40cm	16" 40cm	1008-1616-XXX	☐
	18" 45cm	1008-1618-XXX	☐
	20" 50cm	1008-1620-XXX	☐
	22" 55cm	1008-1622-XXX	☐
17" 43cm	17" 43cm	1008-1717-XXX	☐
	20" 50cm	1008-1720-XXX	☐
	22" 55cm	1008-1722-XXX	☐
18" 45cm	18" 45cm	1008-1818-XXX	☐
	20" 50cm	1008-1820-XXX	☐
	22" 55cm	1008-1822-XXX	☐
20" 50cm	20" 50cm	1008-2020-XXX	☐
	22" 55cm	1008-2022-XXX	☐
22" 55cm	22" 55cm	1008-2222-XXX	☐
Up to 20" 50cm	Custom: _____"	1008-2929-XXX	☐
Up to 25" 63cm	Custom: _____"	1008-2939-XXX	☐



3. CUSTOMISE (IF REQUIRED)

✓	ITEM	NOTES
☐	Leg Length Discrepancy Available in 1" Increments 1008-0292-000	SELECT L R T = _____
☐	Custom backrest cane depth cut-outs Up to 2" additional depth, available in 1" Increments Standard Cut-Out = 3"x3" (8x8cm) 1008-0291-000	S = _____
☐	Additional Cushion retainer bracket (pair) Note: All bases come with a pair of cushion retainers 1008-0250-100	

WHEELCHAIR COMPATIBILITY

Seat Base Width
10" - 20"

Seat Base Width
21" +

Fixed Mount

-0" to +2"

-0" to +1"

Quick Release

-1 to +2"

-1 to +1"

Vigour Standard Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
8" 20cm	8" 20cm	1103-0808-OSC	☐
	10" 25cm	1103-0810-OSC	☐
10" 25cm	10" 25cm	1103-1010-OSC	☐
	12" 30cm	1103-1012-OSC	☐
	XX XX	1103-10XX-OSC	☐
11" 28cm	10" 25cm	1103-1110-OSC	☐
	12" 30cm	1103-1112-OSC	☐
	XX XX	1103-11XX-OSC	☐
12" 30cm	12" 30cm	1103-1212-OSC	☐
	14" 35cm	1103-1214-OSC	☐
	XX XX	1103-12XX-OSC	☐
13" 33cm	12" 30cm	1103-1312-OSC	☐
	14" 35cm	1103-1314-OSC	☐
	XX XX	1103-13XX-OSC	☐
14" 35cm	14" 35cm	1103-1414-OSC	☐
	16" 40cm	1103-1416-OSC	☐
	18" 45cm	1103-1418-OSC	☐
	XX XX	1103-14XX-OSC	☐
15" 38cm	14" 35cm	1103-1514-OSC	☐
	16" 40cm	1103-1516-OSC	☐
	XX XX	1103-15XX-OSC	☐
16" 40cm	16" 40cm	1103-1616-OSC	☐
	18" 45cm	1103-1618-OSC	☐
	20" 50cm	1103-1620-OSC	☐
17" 43cm	XX XX	1103-16XX-OSC	☐
	17" 43cm	1103-1717-OSC	☐
	18" 45cm	1103-1718-OSC	☐
18" 45cm	XX XX	1103-17XX-OSC	☐
	18" 45cm	1103-1818-OSC	☐
	20" 50cm	1103-1820-OSC	☐
19" 48cm	22" 55cm	1103-1822-OSC	☐
	XX XX	1103-18XX-OSC	☐
	18" 45cm	1103-1918-OSC	☐
20" 50cm	20" 50cm	1103-1920-OSC	☐
	XX XX	1103-19XX-OSC	☐
	18" 45cm	1103-2018-OSC	☐
22" 55cm	20" 50cm	1103-2020-OSC	☐
	22" 55cm	1103-2022-OSC	☐
	XX XX	1103-20XX-OSC	☐
24" 60cm	18" 45cm	1103-2218-OSC	☐
	20" 50cm	1103-2220-OSC	☐
	22" 55cm	1103-2222-OSC	☐
Up to 20" 50cm	18" 45cm	1103-2418-OSC	☐
Up to 25" 63cm	20" 50cm	1103-2420-OSC	☐
	22" 55cm	1103-2422-OSC	☐
Up to 20" 50cm	XX = _ _ _ _	1103-2929-OSC	☐
Up to 25" 63cm	XX = _ _ _ _	1103-2939-OSC	☐

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit



TROUGH HEIGHT = 3CM

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover (A) 1123-WWDD-6SC (insert width and depth relevant to cushion size)	
☐	Extra Incontinence Cover (B) 1133-WWDD-7SC (insert width and depth relevant to cushion size)	
☐	Rigidiser Base Panel 1110-WWDD-200 (can be inserted inside cover for a rigid surface)	
☐	5° Wedge 1209-2724-300	

4. CUSTOMISE

Refer to Cushion Customisations (page 11).

Vigour High Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
8" 20cm	8" 20cm	1103-0808-OHC	☐
	10" 25cm	1103-0810-OHC	☐
10" 25cm	10" 25cm	1103-1010-OHC	☐
	12" 30cm	1103-1012-OHC	☐
	XX XX	1103-10XX-OHC	☐
11" 28cm	10" 25cm	1103-1110-OHC	☐
	12" 30cm	1103-1112-OHC	☐
	XX XX	1103-11XX-OHC	☐
12" 30cm	12" 30cm	1103-1212-OHC	☐
	14" 35cm	1103-1214-OHC	☐
	XX XX	1103-12XX-OHC	☐
13" 33cm	12" 30cm	1103-1312-OHC	☐
	14" 35cm	1103-1314-OHC	☐
	XX XX	1103-13XX-OHC	☐
14" 35cm	14" 35cm	1103-1414-OHC	☐
	16" 40cm	1103-1416-OHC	☐
	18" 45cm	1103-1418-OHC	☐
	XX XX	1103-14XX-OHC	☐
15" 38cm	14" 35cm	1103-1514-OHC	☐
	16" 40cm	1103-1516-OHC	☐
	XX XX	1103-15XX-OHC	☐
16" 40cm	16" 40cm	1103-1616-OHC	☐
	18" 45cm	1103-1618-OHC	☐
	20" 50cm	1103-1620-OHC	☐
17" 43cm	XX XX	1103-16XX-OHC	☐
	17" 43cm	1103-1717-OHC	☐
	18" 45cm	1103-1718-OHC	☐
18" 45cm	XX XX	1103-17XX-OHC	☐
	18" 45cm	1103-1818-OHC	☐
	20" 50cm	1103-1820-OHC	☐
19" 48cm	22" 55cm	1103-1822-OHC	☐
	XX XX	1103-18XX-OHC	☐
	18" 45cm	1103-1918-OHC	☐
20" 50cm	20" 50cm	1103-1920-OHC	☐
	XX XX	1103-19XX-OHC	☐
	18" 45cm	1103-2018-OHC	☐
22" 55cm	20" 50cm	1103-2020-OHC	☐
	22" 55cm	1103-2022-OHC	☐
	XX XX	1103-20XX-OHC	☐
24" 60cm	18" 45cm	1103-2218-OHC	☐
	20" 50cm	1103-2220-OHC	☐
	22" 55cm	1103-2222-OHC	☐
Up to 20" 50cm	18" 45cm	1103-2418-OHC	☐
Up to 25" 63cm	20" 50cm	1103-2420-OHC	☐
	22" 55cm	1103-2422-OHC	☐
Up to 20" 50cm	XX = _ _ _ _	1103-2929-OHC	☐
Up to 25" 63cm	XX = _ _ _ _	1103-2939-OHC	☐

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit



TROUGH HEIGHT = 5CM

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover (A) 1123-WWDD-6HC (insert width and depth relevant to cushion size)	
☐	Extra Incontinence Cover (B) 1133-WWDD-7HC (insert width and depth relevant to cushion size)	
☐	Rigidiser Base Panel 1110-WWDD-200 (can be inserted inside cover for a rigid surface)	
☐	5° Wedge 1209-2724-300	

4. CUSTOMISE

Refer to Cushion Customisations (page 11).

Spex Standard Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
8" 20cm	8" 20cm	1105-0808-OSC	☐
	10" 25cm	1105-0810-OSC	☐
10" 25cm	10" 25cm	1105-1010-OSC	☐
	12" 30cm	1105-1012-OSC	☐
	XX XX	1105-10XX-OSC	☐
11" 28cm	10" 25cm	1105-1110-OSC	☐
	12" 30cm	1105-1112-OSC	☐
	XX XX	1105-11XX-OSC	☐
12" 30cm	12" 30cm	1105-1212-OSC	☐
	14" 35cm	1105-1214-OSC	☐
	XX XX	1105-12XX-OSC	☐
13" 33cm	12" 30cm	1105-1312-OSC	☐
	14" 35cm	1105-1314-OSC	☐
	XX XX	1105-13XX-OSC	☐
14" 35cm	14" 35cm	1105-1414-OSC	☐
	16" 40cm	1105-1416-OSC	☐
	18" 45cm	1105-1418-OSC	☐
	XX XX	1105-14XX-OSC	☐
15" 38cm	14" 35cm	1105-1514-OSC	☐
	16" 40cm	1105-1516-OSC	☐
	XX XX	1105-15XX-OSC	☐
16" 40cm	16" 40cm	1105-1616-OSC	☐
	18" 45cm	1105-1618-OSC	☐
	20" 50cm	1105-1620-OSC	☐
17" 43cm	XX XX	1105-16XX-OSC	☐
	17" 43cm	1105-1717-OSC	☐
	18" 45cm	1105-1718-OSC	☐
18" 45cm	XX XX	1105-17XX-OSC	☐
	18" 45cm	1105-1818-OSC	☐
	20" 50cm	1105-1820-OSC	☐
19" 48cm	22" 55cm	1105-1822-OSC	☐
	XX XX	1105-18XX-OSC	☐
	18" 45cm	1105-1918-OSC	☐
20" 50cm	20" 50cm	1105-1920-OSC	☐
	XX XX	1105-19XX-OSC	☐
	18" 45cm	1105-2018-OSC	☐
22" 55cm	20" 50cm	1105-2020-OSC	☐
	22" 55cm	1105-2022-OSC	☐
	XX XX	1105-20XX-OSC	☐
24" 60cm	18" 45cm	1105-2218-OSC	☐
	20" 50cm	1105-2220-OSC	☐
	22" 55cm	1105-2222-OSC	☐
Up to 20" 50cm	18" 45cm	1105-2418-OSC	☐
Up to 25" 63cm	20" 50cm	1105-2420-OSC	☐
	22" 55cm	1105-2422-OSC	☐
Up to 20" 50cm	XX = _ _ _ _	1105-2929-OSC	☐
Up to 25" 63cm	XX = _ _ _ _	1105-2939-OSC	☐

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit
4. Strategic Positioning Base



TROUGH HEIGHT = 3CM

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover (A) 1125-WWDD-6SC (insert width and depth relevant to cushion size)	
☐	Extra Incontinence Cover (B) 1135-WWDD-7SC (insert width and depth relevant to cushion size)	
☐	Rigidiser Base Panel 1110-WWDD-200 (can be inserted inside cover for a rigid surface)	
☐	Extra Positioning Kit 1149-XXXX-062	
☐	5° Wedge 1209-2724-300	

4. CUSTOMISE

Refer to Cushion Customisations (page 11).

Spex High Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
8" 20cm	8" 20cm	1105-0808-OHC	☐
	10" 25cm	1105-0810-OHC	☐
10" 25cm	10" 25cm	1105-1010-OHC	☐
	12" 30cm	1105-1012-OHC	☐
	XX XX	1105-10XX-OHC	☐
11" 28cm	10" 25cm	1105-1110-OHC	☐
	12" 30cm	1105-1112-OHC	☐
	XX XX	1105-11XX-OHC	☐
12" 30cm	12" 30cm	1105-1212-OHC	☐
	14" 35cm	1105-1214-OHC	☐
	XX XX	1105-12XX-OHC	☐
13" 33cm	12" 30cm	1105-1312-OHC	☐
	14" 35cm	1105-1314-OHC	☐
	XX XX	1105-13XX-OHC	☐
14" 35cm	14" 35cm	1105-1414-OHC	☐
	16" 40cm	1105-1416-OHC	☐
	18" 45cm	1105-1418-OHC	☐
	XX XX	1105-14XX-OHC	☐
15" 38cm	14" 35cm	1105-1514-OHC	☐
	16" 40cm	1105-1516-OHC	☐
	XX XX	1105-15XX-OHC	☐
16" 40cm	16" 40cm	1105-1616-OHC	☐
	18" 45cm	1105-1618-OHC	☐
	20" 50cm	1105-1620-OHC	☐
17" 43cm	XX XX	1105-16XX-OHC	☐
	17" 43cm	1105-1717-OHC	☐
	18" 45cm	1105-1718-OHC	☐
18" 45cm	XX XX	1105-17XX-OHC	☐
	18" 45cm	1105-1818-OHC	☐
	20" 50cm	1105-1820-OHC	☐
19" 48cm	22" 55cm	1105-1822-OHC	☐
	XX XX	1105-18XX-OHC	☐
	18" 45cm	1105-1918-OHC	☐
20" 50cm	20" 50cm	1105-1920-OHC	☐
	XX XX	1105-19XX-OHC	☐
	18" 45cm	1105-2018-OHC	☐
22" 55cm	20" 50cm	1105-2020-OHC	☐
	22" 55cm	1105-2022-OHC	☐
	XX XX	1105-20XX-OHC	☐
24" 60cm	18" 45cm	1105-2218-OHC	☐
	20" 50cm	1105-2220-OHC	☐
	22" 55cm	1105-2222-OHC	☐
Up to 20" 50cm	18" 45cm	1105-2418-OHC	☐
Up to 25" 63cm	20" 50cm	1105-2420-OHC	☐
	22" 55cm	1105-2422-OHC	☐
Up to 20" 50cm	XX = _ _ _ _	1105-2929-OHC	☐
Up to 25" 63cm	XX = _ _ _ _	1105-2939-OHC	☐

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit
4. Strategic Positioning Base



TROUGH HEIGHT = 5CM

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover (A) 1125-WWDD-6HC (insert width and depth relevant to cushion size)	
☐	Extra Incontinence Cover (B) 1135-WWDD-7HC (insert width and depth relevant to cushion size)	
☐	Rigidiser Base Panel 1110-WWDD-200 (can be inserted inside cover for a rigid surface)	
☐	Extra Positioning Kit 1149-XXXX-060	
☐	5° Wedge 1209-2724-300	

4. CUSTOMISE

Refer to Cushion Customisations (page 11).

Spex SuperHigh Contour Cushion

1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
8" 20cm	8" 20cm	1105-0808-OSS	☐
	10" 25cm	1105-0810-OSS	☐
10" 25cm	10" 25cm	1105-1010-OSS	☐
	12" 30cm	1105-1012-OSS	☐
	XX XX	1105-10XX-OSS	☐
11" 28cm	10" 25cm	1105-1110-OSS	☐
	12" 30cm	1105-1112-OSS	☐
	XX XX	1105-11XX-OSS	☐
12" 30cm	12" 30cm	1105-1212-OSS	☐
	14" 35cm	1105-1214-OSS	☐
	XX XX	1105-12XX-OSS	☐
13" 33cm	12" 30cm	1105-1312-OSS	☐
	14" 35cm	1105-1314-OSS	☐
	XX XX	1105-13XX-OSS	☐
14" 35cm	14" 35cm	1105-1414-OSS	☐
	16" 40cm	1105-1416-OSS	☐
	18" 45cm	1105-1418-OSS	☐
	XX XX	1105-14XX-OSS	☐
15" 38cm	14" 35cm	1105-1514-OSS	☐
	16" 40cm	1105-1516-OSS	☐
	XX XX	1105-15XX-OSS	☐
16" 40cm	16" 40cm	1105-1616-OSS	☐
	18" 45cm	1105-1618-OSS	☐
	20" 50cm	1105-1620-OSS	☐
	XX XX	1105-16XX-OSS	☐
17" 43cm	17" 43cm	1105-1717-OSS	☐
	18" 45cm	1105-1718-OSS	☐
	XX XX	1105-17XX-OSS	☐
18" 45cm	18" 45cm	1105-1818-OSS	☐
	20" 50cm	1105-1820-OSS	☐
	22" 55cm	1105-1822-OSS	☐
19" 48cm	XX XX	1105-18XX-OSS	☐
	18" 45cm	1105-1918-OSS	☐
	20" 50cm	1105-1920-OSS	☐
20" 50cm	XX XX	1105-19XX-OSS	☐
	18" 45cm	1105-2018-OSS	☐
	20" 50cm	1105-2020-OSS	☐
	22" 55cm	1105-2022-OSS	☐
22" 55cm	XX XX	1105-20XX-OSS	☐
	18" 45cm	1105-2218-OSS	☐
	20" 50cm	1105-2220-OSS	☐
24" 60cm	22" 55cm	1105-2222-OSS	☐
	18" 45cm	1105-2418-OSS	☐
	20" 50cm	1105-2420-OSS	☐
Up to 20" 50cm	22" 55cm	1105-2422-OSS	☐
Up to 25" 63cm	18" 45cm	1105-2929-OSS	☐
	20" 50cm	1105-2939-OSS	☐

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit
4. Strategic Positioning Base



10 - 13" TROUGH HEIGHT = 8CM
14"+ TROUGH HEIGHT = 10CM

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover (A) 1125-WWDD-6SS (insert width and depth relevant to cushion size)	
☐	Extra Incontinence Cover (B) 1135-WWDD-7SS (insert width and depth relevant to cushion size)	
☐	Rigidiser Base Panel 1110-WWDD-200 (can be inserted inside cover for a rigid surface)	
☐	Extra Positioning Kit 1149-XXXX-063	
☐	5° Wedge 1209-2724-300	

4. CUSTOMISE


Refer to Cushion Customisations (page 11).

Spex Constructa Cushion

PACKAGE INCLUDES:

1. Outer Cover
2. Overlay
3. Inner Cover
4. Base Foam:
 - a. Contours
 - b. Ischial Insert
 - c. Positioning Blocks
 - d. Base
5. Adjustment Kit

1. SELECT CUSHION



Includes Strategic Positioning Kit.

Includes Pelvic Adjustment Kit.

-1XXX

-FXXX

With Spare Outer Cover: -XCXX

Without Spare Outer Cover: -X0XX

2. SELECT CONTOURING



-XXSC

-XXHC

-XXSS

3. SELECT SIZE

WIDTH	LENGTH	CODE	✓
14" 35cm	14" 35cm	1102-1414-XXXX	☐
	16" 40cm	1102-1416-XXXX	☐
	18" 45cm	1102-1418-XXXX	☐
16" 40cm	16" 40cm	1102-1616-XXXX	☐
	18" 45cm	1102-1618-XXXX	☐
	20" 50cm	1102-1620-XXXX	☐
18" 45cm	18" 45cm	1102-1818-XXXX	☐
	20" 50cm	1102-1820-XXXX	☐
20" 50cm	20" 50cm	1102-2020-XXXX	☐
	22" 55cm	1102-2022-XXXX	☐

4. SELECT ACCESSORIES

✓	ACCESSORY	IMAGE
☐	Extra Outer Cover, Standard or High Contour Shape: 1102-WWDD-00C-C Flex: 1102-WWDD-F0C-C (insert width relevant to cushion size)	
☐	Extra Outer Cover, Superhigh Contour Shape: 1102-WWDD-00S-C Flex: 1102-WWDD-F0S-C (insert width relevant to cushion size)	
☐	Flexy Contour 1102-WW00-00F (insert width relevant to cushion size)	
☐	Constructa Flex Extra Pelvic Adjustment Kit 1102-WW00-F00 (insert width relevant to cushion size)	
☐	Constructa Shape Extra Strategic Positioning Kit 1102-WW00-000 (insert width relevant to cushion size)	
☐	Rigidiser Base Panel (can be inserted inside cover for a rigid surface) 1110-WWDD-200	
☐	5° Wedge 1209-2724-300	
☐	Extra 7.5 Degree Wedge (Flex 1 leg only) 1102-WWDD-321	

5. CUSTOMISE

Customise your Constructa Cushion with extra components available from our **Constructa Cushion Spare Parts Script**.

Spex Wonderseat Constructa Cushion

PACKAGE INCLUDES:

1. Outer Cover
2. Inner Cover
3. Overlay
4. Base Foam:
 - a. Contours
 - b. Ischial Insert
 - c. Positioning Blocks
 - d. Base
5. Adjustment Kit (not pictured)



2. SELECT COLOUR

Your colour selection applies to any extra covers ordered



1. SELECT SIZE

WIDTH		LENGTH		CODE	✓
10"	25cm	12"	30cm	1102-1012-WSC	☐
12"	30cm	16"	40cm	1102-1216-WSC	☐

3. SELECT ACCESSORIES

✓ ACCESSORY

IMAGE

- Extra Outer Cover
- ☐ 10x12": 1102-1012-WSC-C
- ☐ 12x16": 1102-1216-WSC-C



XLella Cushion

PACKAGE INCLUDES:

1. Outer Cover
2. Incontinence Cover
3. Cushion Foam Unit
4. Strategic Positioning Base (Spex Classic Cushions only)



1. SELECT XLELLA CUSHION STYLE

Vigour
Standard Contour



1103-WWDD-OSC

Spex Classic
Standard Contour




1105-WWDD-OSC

Spex Classic
High Contour



1105-WWDD-OHC



2. SELECT SIZING

To order replace
WW and DD in code

		Cushion Width (WW) (inches)					
		20	22	24	26	28	30
Cushion Depth (DD) (inches)	18						
	19						
	20						
	21						
	22						
	23						
	24						
	25						
	26						

XLella Cushion Size Available

XLella Cushion Size Made
To Order - Cushion Depth
can be specified to suit

3. SELECT COLOUR

Your colour selection applies to
any extra covers ordered

Black
(default)

Chilli
1209-2491-022

Mango
1209-2491-030

Bumblebee
1209-2491-044

Calypso
1209-2491-028

Royal
1209-2491-077

Ocean
1209-2491-060

Leaf
1209-2491-054

Granite
1209-2491-090

4. SELECT ACCESSORIES

✓ ACCESSORY

IMAGE

Extra Outer Cover (A)


- ☐ Vigour Std: 1123-WWDD-6SC
- ☐ Classic Std: 1125-WWDD-6SC
- ☐ Classic High: 1125-WWDD-6HC
(insert width & depth as relevant)

Extra Incontinence Cover (B)

- ☐ Vigour Std: 1133-WWDD-7SC
- ☐ Classic Std: 1135-WWDD-7SC
- ☐ Classic High: 1135-WWDD-7HC
(insert width & depth as relevant)

Extra Positioning Kit

(Spex Classic Cushions only)
1149-XXXX-062



5. CUSTOMISE


✓ CUSTOMISATION

IMAGE

Remove Medial Thigh Contouring

1119-0282-000


- ☐ Achieve flat front surface without any contouring. Ideal relief for excess tissue and to facilitate easy transfers.



Front Surface Cut-Back

M = _____

- ☐ 1119-0290-000
Sloping front for users with excess tissue or to provide relief when transferring and/or foot propelling.



Cushion Customisations

☐ Medial Thigh Support - Specified measurements 1119-0281-000	Measurements (cm) E _____ F _____ G _____ H _____			Reference
	Pre-set (cm) Standard ☐ -0283- High ☐ -0284- SuperHigh ☐ -0285-			

Leg Length Discrepancy 1119-0293-000	Users Left (cm) L _____	Users Right (cm) L _____	Reference
--	----------------------------	-----------------------------	----------------------

Contacture Undercut - Front of Cushion 1119-0290-000	Users Left (cm) M _____	Users Right (cm) M _____	Reference
--	----------------------------	-----------------------------	----------------------

☐ Hip Belt Cut Out 1119-0294-000	Users Left (cm) ☐ N _____ O _____	Users Right (cm) ☐ N _____ O _____	Reference
--	---	--	----------------------

☐ Angled Shelf Cushion Base 45° shelf starts where ischial well stops. 1119-0261-000	Measurements (cm) R _____	Reference
---	------------------------------	----------------------

½" Gel Foam Overlay
(N/A for cushions 22" x 22" and greater).

The soft upper layer will be replaced with this overlay.
Max Size = 22" x 32"
1119-0391-000

½" Soft Pressure Relieving Foam Overlay

½" Soft Pressure Relieving Foam.
Max Size = 32" x 32"
1119-0390-000


1" Soft Gel Foam Ischial Well Insert

1119-0351-000


Universal Positioning Base

1. SELECT BASE


CUSHION SPACER KIT



ISCHIAL CONTOUR CUSHION SPACER KIT



HIP FLEXION CONTOUR SPACER KIT



PACKAGE INCLUDES:

10mm spacer, 20mm spacer
8 x obliquity spacer squares
Raises cushion for extra long leg lengths,
or as a base for attaching contouring
support wedges.

10mm spacer base
25mm front ischial ledge
8 x obliquity spacer squares
Helps to prevent posterior pelvic sliding.

Base with 25mm high ischial ledge
8 x obliquity spacer squares
Helps to prevent posterior pelvic sliding.
Allows for 10° Hip Flexion relief.

2. SELECT SIZE

	WIDTH	LENGTH	CODE	✓
Cushion spacer kit	14" 35cm	14" 35cm	1113-01-1414	☐
		16" 40cm	1113-01-1416	☐
		18" 45cm	1113-01-1418	☐
	16" 40cm	16" 40cm	1113-01-1616	☐
		18" 45cm	1113-01-1618	☐
		20" 50cm	1113-01-1620	☐
	18" 45cm	18" 45cm	1113-01-1818	☐
		20" 50cm	1113-01-1820	☐
	20" 50cm	18" 45cm	1113-01-2018	☐
Ischial contour cushion Spacer kit	14" 35cm	14" 35cm	1113-02-1414	☐
		16" 40cm	1113-02-1416	☐
		18" 45cm	1113-02-1418	☐
	16" 40cm	16" 40cm	1113-02-1616	☐
		18" 45cm	1113-02-1618	☐
		20" 50cm	1113-02-1620	☐
	18" 45cm	18" 45cm	1113-02-1818	☐
		20" 50cm	1113-02-1820	☐
	20" 50cm	18" 45cm	1113-02-2018	☐
Hip flexion contour Spacer kit	14" 35cm	14" 35cm	1113-03-1414	☐
		16" 40cm	1113-03-1416	☐
		18" 45cm	1113-03-1418	☐
	16" 40cm	16" 40cm	1113-03-1616	☐
		18" 45cm	1113-03-1618	☐
		20" 50cm	1113-03-1620	☐
	18" 45cm	18" 45cm	1113-03-1818	☐
		20" 50cm	1113-03-1820	☐
	20" 50cm	18" 45cm	1113-03-2018	☐
		20" 50cm	1113-03-2020	☐


3. SELECT CONTOUR KIT

LENGTH	CODE	✓
14" 35cm	1113-05-0014	☐
16" 40cm	1113-05-0016	☐
18" 45cm	1113-05-0018	☐
20" 50cm	1113-05-0020	☐
XX = ----	1113-05-0029	☐



4. SELECT COVER COMBINATION

WIDTH	LENGTH	CODE	✓
14" 35cm	14" 35cm	1113-07-1414	☐
	16" 40cm	1113-07-1416	☐
	18" 45cm	1113-07-1418	☐
16" 40cm	16" 40cm	1113-07-1616	☐
	18" 45cm	1113-07-1618	☐
	20" 50cm	1113-07-1620	☐
18" 45cm	18" 45cm	1113-07-1818	☐
	20" 50cm	1113-07-1820	☐
20" 50cm	18" 45cm	1113-07-2018	☐
	20" 50cm	1113-07-2020	☐
XX = ----	XX = ----	1113-07-2929	☐



Zygo Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
 2. Foam System
 3. Outer Cover
 4. Foam Wedges
 5. Cover Extension
- (Note: hardware must be selected)

1. SELECT ZYGO CONTOUR OPTION



Active Lo
1251-WWHH-OAL



Active Mid
1251-WWHH-OAM



Deep Lo
1251-WWHH-ODL



Deep Mid
1251-WWHH-ODM

2. SELECT SIZING

To order replace
WW and HH in code

		Wheelchair Width (inches)					
		10	12	14	16	18	20
Back Support Height (inches)	7						
	8						
	9						
	10						
	11						
	12						
	13						
	14						
	15						
	16						
	17						
	18						

Lo Back Support

Mid Back Support

3. SELECT HARDWARE - DROP MOUNT

		Back Support Width (inches)					
		10	12	14	16	18	20
Back Support Height (inches)	7		P				
	8			S/M			
	9	P			S/M	S/M	
	10		P				
	11	P					M/L
	12		P				
	13		P	S/M			
	14				S/M		
	15			S/M		M/L	M/L
	16			S/M	S/M		
	17				S/M	M/L	
	18					M/L	M/L



Paediatric
1259-1400-OOP



Small
1259-1400-OOS



Medium
1259-1400-00M



Large
1259-1400-00L

4. ADD ACCESSORIES




TSS Lite contouring & positioning kit
1250-3955-000



Harness Guides
1400-6692-000

5. BUILD CODE

REFER TO THE ZYGO BACK SUPPORT CODES ABOVE



Vigour Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
 2. Foam System
 3. Outer Cover
 4. Posterior Pelvic Support & Positioning Wedges
- (Note: hardware must be selected)



1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
10" 25cm	12" 30cm	1252-1012-00H	☐
	14" 35cm	1252-1014-00H	☐
	16" 40cm	1252-1016-00H	☐
11" 28cm	12" 30cm	1252-1112-00H	☐
	14" 35cm	1252-1114-00H	☐
	16" 40cm	1252-1116-00H	☐
12" 30cm	14" 35cm	1252-1214-00H	☐
	16" 40cm	1252-1216-00H	☐
	18" 45cm	1252-1218-00H	☐
13" 33cm	14" 35cm	1252-1314-00H	☐
	16" 40cm	1252-1316-00H	☐
	18" 45cm	1252-1318-00H	☐
14" 35cm	16" 40cm	1252-1416-00H	☐
	18" 45cm	1252-1418-00H	☐
	20" 50cm	1252-1420-00H	☐
15" 38cm	16" 40cm	1252-1516-00H	☐
	18" 45cm	1252-1518-00H	☐
	20" 50cm	1252-1520-00H	☐
16" 40cm	18" 45cm	1252-1618-00H	☐
	20" 50cm	1252-1620-00H	☐
	23" 58cm	1252-1623-00H	☐
17" 43cm	18" 45cm	1252-1718-00H	☐
	20" 50cm	1252-1720-00H	☐
	23" 58cm	1252-1723-00H	☐
18" 45cm	18" 45cm	1252-1818-00H	☐
	20" 50cm	1252-1820-00H	☐
	23" 58cm	1252-1823-00H	☐
19" 48cm	18" 45cm	1252-1918-00H	☐
	20" 50cm	1252-1920-00H	☐
	23" 58cm	1252-1923-00H	☐
20" 50cm	18" 45cm	1252-2018-00H	☐
	20" 50cm	1252-2020-00H	☐
	23" 58cm	1252-2023-00H	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black
(default)

Chilli
1209-2491-022

Mango
1209-2491-030

Bumblebee
1209-2491-044

Calypso
1209-2491-028


Royal
1209-2491-077

Ocean
1209-2491-060


Leaf
1209-2491-054

Granite
1209-2491-090


3. SELECT HARDWARE




Single Mount
1259-2100-100




Dual Mount
1259-2200-100




Quad Mount
1259-2400-100




Quad Mount,
Heavy Duty,
up to 23" wide
1259-2405-100




Quad Mount,
Heavy Duty, 23"
wide +
1259-2406-100



Direct Mount
1259-2460-100




Drop Mount
Paediatric:
1259-1400-00P
Small:
1259-1400-00S
Medium:
1259-1400-00M
Large:
1259-1400-00L



Quad Extension
Brackets
1259-2450-100


4. ADD ACCESSORIES




Additional backrest
lumber pad and lateral
stability wedges
1209-2713-300



Asymmetric narrow
contour wedge
(200mm x 100mm x
40mm)
1209-2722-300




Asymmetric bi-angular
contour wedge
(200mm x 100mm x
40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300



Extra Outer Cover (A)
Extra Incontinence
Cover (B)
1209-2713-300



T-shape
modification
1209-2931-000
Surface / I Shape
modification
1209-2936-000

Vigour Mid Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
2. Foam System
3. Outer Cover
4. Lumbar Pad and Positioning Wedges


(Note: hardware must be selected)




1. SELECT SIZE

WIDTH		LENGTH		CODE	✓
10"	25cm	10"	25cm	1252-1010-00M	☐
12"	30cm	12"	30cm	1252-1212-00M	☐
14"	35cm	14"	35cm	1252-1414-00M	☐
16"	40cm	16"	40cm	1252-1616-00M	☐
17"	43cm	17" 43cm		1252-1717-00M	☐
18"	45cm			1252-1817-00M	☐
20"	50cm			1252-2017-00M	☐
User-matched size up to 20" wide				1252-2929-00M	☐


2. SELECT HARDWARE




Single Mount
1259-2100-100




Dual Mount
1259-2200-100




Quad Mount
1259-2400-100




Quad Mount, Heavy Duty, up to 23" wide
1259-2405-100




Quad Mount, Heavy Duty, 23" wide +
1259-2406-100




Direct Mount
1259-2460-100



Drop Mount
Paediatric:
1259-1400-00P
Small:
1259-1400-00S
Medium:
1259-1400-00M
Large:
1259-1400-00L




Quad Extension Brackets
1259-2450-100




Width Adjustment Bracket Set
1259-2490-100


3. ADD ACCESSORIES




Additional backrest lumber pad and lateral stability wedges
1209-2713-300



Asymmetric narrow contour wedge
(200mm x 100mm x 40mm)
1209-2722-300



Asymmetric bi-angular contour wedge
(200mm x 100mm x 40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300



Extra Outer Cover (A)
1252-WWLL-00M-I
Extra Incontinence Cover (B)
1252-WWLL-00M-C

Growth-Adjust Back Support

GROWTH-ADJUST



SPEX CLASSIC GROWTH-ADJUST



SPEX CLASSIC T SHAPE GROWTH-ADJUST



PACKAGE INCLUDES:

1. Growth-Adjust Aluminium Shell
2. Foam System
3. Outer Cover
4. Growth Kit of foam parts
5. Posterior Pelvic Support & Positioning Wedges

(Note: hardware must be selected)



1. SELECT SIZE

GROWTH-ADJUST

WIDTH	LENGTH	CODE	✓
12"	11"-15" 28-38cm	1252-1211-OGA	☐
30cm	14"-18" 35-45cm	1252-1214-OGA	☐
14"	12"-17" 30-43cm	1252-1412-OGA	☐
35cm	16"-21" 40-53cm	1252-1416-OGA	☐
16"	14"-20" 35-50cm	1252-1614-OGA	☐
40cm	18"-24" 45-60cm	1252-1618-OGA	☐

SPEX CLASSIC GROWTH-ADJUST

WIDTH	LENGTH	CODE	✓
12"	11"-15" 28-38cm	1255-1211-OGA	☐
30cm	14"-18" 35-45cm	1255-1214-OGA	☐
14"	12"-17" 30-43cm	1255-1412-OGA	☐
35cm	16"-21" 40-53cm	1255-1416-OGA	☐
16"	14"-20" 35-50cm	1255-1614-OGA	☐
40cm	18"-24" 45-60cm	1255-1618-OGA	☐

SPEX CLASSIC T SHAPE GROWTH-ADJUST

WIDTH	LENGTH	CODE	✓
12"	11"-15" 28-38cm	1255-1211-TGA	☐
30cm	14"-18" 35-45cm	1255-1214-TGA	☐
14"	12"-17" 30-43cm	1255-1412-TGA	☐
35cm	16"-21" 40-53cm	1255-1416-TGA	☐
16"	14"-20" 35-50cm	1255-1614-TGA	☐
40cm	18"-24" 45-60cm	1255-1618-TGA	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT HARDWARE

Single Mount 1259-2100-100	Dual Mount 1259-2200-100	Quad Mount 1259-2400-100	Drop Mount ☐ Paediatric: 1259-1400-00P ☐ Small: 1259-1400-00S ☐ Medium: 1259-1400-00M ☐ Large: 1259-1400-00L

4. ADD ACCESSORIES

Manta Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
2. Foam System
3. Outer Cover
4. Posterior Pelvic Support & Positioning Wedges

(Note: hardware must be selected)



1. SELECT SIZE


WIDTH	LENGTH	CODE	✓
12"	14" 35cm	1254-1214-000	☐
30cm	16" 40cm	1254-1216-000	☐
13"	14" 35cm	1254-1314-000	☐
33cm	16" 40cm	1254-1316-000	☐
14"	16" 40cm	1254-1416-000	☐
35cm	18" 45cm	1254-1418-000	☐
15"	16" 40cm	1254-1516-000	☐
38cm	18" 45cm	1254-1518-000	☐
16"	18" 45cm	1254-1618-000	☐
40cm	20" 50cm	1254-1620-000	☐
17"	18" 45cm	1254-1718-000	☐
43cm	20" 50cm	1254-1720-000	☐
18"	18" 45cm	1254-1818-000	☐
45cm	20" 50cm	1254-1820-000	☐
19"	18" 45cm	1254-1918-000	☐
48cm	20" 50cm	1254-1920-000	☐
20"	18" 45cm	1254-2018-000	☐
50cm	20" 50cm	1254-2020-000	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT HARDWARE




Single Mount
1259-2100-100

Dual Mount
1259-2200-100

Quad Mount
1259-2400-100

Quad Mount, Heavy Duty, up to 23" wide
1259-2405-100




Quad Mount, Heavy Duty, 23" wide +
1259-2406-100

Direct Mount
1259-2460-100

Drop Mount
Paediatric:
1259-1400-00P
Small:
1259-1400-00S
Medium:
1259-1400-00M
Large:
1259-1400-00L

Quad Extension Brackets
1259-2450-100


4. ADD ACCESSORIES



Additional backrest lumber pad and lateral stability wedges
1209-2713-300

Asymmetric narrow contour wedge
(200mm x 100mm x 40mm)
1209-2722-300

Asymmetric bi-angular contour wedge
(200mm x 100mm x 40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300

Extra Outer Cover (A) Extra Incontinence Cover (B)

Adapta Back Support

1. SELECT BACK SUPPORT




Adapta²
1257-WWHH-200




Adapta³
1257-WWHH-300

2. SELECT FOAM LENGTH - ADAPTA³ ONLY



FRONT

Standard
1257-WWHH-300



FRONT

Kyphotic
1257-WWHH-30K

3. SELECT SIZING

To order replace
WW and HH in code

		Wheelchair Width (inches)					
		12	14	16	18	20	22
Back Support Height (inches)	12						
	13						
	14						
	15						
	16						
	17		+2"				
	18		+2"				
	19		+2"	+2"			
	20						
	21			+2"	+2"	+2"	
	22						
	23				+2"	+2"	

Adapta² Back Support Adapta³ Back Support

ADAPTA³ HAS 2" OF HEIGHT GROWTH

4. SELECT COLOUR

Your colour selection applies to
any extra covers ordered

Black
(default)

Chilli
1209-2491-022

Mango
1209-2491-030

Bumblebee
1209-2491-044

Calypso
1209-2491-028


Royal
1209-2491-077

Ocean
1209-2491-060


Leaf
1209-2491-054

Granite
1209-2491-090


5. SELECT HARDWARE - FOR ADAPTA




Forward Dual Mount
(Default for Adapta²)
1259-2210-100



Dual Mount
(Default for Adapta³)
1259-2205-100




Forward Extended Dual Mount
1259-2215-100



Single Mount
1259-2100-100

6. SELECT UPGRADES & ACCESSORIES



Depth Reduced Axial Swing-Away Lateral Support
1266-2331-100




TSS Lite Contouring Kit
1250-3955-000

Adapta Shoulder Harness Guide Hardware
1400-6693-000

15mm Thick Low Resilient Foam to Front of Backrest
1209-2348-300

SELECT UPGRADES - FOR ADAPTA³ ONLY



Flipback Upgrade for Top Section
1257-WWHH-XFX




Reduced Length Stem - Stylo 160
(mandatory for head support use)
1274-1255-000



Reduced Length Stem - Stylo 260
(mandatory for head support use)
1274-1285-000

7. BUILD CODE

REFER TO THE ADAPTA BACK SUPPORT CODES ABOVE



Comfi Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
2. Foam System
3. Outer Cover
4. Posterior Pelvic Support & Positioning Wedges

(Note: hardware must be selected)



1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
10" 25cm	12" 30cm	1253-1012-00H	☐
	14" 35cm	1253-1014-00H	☐
	16" 40cm	1253-1016-00H	☐
11" 28cm	12" 30cm	1253-1112-00H	☐
	14" 35cm	1253-1114-00H	☐
	16" 40cm	1253-1116-00H	☐
12" 30cm	14" 35cm	1253-1214-00H	☐
	16" 40cm	1253-1216-00H	☐
	18" 45cm	1253-1218-00H	☐
13" 33cm	14" 35cm	1253-1314-00H	☐
	16" 40cm	1253-1316-00H	☐
	18" 45cm	1253-1318-00H	☐
14" 35cm	16" 40cm	1253-1416-00H	☐
	18" 45cm	1253-1418-00H	☐
	20" 50cm	1253-1420-00H	☐
15" 38cm	16" 40cm	1253-1516-00H	☐
	18" 45cm	1253-1518-00H	☐
	20" 50cm	1253-1520-00H	☐
16" 40cm	18" 45cm	1253-1618-00H	☐
	20" 50cm	1253-1620-00H	☐
	23" 58cm	1253-1623-00H	☐
17" 43cm	18" 45cm	1253-1718-00H	☐
	20" 50cm	1253-1720-00H	☐
	23" 58cm	1253-1723-00H	☐
18" 45cm	18" 45cm	1253-1818-00H	☐
	20" 50cm	1253-1820-00H	☐
	23" 58cm	1253-1823-00H	☐
19" 48cm	18" 45cm	1253-1918-00H	☐
	20" 50cm	1253-1920-00H	☐
	23" 58cm	1253-1923-00H	☐
20" 50cm	18" 45cm	1253-2018-00H	☐
	20" 50cm	1253-2020-00H	☐
	23" 58cm	1253-2023-00H	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT HARDWARE

Single Mount 1259-2100-100	Dual Mount 1259-2200-100	Quad Mount 1259-2400-100	Quad Mount, Heavy Duty, up to 23" wide 1259-2405-100
Quad Mount, Heavy Duty, 23" wide + 1259-2406-100	Direct Mount 1259-2460-100	Drop Mount - Paediatric: 1259-1400-00P - Small: 1259-1400-00S - Medium: 1259-1400-00M - Large: 1259-1400-00L	Quad Extension Brackets 1259-2450-100

4. ADD ACCESSORIES

Additional backrest lumber pad and lateral stability wedges 1209-2713-300	Asymmetric narrow contour wedge (200mm x 100mm x 40mm) 1209-2722-300	Asymmetric bi-angular contour wedge (200mm x 100mm x 40mm) 1209-2723-300
Asymmetric contour kit 1209-2720-300	Extra Outer Cover (A) Extra Incontinence Cover (B)	

Spex Classic Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
 2. Foam System
 3. Outer Cover
 4. Spare Pack of foam cubes
 5. Posterior Pelvic Support & Positioning Wedges
- (Note: hardware must be selected)



1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
10" 25cm	12" 30cm	1255-1012-OSS	☐
	14" 35cm	1255-1014-OSS	☐
	16" 40cm	1255-1016-OSS	☐
11" 28cm	12" 30cm	1255-1112-OSS	☐
	14" 35cm	1255-1114-OSS	☐
	16" 40cm	1255-1116-OSS	☐
12" 30cm	14" 35cm	1255-1214-OSS	☐
	16" 40cm	1255-1216-OSS	☐
	18" 45cm	1255-1218-OSS	☐
13" 33cm	14" 35cm	1255-1314-OSS	☐
	16" 40cm	1255-1316-OSS	☐
	18" 45cm	1255-1318-OSS	☐
14" 35cm	16" 40cm	1255-1416-OSS	☐
	18" 45cm	1255-1418-OSS	☐
	20" 50cm	1255-1420-OSS	☐
15" 38cm	16" 40cm	1255-1516-OSS	☐
	18" 45cm	1255-1518-OSS	☐
	20" 50cm	1255-1520-OSS	☐
16" 40cm	18" 45cm	1255-1618-OSS	☐
	20" 50cm	1255-1620-OSS	☐
	23" 58cm	1255-1623-OSS	☐
17" 43cm	18" 45cm	1255-1718-OSS	☐
	20" 50cm	1255-1720-OSS	☐
	23" 58cm	1255-1723-OSS	☐
18" 45cm	18" 45cm	1255-1818-OSS	☐
	20" 50cm	1255-1820-OSS	☐
	23" 58cm	1255-1823-OSS	☐
19" 48cm	18" 45cm	1255-1918-OSS	☐
	20" 50cm	1255-1920-OSS	☐
	23" 58cm	1255-1923-OSS	☐
20" 50cm	18" 45cm	1255-2018-OSS	☐
	20" 50cm	1255-2020-OSS	☐
	23" 58cm	1255-2023-OSS	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black (default)	Chilli 1209-2491-022	Mango 1209-2491-030
Bumblebee 1209-2491-044	Calypso 1209-2491-028	Royal 1209-2491-077
Ocean 1209-2491-060	Leaf 1209-2491-054	Granite 1209-2491-090

3. SELECT HARDWARE

Single Mount 1259-2100-100	Dual Mount 1259-2200-100	Quad Mount 1259-2400-100	Quad Mount, Heavy Duty, up to 23" wide 1259-2405-100
Quad Mount, Heavy Duty, 23" wide + 1259-2406-100	Direct Mount 1259-2460-100	Drop Mount - Paediatric: 1259-1400-00P - Small: 1259-1400-00S - Medium: 1259-1400-00M - Large: 1259-1400-00L	Quad Extension Brackets 1259-2450-100

4. ADD ACCESSORIES

Additional backrest lumber pad and lateral stability wedges 1209-2713-300	Asymmetric narrow contour wedge (200mm x 100mm x 40mm) 1209-2722-300	Asymmetric bi-angular contour wedge (200mm x 100mm x 40mm) 1209-2723-300
Asymmetric contour kit 1209-2720-300	Extra Outer Cover (A) Extra Incontinence Cover (B)	Extra Positioning Kit (pack 10 cubes) 1255-XXXX-300

Spex Classic Back Support with WheelAir

PACKAGE INCLUDES:

1. Spex Classic backrest squab with unique contouring cube system, outer cover and rectangular backrest shell.
2. Supracore Front Pocket to allow for greater airflow.
3. WheelAir Ventilation System, includes: Fanbox, Remote Controller, Nozzle Tube Magnets, Nozzle Channels 90°, Flexible Tubes, 90° elbows, mounting bracket and Air Channels.

(Note: hardware must be selected)



1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
14" 35cm	16" 40cm	1255-1416-OSS-WAIR	☐
	18" 45cm	1255-1418-OSS-WAIR	☐
	20" 50cm	1255-1420-OSS-WAIR	☐
15" 38cm	16" 40cm	1255-1516-OSS-WAIR	☐
	18" 45cm	1255-1518-OSS-WAIR	☐
	20" 50cm	1255-1520-OSS-WAIR	☐
16" 40cm	18" 45cm	1255-1618-OSS-WAIR	☐
	20" 50cm	1255-1620-OSS-WAIR	☐
	23" 58cm	1255-1623-OSS-WAIR	☐
17" 43cm	18" 45cm	1255-1718-OSS-WAIR	☐
	20" 50cm	1255-1720-OSS-WAIR	☐
	23" 58cm	1255-1723-OSS-WAIR	☐
18" 45cm	18" 45cm	1255-1818-OSS-WAIR	☐
	20" 50cm	1255-1820-OSS-WAIR	☐
	23" 58cm	1255-1820-OSS-WAIR	☐
19" 48cm	18" 45cm	1255-1918-OSS-WAIR	☐
	20" 50cm	1255-1920-OSS-WAIR	☐
	23" 58cm	1255-1923-OSS-WAIR	☐
20" 50cm	18" 45cm	1255-2018-OSS-WAIR	☐
	20" 50cm	1255-2020-OSS-WAIR	☐
	23" 58cm	1255-2023-OSS-WAIR	☐
21" 53cm	18" 45cm	1255-2118-OSS-WAIR	☐
	20" 50cm	1255-2120-OSS-WAIR	☐
	23" 58cm	1255-2123-OSS-WAIR	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black
(default)

Chilli
1209-2491-022

Mango
1209-2491-030

Bumblebee
1209-2491-044

Calypso
1209-2491-028


Royal
1209-2491-077

Ocean
1209-2491-060


Leaf
1209-2491-054

Granite
1209-2491-090


3. SELECT HARDWARE




Single Mount
1259-2100-100




Dual Mount
1259-2200-100




Quad Mount
1259-2400-100




Quad Mount,
Heavy Duty,
up to 23" wide
1259-2405-100




Quad Mount,
Heavy Duty, 23"
wide +
1259-2406-100



Direct Mount
1259-2460-100




Drop Mount
☐ Paediatric:
1259-1400-00P
☐ Small:
1259-1400-00S
☐ Medium:
1259-1400-00M
☐ Large:
1259-1400-00L



Quad Extension
Brackets
1259-2450-100


4. ADD ACCESSORIES




Additional backrest
lumber pad and lateral
stability wedges
1209-2713-300



Asymmetric narrow
contour wedge
(200mm x 100mm x
40mm)
1209-2722-300




Asymmetric bi-angular
contour wedge
(200mm x 100mm x
40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300



Extra Outer Cover (A)
Extra Incontinence
Cover (B)



Extra Positioning Kit
(pack 10 cubes)
1255-XXXX-300

Spex T-Shape Back Support

PACKAGE INCLUDES:

1. T-Shape Aluminium Shell
 2. Foam System
 3. Outer Cover
 4. Spare Pack of foam cubes
 5. Posterior Pelvic Support & Positioning Wedges
- (Note: hardware must be selected)



1. SELECT SIZE

WIDTH	LENGTH	CODE	✓
10" 25cm	12" 30cm	1255-1012-OTS	☐
	14" 35cm	1255-1014-OTS	☐
	16" 40cm	1255-1016-OTS	☐
11" 28cm	12" 30cm	1255-1112-OTS	☐
	14" 35cm	1255-1114-OTS	☐
	16" 40cm	1255-1116-OTS	☐
12" 30cm	14" 35cm	1255-1214-OTS	☐
	16" 40cm	1255-1216-OTS	☐
	18" 45cm	1255-1218-OTS	☐
13" 33cm	14" 35cm	1255-1314-OTS	☐
	16" 40cm	1255-1316-OTS	☐
	18" 45cm	1255-1318-OTS	☐
14" 35cm	16" 40cm	1255-1416-OTS	☐
	18" 45cm	1255-1418-OTS	☐
	20" 50cm	1255-1420-OTS	☐
15" 38cm	16" 40cm	1255-1516-OTS	☐
	18" 45cm	1255-1518-OTS	☐
	20" 50cm	1255-1520-OTS	☐
16" 40cm	18" 45cm	1255-1618-OTS	☐
	20" 50cm	1255-1620-OTS	☐
	23" 58cm	1255-1623-OTS	☐
17" 43cm	18" 45cm	1255-1718-OTS	☐
	20" 50cm	1255-1720-OTS	☐
	23" 58cm	1255-1723-OTS	☐
18" 45cm	18" 45cm	1255-1818-OTS	☐
	20" 50cm	1255-1820-OTS	☐
	23" 58cm	1255-1823-OTS	☐
19" 48cm	18" 45cm	1255-1918-OTS	☐
	20" 50cm	1255-1920-OTS	☐
	23" 58cm	1255-1923-OTS	☐
20" 50cm	18" 45cm	1255-2018-OTS	☐
	20" 50cm	1255-2020-OTS	☐
	23" 58cm	1255-2023-OTS	☐

2. SELECT COLOUR

Your colour selection applies to any extra covers ordered

Black
(default)

Chilli
1209-2491-022

Mango
1209-2491-030

Bumblebee
1209-2491-044

Calypso
1209-2491-028


Royal
1209-2491-077

Ocean
1209-2491-060


Leaf
1209-2491-054

Granite
1209-2491-090


3. SELECT HARDWARE




Single Mount
1259-2100-100




Dual Mount
1259-2200-100




Quad Mount
1259-2400-100




Quad Mount,
Heavy Duty,
up to 23" wide
1259-2405-100




Quad Mount,
Heavy Duty, 23"
wide +
1259-2406-100



Direct Mount
1259-2460-100




Drop Mount
☐ Paediatric:
1259-1400-00P
☐ Small:
1259-1400-00S
☐ Medium:
1259-1400-00M
☐ Large:
1259-1400-00L



Quad Extension
Brackets
1259-2450-100


4. ADD ACCESSORIES




Additional backrest
lumber pad and lateral
stability wedges
1209-2713-300



Asymmetric narrow
contour wedge
(200mm x 100mm x
40mm)
1209-2722-300



Asymmetric bi-angular
contour wedge
(200mm x 100mm x
40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300



Extra Outer Cover (A)
Extra Incontinence
Cover (B)



Extra Positioning Kit
(pack 10 cubes)
1255-XXXX-300

Spex SuperShape Back Support

PACKAGE INCLUDES:

1. Aluminium Shell
2. Foam System
3. Outer Cover
4. Spare Pack of foam cubes including firm EVA Posterior Pelvic Support & Positioning Wedges


(Note: hardware must be selected)



1256-WWDD-QSS


1. SELECT SIZING

To order replace **WW** and **HH** in code




2. SELECT COLOUR


Your colour selection applies to any extra covers ordered




3. SELECT HARDWARE




Single Mount
1259-2100-100




Dual Mount
1259-2200-100




Quad Mount
1259-2400-100




Quad Mount, Heavy Duty, up to 23" wide
1259-2405-100




Quad Mount, Heavy Duty, 23" wide +
1259-2406-100



Direct Mount
1259-2460-100




Drop Mount
☐ **Paediatric:**
1259-1400-00P
☐ **Small:**
1259-1400-00S
☐ **Medium:**
1259-1400-00M
☐ **Large:**
1259-1400-00L




Quad Extension Brackets
1259-2450-100


4. ADD ACCESSORIES




Additional backrest lumber pad and lateral stability wedges
1209-2713-300



Asymmetric narrow contour wedge
(200mm x 100mm x 40mm)
1209-2722-300



Asymmetric bi-angular contour wedge
(200mm x 100mm x 40mm)
1209-2723-300




Asymmetric contour kit
1209-2720-300



Extra Outer Cover (A)
Extra Incontinence Cover (B)




Extra Positioning Kit
(pack 10 cubes)
1255-XXXX-300



T-shape modification
1209-2933-000

5. BUILD CODE

REFER TO THE SUPERSHAPE BACK SUPPORT CODES ABOVE



XLella Back Support

1. SELECT BACK SUPPORT



XLella, Vigour Padding
1258-WWHH-00H



XLella, Classic Padding
1258-WWHH-OSS



XLella2, Biangular
1258-WWHH-200

2. SELECT SIZING

To order replace
WW and HH in code


		Wheelchair Width (inches)					
		20	22	24	26	28	30
Back Support Height (inches)	12						
	13						
	14						
	15						
	16						
	17						
	18						
	19						
	20						
	21						
	22						
	23						
	24						
	25						

XLella Vigour or Classic


XLella* Biangular

3. SELECT COLOUR


Your colour selection applies to
any extra covers ordered



Black
(default)




Chilli
1209-2491-022




Mango
1209-2491-030



Bumblebee
1209-2491-044



Calypso
1209-2491-028



Royal
1209-2491-077



Ocean
1209-2491-060




Leaf
1209-2491-054




Granite
1209-2491-090

4. SELECT HARDWARE - FOR XLELLA




Fixed Bariatric Hardware (default)
1259-2500-100




Quick-Release Bariatric Hardware
1259-2505-100

5. SELECT UPGRADES & ACCESSORIES




Cover Extension
☐ 22"
1258-2200-500
☐ 24"
1258-2400-500
☐ 26"
1258-2600-500
☐ 28"
1258-2800-500




Reinforcement Plate
for Lateral Supports
(forceful side leaning
or users over 100kg)
1266-2450-100




Additional backrest
lumber pad and lateral
stability wedges
1209-2713-300



Asymmetric narrow
contour wedge
(200mm x 100mm x
40mm)
1209-2722-300



Asymmetric bi-angular
contour wedge
(200mm x 100mm x
40mm)
1209-2723-300



Asymmetric contour kit
1209-2720-300




Extra Outer Cover (A)
Extra Incontinence
Cover (B)



Heavy-Duty kit for
Stylo Head Support
1274-0000-HDK

6. BUILD CODE

REFER TO THE XLELLA BACK SUPPORT CODES ABOVE



Back Support Custom Sizes

1. SELECT BACK SUPPORT

(cover & back support can be selected simultaneously)

UP TO 20"

	COVER	BACK SUPPORT
Custom Vigour	☐	☐
Vigour	☐	☐
Custom Manta	☐	☐
Adapta	☐	☐
Custom Comfi	☐	☐
Custom Spex Classic	☐	☐
Custom Spex T-Shape	☐	☐
Spex Growth-Adjust	☐	☐
Spex Supershapes	☐	☐
Xlella Vigour	☐	☐
Xlella Classic	☐	☐
Xlella ²	☐	☐

UP TO 25"

SIZES ABOVE 21", REFER XLELLA

	COVER	BACK SUPPORT
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐
☐	☐	☐

SIZES ABOVE 21", REFER XLELLA

IMPORTANT

When ordering custom size back supports, lead time can be extended by up to 12 working days.

2. ENTER MEASUREMENTS

1

WHEELCHAIR WIDTH

XX
 1 = ____


2

SHELL HEIGHT

XX
 2 = ____

3

For standard mounting using Single/Dual/Quad mounts the distance between edge of shell and inside of cane needs to be a minimum of 20mm and maximum of 45mm (each side).



3. SPECIFY PADDING

Default Padding
(wheelchair width -1", for Xlella backs, wheelchair width +1")
Default

Slim Fit Padding
(same width as shell)
1249-2339-300

Custom Fit Padding
1249-2349-30

Default

1249-2339-300

1249-2349-300

LENGTH
 L = ____

WIDTH
 W = ____

Back Support Accessories

TESSELLATED SHAPE SYSTEM KIT

WIDTH	HEIGHT	ITEM CODE	✓
10" 25cm	12" 30cm	1250-1012-000	☐
	14" 35cm	1250-1014-000	☐
	16" 40cm	1250-1016-000	☐
12" 30cm	14" 35cm	1250-1214-000	☐
	16" 40cm	1250-1216-000	☐
	18" 45cm	1250-1218-000	☐
14" 35cm	16" 40cm	1250-1416-000	☐
	18" 45cm	1250-1418-000	☐
	20" 50cm	1250-1420-000	☐
16" 40cm	18" 45cm	1250-1618-000	☐
	20" 50cm	1250-1620-000	☐
	23" 58cm	1250-1623-000	☐
18" 45cm	18" 45cm	1250-1818-000	☐
	20" 50cm	1250-1820-000	☐
	23" 58cm	1250-1823-000	☐
20" 50cm	18" 45cm	1250-2018-000	☐
	20" 50cm	1250-2020-000	☐
	23" 58cm	1250-2023-000	☐
Up to 20" (50cm)	XX XX	1250-2929-000	☐
Up to 25" (63cm)	XX XX	1250-2939-000	☐
TSS to suit back ordered on Script (inc. shape mods)		1250-3949-000	☐

STANDARD TESSELLATED SHAPE SYSTEM KIT




ELEMENTS BUILDING KIT ONLY

TSS Elements building kit only
(no backboard).
Use when mounting shape
elements directly to back shell.
1250-3959-000




EXTRA POSITIONING WEDGES




**Additional backrest
lumber pad and lateral
stability wedges**
1209-2713-300



**Asymmetric narrow
contour wedge**
(20cm x 10cm x 4cm)
1209-2722-300




**Asymmetric bi-angular
contour wedge**
(20cm x 10cm x 4cm)
1209-2723-300




Asymmetric contour kit
1209-2720-300


EXTRA HARDWARE



Single Mount
for second
wheelchair
1259-2120-100




Dual Mount
for second
wheelchair
1259-2220-100



Quad Mount
for second
wheelchair
1259-2420-100

ALTERNATIVE MOUNTING OPTIONS



**Spex Power
Recline Adapter
Mount**
1259-2810-200




**Spex Standing
Powerchair
Adapter Mount**
1259-2610-200

SUPRACOR® HONEYCOMB

**Envelope Cover only for
Supracor® insert**
1208-2951-300

**Envelope Cover with
Supracor insert**
1208-2952-300




Supracor Honeycomb covers add
13mm to cushion height and provide
heat & sweat reduction.

Supracor Honeycomb covers must be
ordered with back support or spare cover.

Lateral Trunk Support


1. SELECT HARDWARE

(Laterals are scripted as each, not in pairs)




Fixed Hardware

- ☐ Fixed Bracket: 1266-2101-100
- ☐ Paediatric Fixed Bracket: 1266-2111-100




Swing Away Hardware

- ☐ Straight: 1266-2203-100
- ☐ Extended Straight: 1266-2202-100




Offset Swing Away

- ☐ Offset: 1266-2101-100
- ☐ Extended Offset: 1266-2204-100



Axial Swing Away Hardware

- ☐ Axial: 1266-2301-100
- ☐ Depth-Reduced Axial: 1266-2331-100
- ☐ Paediatric Axial
- ☐ Lateral: 1266-2613-100-[L/R]
- ☐ Paediatric Swing Away: 1266-2623-400
- ☐ Paediatric Swing Away Cane Mount: 1266-2623-300-L/R
- ☐ Paediatric Axial Forward Mount: 1266-2623-200-L/R



Quick-Release

- ☐ Back Mount Fixed Angle: 1266-2601-100-L/R
- ☐ Cane Mount Fixed Angle: 1266-2601-200-L/R
- ☐ Back Mount, Swing-Away, Angle Adjustment: 1266-2621-100-L/R
- ☐ Cane Mount, Swing-Away, Angle Adjustment: 1266-2621-200-L/R

2. SELECT PAD & SIZE

TEARDROP



SIZES AVAILABLE

0 (140mm x 80mm)

-XX51-

-XX52-

SELECT

Left

Right

NARROW STRAIGHT



SIZES AVAILABLE

0 (85mm x 115mm)

-XX10-

1 (95mm x 145mm)

-XX11-

2 (95mm x 180mm)

-XX12-

3 (115mm x 210mm)

-XX13-

4 (115mm x 240mm)

-XX14-

SELECT

L R

L R

L R

L R

L R

L R

L R

L R

NARROW CURVED



SIZES AVAILABLE

1 (95mm x 145mm)

-XX21-

2 (95mm x 180mm)

-XX22-

3 (115mm x 210mm)

-XX23-

4 (115mm x 240mm)

-XX24-

SELECT

L R

L R

L R

L R

L R

L R

L R

BROAD



SIZES AVAILABLE

1 (105mm x 130mm)

-XX31-

2 (125mm x 140mm)

-XX32-

3 (140mm x 150mm)

-XX33-

4 (160mm x 155mm)

-XX34-

5 (195mm x 165mm)

-XX35-

6 (230mm x 175mm)

-XX36-

SELECT

L R

L R

L R

L R

L R

L R

L R

L R

L R

L R

3. SELECT COVER



Standard (-400)



Comfi (-800)

Incompatible with Quick-Release

4. CUSTOMISE



Left Cane Mount

(attachment directly to wheelchair canes)
1266-2411-100




Right Cane Mount

(attachment directly to wheelchair canes)
1266-2412-100



Reinforcement plate

(forceful side leaning or clients over 100kg)
1266-2450-100




Deep Contour Bracket

(extends depth by 30mm)
1266-2420-100



Deep Contour Extended Bracket Kit

(extends depth up to 115mm)
1266-2430-100



Extender Kit for Paediatric Swing-Away Lateral

(extends depth up to 40mm)
1266-2440-100

5. BUILD CODE FOR LATERAL PAD


REFER TO THE LATERAL TRUNK SUPPORT CODES ABOVE

1	2	6	6	-	3	0						
						ADD PAD CODES			ADD COVER CODES			

Axial Biangular Lateral Trunk Support

1. SELECT HARDWARE


(Laterals are scripted as each, not in pairs)



L R
1266-2501-100

2. SELECT BIANGULAR PAD SIZE

AXIAL PAD



SIZES AVAILABLE

4 (165mm x 175mm)

-XX44-

6 (215mm x 195mm)

-XX46-

8 (265mm x 195mm)

-XX48-

SELECT

L

R

L

R

L

R


5. CUSTOMISE




Left Cane Mount
(attachment directly to wheelchair canes)
1266-2411-100



Right Cane Mount
(attachment directly to wheelchair canes)
1266-2412-100



Reinforcement plate
(forceful side leaning or clients over 100kg)
1266-2450-100




Deep Contour Bracket
(extends depth by 30mm)
1266-2420-100



Deep Contour Extended Bracket Kit
(extends depth up to 115mm)
1266-2430-100

3. SELECT PADDING



12mm Pressure-Reduction
Overall Pad Thickness 30mm
-X00



24mm Pressure-Reduction
Overall Pad Thickness 45mm
-X10

4. SELECT COVER




Standard (-4XX)



Comfi (-8XX)


6. BUILD CODE FOR LATERAL PAD

REFER TO THE LATERAL TRUNK SUPPORT CODES ABOVE




Arm Supports

1. SELECT ARM SUPPORT HARDWARE




- ☐ **Fixed**
-01XX-
- ☐ **Fixed Paediatric**
-02XX-




- ☐ **Angle Adjust**
-03XX-

2. SELECT ARM SUPPORT MOUNTING




- ☐ **Seat Rail**
-2XX-




- ☐ **Seat Pan**
-1XX-

3. SELECT ARM SUPPORT MECHANISM



- ☐ **Quick Release**
-X1X-



- ☐ **Swing-Down**
-X2X-

4. SELECT COMFI ARM SUPPORT PADS



SIZES AVAILABLE

1 NARROW	1156-5811-410	40mm x 250mm	☐
2 NARROW	1156-5812-410	50mm x 300mm	☐
3 NARROW	1156-5813-410	60mm x 350mm	☐
1 WIDE	1156-5821-410	70mm x 250mm	☐
2 WIDE	1156-5822-410	80mm x 300mm	☐
3 WIDE	1156-5823-410	90mm x 350mm	☐

CUSTOM


1156-5829-410

XX
A = ____

XX
B = ____


5. CUSTOMISE

(tick any customisations that apply)



Controller Cutout
1300-5839-400

XX
C = ____ D = ____



☐ **Custom Thickness**
1300-5842-400


XX
E = ____



☐ **½" Gel-Foam overlay**
(custom fabric on arm support pad)
1300-5847-400




Spextex waterproof cover
(knitted fabric surface)
1300-5849-219




Rear elbow pad & bracket
1300-5844-400

XX
F = ____ G = ____



Side elbow pad
1300-5845-400

XX
I = ____ J = ____



Arm troughing
1300-5835-400

XX
A = ____ B = ____
H = ____

6. BUILD CODE

REFER TO THE ARM SUPPORT CODES ABOVE

1 1 5 6 - 0 0 - 0


ADD HARDWARE CODES

ADD MOUNTING CODES


ADD MECHANISM CODES

Thigh Supports

1. SELECT THIGH SUPPORT HARDWARE




☐ Fixed
-0002-




☐ Axial
-0001-

2. SELECT THIGH SUPPORT MOUNTING



☐ Seat Rail
-2XX



☐ Seat Pan
-1XX

3. SELECT THIGH SUPPORT MECHANISM



☐ Quick Release
-X1X



☐ Swing-Down
-X2X

4) a. SELECT THIGH SUPPORT PADS

COMFY THIGH SUPPORT PADS



SIZES AVAILABLE

10" (25cm)	1156-1010-400	90 x 265 x 35mm	☐
12" (30cm)	1156-1012-400	95 x 310 x 35mm	☐
14" (35cm)	1156-1014-400	130 x 360 x 35mm	☐
16" (40cm)	1156-1016-400	135 x 410 x 35mm	☐
18" (45cm)	1156-1018-400	145 x 460 x 35mm	☐
20" (50cm)	1156-1020-400	160 x 510 x 35mm	☐

CUSTOM

1156-1029-400 ☐ A = XX B = XX E = XX

4) b. SELECT THIGH SUPPORT PADS

WONDERSEAT THIGH SUPPORT PADS - NON-CUSTOMISABLE



☐ 5" Deep, 11.5 x 8.5cm
1156-1105-400



☐ Snug - 8" Deep, 19 x 8.5cm
1156-1208-400-[L/R]



☐ 8" Deep, 19 x 8.5cm
1156-1108-400



☐ 10" Deep, 26.5 x 9cm
1156-1010-400


5. CUSTOMISE

(tick any customisations that apply)
(N/A Wonderseat Thigh Support Pads)

SELECT	ITEM	IMAGE
☐ L ☐ R	Custom Cut-out 1153-1139-400 C = XX D = XX	
☐	Custom Thickness 1153-1145-400 E = XX	
☐	1/2" Gel-Foam overlay (custom fabric on thigh support pad) 1153-1147-400	
☐	Spextex waterproof cover (wipe-down surface) 1153-1149-209	


6. BUILD CODE

REFER TO THE THIGH SUPPORT CODES ABOVE




Combo Arm/Thigh Hardware

1. SELECT SEAT PAN MOUNTING



Seat Rail


☐ -2XX



Seat Pan


☐ -1XX

2. SELECT MECHANISM



Quick Release

☐ -X1X



Swing-Down

☐ -X2X

3. SELECT CONFIGURATION

SIZES AVAILABLE

ANGLE-ADJUST ARM + AXIAL THIGH XXXX-0301-XXX ☐

ANGLE-ADJUST ARM + FIXED THIGH XXXX-0302-XXX ☐

FIXED ARM + AXIAL THIGH XXXX-0101-XXX ☐

PAEDIATRIC FIXED ARM + AXIAL THIGH XXXX-0201-XXX ☐

FIXED ARM + FIXED THIGH XXXX-0102-XXX ☐

PAEDIATRIC FIXED ARM + FIXED THIGH XXXX-0202-XXX ☐

4. SELECT COMFI ARM SUPPORT PADS



SIZES AVAILABLE

1 NARROW 1156-5811-410 40mm x 250mm ☐

2 NARROW 1156-5812-410 50mm x 300mm ☐

3 NARROW 1156-5813-410 60mm x 350mm ☐

1 WIDE 1156-5821-410 70mm x 250mm ☐

2 WIDE 1156-5822-410 80mm x 300mm ☐

3 WIDE 1156-5823-410 90mm x 350mm ☐

CUSTOM

1156-5829-410

A = ____

B = ____

5. SELECT THIGH SUPPORT PADS

COMFI THIGH PADS



SIZES AVAILABLE

10" (25cm) 1156-1010-400 90 x 265 x 35mm ☐

12" (30cm) 1156-1012-400 95 x 310 x 35mm ☐

14" (35cm) 1156-1014-400 130 x 360 x 35mm ☐

16" (40cm) 1156-1016-400 135 x 410 x 35mm ☐

18" (45cm) 1156-1018-400 145 x 460 x 35mm ☐

20" (50cm) 1156-1020-400 160 x 510 x 35mm ☐

CUSTOM

1156-1029-400

A = ____

B = ____

E = ____

WONDERSEAT THIGH PADS



☐ 5" Deep, 11.5 x 8.5cm
1156-1105-400



☐ Snug - 8" Deep, 19 x 8.5cm
1156-1208-400-[L/R]



☐ 8" Deep, 19 x 8.5cm
1156-1108-400



☐ 10" Deep, 26.5 x 9cm
1156-1010-400

4. BUILD CODE

REFER TO THE ARM/THIGH SUPPORT CODES ABOVE

1 1 5 6 - - - - - 0

ADD
CONFIGURATION
CODES

ADD MOUNTING
CODES

ADD MECHANISM
CODES

Medial Thigh Support


1. SELECT HARDWARE



☐ Standard
-150X-



☐ Axial
-160X-



2. SELECT MEDIAL THIGH SUPPORT SIZE

STANDARD PAD

☐	SIZE 1	XXXX-1501-XXX	A = 60mm	B = 145mm	C = 95mm	D = 40mm	E = 40mm	F = 65mm
☐	SIZE 2	XXXX-1502-XXX	A = 60mm	B = 155mm	C = 100mm	D = 60mm	E = 50mm	F = 65mm
☐	SIZE 3	XXXX-1503-XXX	A = 60mm	B = 185mm	C = 105mm	D = 60mm	E = 60mm	F = 75mm
☐	CUSTOM	XXXX-1509-XXX	XX A = ____	XX B = ____	XX C = ____	XX D = ____	XX E = ____	XX F = ____

AXIAL ADJUSTABLE PAD

☐	SIZE 1	XXXX-1601-XXX	A = 70-140mm	B = 145mm	C = 95mm	D = 50-120mm	E = 40mm	F = 65mm
☐	SIZE 2	XXXX-1602-XXX	A = 70-140mm	B = 155mm	C = 100mm	D = 60-140mm	E = 50mm	F = 65mm
☐	SIZE 3	XXXX-1603-XXX	A = 70-160mm	B = 185mm	C = 105mm	D = 60-140mm	E = 60mm	F = 75mm
☐	CUSTOM	XXXX-1609-XXX	XX A = ____	XX B = ____	XX C = ____	XX D = ____	XX E = ____	XX F = ____

3. CUSTOMISE

(tick any customisations that apply)

SELECT ITEM

IMAGE

☐	½" Pressure Relieving foam overlay (Custom fabric on medial thigh support cover) 1160-1576-000	
☐	½" Gel-Foam overlay (Custom fabric on medial thigh support cover) 1153-1147-400	
☐	Spextex waterproof cover (Wipe-down surface) 1153-1149-209	
☐	Spextex waterproof cover (Knitted fabric surface) 1153-1149-209	

4. BUILD CODE

REFER TO THE MEDIAL THIGH SUPPORT CODES ABOVE

1 1 6 2 - - - - - 0 0 0

ADD HARDWARE
CODES

ADD SIZE CODES

Stylo Head Support

1. SELECT HARDWARE



Stylo 130
1274-1232-000

Stylo 160
1274-1252-000

Stylo 260
1274-1282-000




Stylo Flip-Back 160
1274-2252-000

Stylo Flip-Back 260
1274-2282-000

2. SELECT PAD & SIZE

(Width x Height x Depth)

COMFORT




250 1275-0031-000

260mm x 125mm
x 90mm

350 1275-0032-000

340mm x 165mm
x 120mm

CIRCLE




150 1275-0011-000

175mm x 145mm
x 55mm

250 1275-0012-000

220mm x 185mm
x 70mm

CONTOUR




250 1275-0021-000

280mm x 165mm
x 140mm

350 1275-0022-000

350mm x 210mm
x 170mm

ADJUSTABLE LATERAL




150 1275-0061-000

175mm x 145mm
x 55mm

250 1275-0062-000

220mm x 185mm
x 70 mm

STANDARD AND EXTENDED LATERAL



Standard Paediatric 1275-0041-000

220mm x 150mm
x 145mm

Standard 1275-0043-000

250mm x 170mm
x 180mm

Extended Paediatric 1275-0042-000


230mm x 180mm
x 160mm

Extended 1275-0044-000


260mm x 190mm
x 180mm

3. SELECT COLOUR

Your colour selection applies to any extra covers ordered



4. CUSTOMISE



1/2" soft pressure-relieving foam overlay
1271-0019-000

Spextex waterproof cover - wipe-down surface
1271-0039-209

Spextex waterproof cover - knitted fabric surface
1271-0039-219




Neoprene cover
1271-0039-291

Stylo Heavy-Duty Compatibility Kit
(Allows attachment to other backs and reinforces connection, reducing shell flex and ball slip)
1274-0000-HDK

Stylo Slim Mount
(For use where space is limited, e.g. power wheelchairs)
1274-4140-200

CUSTOMISE (ADJUSTABLE LATERAL HEADREST ONLY)



Custom wing headrest pads
1270-4149-300


Swing-away lateral mech upgrade
1272-4250-100

XX
L = ____ H = ____


Harnesses

1. SELECT HARNESS

Tick option that applies




H-HARNESS
-664X-



RETRACTOR
-662X-



VEST
-661X-



CENTREPOINT
-666X-

2. SELECT SIZE




4. SELECT ACCESSORIES



Shoulder Strap Guides
(Attachment Hardware Included)
1400-6692-000


3. CUSTOMISE

(tick any customisations that apply)




Removable neoprene sleeves
1409-6810-000


XX
L = ____



Buckles
(Fitted to lower shoulder harness straps)
1409-6820-000



Custom shape alteration
(More information required with order)
1409-6834-000




Extra-long webbing straps
1409-6815-000


XX
L = ____



Change buckle to centre-release style
(H-Harness & Centrepoint only)
1409-6828-000



Change buckle to anti-escape style
1409-6826-000



Convert shoulder harness to rear-pull style
1409-6841-000

5. BUILD CODE

REFER TO THE HARNESS CODES ABOVE

1 4 0 9 - - - - - 0 1 7

ADD HARNESS
CODES

ADD
SIZE CODES

IMPORTANT

- S, M, L & XL available only for retractor harness.
- XS, S, M & L available only for vest harness.

Chest Straps

SELECT SIZE



- ☐ **Small**
(35cm W x 7.5cm H)
1409-6821-017
- ☐ **Medium**
(40cm W x 8cm H)
1409-6822-017
- ☐ **Large**
(55cm W x 10cm H)
1409-6823-017

Hip Belts

1. SELECT HIP BELT

Tick option that applies



✓	SIZE	2 POINT CENTRE PULL	2 POINT DUAL CENTRE PULL	4 POINT SIDE PULL	4 POINT CENTRE PULL	4 POINT DUAL CENTRE PULL
☐	X-SMALL	XXX5-6210-XXX	XXX5-6220-XXX	XXX6-6400-XXX	XXX6-6410-XXX	XXX6-6420-027
☐	SMALL	XXX5-6211-XXX	XXX5-6221-XXX	XXX6-6401-XXX	XXX6-6411-XXX	XXX6-6421-027
☐	MEDIUM	XXX5-6212-XXX	XXX5-6222-XXX	XXX6-6402-XXX	XXX6-6412-XXX	XXX6-6422-027
☐	LARGE	XXX5-6213-XXX	XXX5-6223-XXX	XXX6-6403-XXX	XXX6-6413-XXX	XXX6-6423-027

2. SELECT BUCKLE

Tick option that applies



SIDE RELEASE
-027



CENTRE RELEASE
-028



ANTI-ESCAPE
-026



PEN-POPPER
-025

3. SELECT OPTIONAL ADD-ON

Pressure reducing gel-foam slip on sleeves for all hip belts.
Retrofits to hip belt pads.

☐ 1409-6812-000



4. BUILD CODE

REFER TO THE HIP BELT & BUCKLE CODES ABOVE

1 4 0 - - - - -

ADD HIP BELT CODES

ADD BUCKLE CODES

IMPORTANT

- The 4 point dual centre-pull can only be used with the side release buckle.
- The 2 point dual centre-pull can only be used with the side release buckle.
- The anti-escape buckle is not available for x-small hip belts.

Limb Supports

PADDED FOREARM STRAPS

SIZE	SET OF FOUR	✓	ONE ONLY	QTY
SMALL	1408-6511-017	○	1408-6521-017	
MEDIUM	1408-6512-017	○	1408-6522-017	
LARGE	1408-6513-017	○	1408-6523-017	

PADDED FOREARM STRAPS



HEEL LOOPS

SIZE	LENGTH	PAIR	✓
SMALL	230mm	1403-6181-017	○
MEDIUM	255mm	1403-6182-017	○
LARGE	275mm	1403-6183-017	○

HEEL LOOPS



FOOTFAST STABILISERS

SIZE	CIRCUMFERENCE	STANDARD PAIR	✓	WITH REAR RETAINER STRAP	✓
SMALL	15-23cm (Inner)	1401-6161-017	○	1401-6165-017	○
MEDIUM	22-30cm (Inner)	1401-6162-017	○	1401-6166-017	○
LARGE	26-38cm (Inner)	1401-6163-017	○	1401-6167-017	○

FOOTFAST STABILISERS



REAR RETAINER STRAP



COMFORT PADDED & NEOPRENE CALF STRAPS

W/C WIDTH	HEIGHT	COMFORT PADDED CODE	✓	NEOPRENE CODE	✓
10"	5cm	1402-6170-017	○	1402-6070-017	○
12"	5cm	1402-6171-017	○	1402-6071-017	○
14"	7cm	1402-6172-017	○	1402-6072-017	○
16"	9cm	1402-6173-017	○	1402-6073-017	○
17"	9cm	1402-6174-017	○	1402-6074-017	○
18"	11cm	1402-6175-017	○	1402-6075-017	○
20"	13cm	1402-6176-017	○	1402-6076-017	○
XX	XX	1402-6177-017	○	1402-6077-017	○

PADDED CALF STRAPS



NEOPRENE CALF STRAPS



CALF PANELS

W/C WIDTH	WIDTH	HEIGHT	CODE	✓
14-15"	14"	8"	1402-6314-017	○
16-17"	16"	8"	1402-6316-017	○
18-19"	18"	10"	1402-6318-017	○
20"	20"	10"	1402-6320-017	○
Custom	XX	XX	1402-6329-017	○

CALF PANELS



SIZE	WIDTH	DEPTH	CUT-OUT DEPTH	CODE	✓
SMALL	50cm	8"	15cm	1311-7441-209	☐
MEDIUM	60cm	8"	15cm	1311-7442-209	☐
LARGE	68cm	10"	20cm	1311-7443-209	☐



Spex Notes

[illegible][illegible]



Manufactured by Spex Ltd

solutions@spexseating.com

spexseating.com

