

Starlight-NXT




CAR SEAT

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Seat Width	
B Seat Depth	
C Shoulder Strap Height	
D Foot Support Height	
E Chest Width	
F Axilla Height	
User Weight	



Basic Packages

Includes:			
- Adjustable head support - Back support infill padding - Slim lateral supports 20° tilt - Base board accessory - Clamp Bow Stabiliser - Swivel base system - (Left turning) to enable turning of the car seat towards the doorway - Grey/Black upholstery	Left Turning 4141-0011-090		

Specifications

	Seat Width (pelvis)	Seat Width (knee)	Seat Depth	Back Height	Crown Height	Lateral Height Adjustment	Lateral Width Adjustment	Product Weight	Weight Capacity
Standard	30cm	35cm	26cm	32–42cm	51–62cm	17–27cm	14–19cm	8.5kg	9–36kg
XL	30cm	35cm	32cm	44–53cm	67–76cm	25–38cm	14–19cm	8.5kg	9–36kg

XL Upgrade Kit

Includes:			
- Extended Harness - Longer back support padding - Additional mounting plates (Harness upgrade only possible when child is over 15kg)	4141-0140-000		

Accessories

ISOfix Adapter

4170-6911-100



Foot Support

Standard
4141-5240-000



Abductor Block

Removable
4141-1510-000



Tray Table

Padded
4141-5530-000



Seat Protection Mat

4141-0190-000



Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other markings on the paper.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Starlight-NXT_20250716_NZ