

**ORDER FORM - Medifab**

3- Shipping information

Ship to:*

Name: _____

Address:

City: _____ State: _____

Zip Code: _____

Phone: _____

Phone: _____

Phone: _____

Phone: _____

Phone: _____

Phone: _____

Phone: _____

MA083	MA084	MA085	MA086	MA087	MA088	MA089
<i>English Manual (STD)</i>	<i>German Manual</i>	<i>Italian Manual</i>	<i>Dutch Manual</i>	<i>Spanish Manual</i>	<i>French Manual</i>	<i>Portuguese Manual</i>

MODEL		English Manual (179)	German Manual (179)	Italian Manual (179)	Spanish Manual (179)	French Manual (179)	Portuguese Manual (179)		
WC20	☐	Ultralight VELOCE Wheelchair (120 kg - 265 lb weight limit)							
WC09	☐	HD Heavy Duty package (159 kg - 350 lb weight capacity)							
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber							
Frame									
SEAT WIDTH ³ - cm (inches)									
SW	☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ² ☐	48.26 (19") ² ☐		
	☐	50.80 (20") ² ☐				53.34 (21") ⁷ ☐	55.88 (22") ⁷ ☐		
SEAT DEPTH - cm (inches)									
SD	☐	35.56 (14") ¹ ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ☐	48.26 (19") ☐		
	☐	50.80 (20") ☐							
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴									
SP01	☐	Standard frame length configuration ⁵							
SP02	☐	Frame length 2.54 cm (1") longer than standard configuration							
SP03	☐	Frame length 5.08 cm (2") longer than standard configuration							
SP04	☐	Frame length 2.54 cm (1") shorter than standard configuration							
SP05	☐	Frame length 5.08 cm (2") shorter than standard configuration							
TRANSIT TIE DOWN									
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶			Short belt 50-104 cm (20" - 41")				
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶			Long belt 50-116 cm (20" - 46")				
AC35	☐	WC-19 Transit Tie Down without belt ⁶							

1- Only available with frame length 1" and 2" longer; 2- Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;

3. Width calculator


4- See chart on page 7

5- 4.375" from seat sling to front frame:

6- May interfere with swing away armrest and Clamp on height adjustable push handle: Not available with HD option: 7- Only available with HD option and the one piece adjustable angle flip-up

FRONT Seat-to-floor height

Caster sizes		Height available +/- 1/4"	
☐	3" Caster	35.56cm to 53.34cm (14" to 21")	Specify desired height: _____ (inches)
☐	4" Caster	36.83cm to 53.34cm (14.5" to 21")	
☐	5" Caster	39.37cm to 53.34cm (15.5" to 21")	
☐	6" Caster	41.91cm to 53.34cm (16.5" to 21") (N/A with 90° front frame)	




CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		
CA10	☐	Pneumatic casters (1 ¼")		6"		
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks <i>(Not available with caster pin locks)</i>				


REAR
 Seat-to-floor Height

WHEEL SIZES		Height available	Specify desired height: (inches)
☐	22" (501) wheel	34.29cm to 46.35cm (13.5" to 18.25")	
☐	24" (540) wheel	36.83cm to 48.89cm (14.5" to 19.25")	
☐	25" (559) wheel	38.10cm to 50.16cm (15" to 19.75")	
☐	26" (590) wheel	39.37cm to 51.43cm (15.5" to 20.25")	


Configur

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (½") of selected Seat to Floor height"

FRONT FRAME ANGLE						
FA70	☐	70° Front frame angle				
FA75	☐	75° Front frame angle				
FA80	☐	80° Front frame angle				
FA85	☐	85° Front frame angle				
FA90	☐	90° Front frame angle				

SEAT TO FOOTREST LENGTH	
Set Length to ¹ :	Length must be at least 4,5 cm (1 ¾") shorter than front seat-to-floor height
1 - +/- 1.27 cm (1/2")	
	

Footrest width / Taper selection					
Seat Width		Straight ¹		1" Taper	
Outside Measurement		6,4 cm Narrower than Seat Width		11,4 cm Narrower than Seat Width	
35.56	14"	29.21 cm	11.5"	24.13 cm	9.5"
38.10	15"	31.75 cm	12.5"	26.67 cm	10.5"
40.64	16"	34.29 cm	13.5"	29.21 cm	11.5"
43.18	17"	36.83 cm	14.5"	31.75 cm	12.5"
45.72	18"	39.37 cm	15.5"	34.29 cm	13.5"
48.26	19"	41.91 cm	16.5"	36.83 cm	14.5"
50.80	20"	44.45 cm	17.5"	39.37 cm	15.5"
53.34	21"	46.99 cm	18.5"	41.91 cm	16.5"
55.88	22"	49.53 cm	19.5"	44.45 cm	17.5"



1- Not compatible with auto-folding footplate

FOOTPLATES						
FP02	☐	Newton adjustable angle footplate - Composites <i>(Two piece)</i>				
FP03	☐	Newton adjustable angle - Aluminum <i>(Two piece)</i>				
FP04	☐	Newton adjustable angle footplate - Carbon fiber with rubber inserts <i>(Two piece)</i>				
FP06	☐	One piece adjustable angle flip-up - Aluminum				
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				
FP28	☐	One piece auto-folding adjustable angle - Aluminum				
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop <i>(not available with one piece footplate)</i>				
FP09	☐	Calf strap				
FP30	☐	Bodypoint Aeromesh Padded calf strap				

Backrest						
FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS						
BC03-9		Ultralight carbon fiber back cane	22.86cm (9")	*Temporarily unavailable		
BC03-12		Ultralight carbon fiber back cane	30.48cm (12")	*Temporarily unavailable		
BC03-15		Ultralight carbon fiber back cane	38.10cm (15")	*Temporarily unavailable		
BC03-18		Ultralight carbon fiber back cane	45.72cm (18")	*Temporarily unavailable		
BC03-21		Ultralight carbon fiber back cane	53.34cm (21")	*Temporarily unavailable		
BC04-13		8° Bend Carbon fiber back cane	34.29cm (13 ½)	*Temporarily unavailable		
BC04-15		8° Bend Carbon fiber back cane	38.10cm (15")	*Temporarily unavailable		
BC04-18		8° Bend Carbon fiber back cane	45.72cm (18")	*Temporarily unavailable		

* Replacement options below: 1-Carbon fiber height adjustable; 2- Fixed height aluminum; 3- carbon fiber height and angle adjustable

HEIGHT ADJUSTABLE WITH PUSH HANDLES (MUST select ONE push handle options)						
BC06-912	☐	Straight Carbon fiber back cane	22.86 - 30.48cm (9-12")	increment of 1/2"	Set at ² :	
BC06-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at ² :	
BC06-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at ² :	
BC06-1821	☐	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at ² :	
BC08-1315	☐	8° Bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	increment of 1/2"	Set at ² :	
BC08-1518	☐	8° Bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at ² :	
BC08-1821	☐	8° Bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at ² :	

FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS						
BC30	☐	Straight aluminum back cane	9"			
BC36	☐	Straight aluminum back cane	12"			
BC42	☐	Straight aluminum back cane	15"			
BC63	☐	8° bend aluminum back cane	13 ½"			
BC64	☐	8° bend aluminum back cane	15"			

ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC09-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at ² :	
BC09-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at ² :	
BC09-1821	☐	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at ² :	
BC10-1315	☐	8° bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	increment of 1/2"	Set at ² :	
BC10-1518	☐	8° bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at ² :	
BC10-1821	☐	8° bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at ² :	
		Select back cane angle:	85° ¹ ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐			

PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	☐	Push handles				
PH05	☐	Fold down push handles ^{3,5}				
PH06-920	☐	Clamp on height adjustable push handles 22.86 - 50.80cm (9" - 20") ^{3,4,5}				
PH06-923	☐	Clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3,4,5}				
PH02	☐	Omit push handles				

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 120 kg/265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes; 5- Not compatible with headrest and attendant lock

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				
BH08	☐	Adjustable tension Active-Mesh back upholstery				
BH99	☐	Omit Back upholstery				

SEAT CUSHION						
SC01	☐	5.08 cm (2") Cushion				
SC02	☐	7.62 cm (3") Cushion				
SC99	☐	Omit seat cushion				

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; **Only available with 24" wheel size; Not available with 21" and 22" wide**; 2- Only standard handrim position available; 3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels

Handrim					Available	
HR01	☐	Aluminum anodized			22", 24", 25", 26"	
HR61	☐	Aluminum black hard anodized ²			22", 24", 25"	
HR02	☐	Plastic coated ¹			22", 24", 25", 26"	
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"	
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See page 6 for the color selection and price		22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	
HR05	☐	Surge LT ^{1,2,3} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	
HR04	☐	Natural Fit ^{1,2,3}			22", 24", 25", 26"	
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"	
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"	
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"	
HR99	☐	Omit handrim				

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

Handrim option					
ASSEMBLY POSITION					
HR29	☐	Standard Handrim position	Full tab		
HR030	☐	Narrow Handrim position no cut	Full tab		
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions		
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐		
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	☐	Standard grip			
HR27	☐	Super grip			
HR099	☐	Omit thumb piece			


3- Only available with Aluminum and Air Grip handrims;

AXLE ¹					
AX01	☐	Quick release axle			
AX03	☐	Quad release axle			

1- Not required with one arm drive

REAR WHEEL CAMBER				
CM00	☐	0° camber		
CM03	☐	3° camber		
CM06	☐	6° camber		

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												
2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"	10,2 4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 3,8 cm/1.5"

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2.5").

Possible conflict between a 26" rear wheel and a center of gravity over 6,3 cm (2.5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK				
WL08	☐	Newton Wheel lock - push to lock		
WL16	☐	Aluminum push to lock with extensions		
WL02	☐	Pull to lock		
WL04	☐	Newton grade aid push to lock (<i>works best with pneumatic or high pressure tire</i>)		
WL13 WL14	☐	6" removable extension handle ¹		
WL05	☐	Scissors Lock		
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐
WL99	☐	Omit wheel locks		

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Handrims				
HEIGHT ADJUSTABLE FLIP BACK ARMREST (<i>convertible to "T" armrest</i>) ¹				
AR01-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐ Unlocked STD ☐
AR01-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐ Unlocked STD ☐
AR02-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐ Unlocked STD ☐
AR02-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐ Unlocked STD ☐
HEIGHT ADJUSTABLE "T" ARMREST ¹				
AR13-610	☐	25.4cm (10") Desk length armpad ⁴	15.25 - 25cm (6" - 10")	
AR03-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	
AR03-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	
AR14-610	☐	35.56cm (14") Full length armpad ⁴	15.25 - 25cm (6" - 10")	
AR04-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	
AR04-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²				
AR05	☐	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")	
AR23	☐	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")	
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")	
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")	
ARMREST OPTIONS				
AR12	☐	35.56cm (14") Full Length Armpad - Gel Ovation (<i>pair</i>) ³		
AR99	☐	Omit armrest		

1- May have a 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid back rests; Not available with fixed carbon fiber side guards

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (<i>not available with Flip Back and "T" Armrest</i>)				
SG01	☐	Plastic side guards (<i>removable</i>)		
SG02	☐	Carbon fiber side guards (<i>fixed</i>)		
SG03	☐	Carbon fiber side guards (<i>removable</i>)		
SG06	☐	Plastic fender side guards - straight bracket (<i>removable</i>)		
SG08	☐	Plastic fender side guards - offset bracket (<i>removable</i>) ¹		
SG04	☐	Carbon fiber fender side guards (<i>removable</i>)		
SG99	☐	Omit side guard		
SIDE GUARDS OPTION				
SG07	☐	Side guards clip (<i>pair</i>) ²		

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - Small 88.9cm (35") ¹					
SB01-L	☐	Velcro adjustable - Long 152.40cm (60")					
SB02-S	☐	Auto buckle - Small 88.9cm (35") ¹					
SB02-M	☐	Auto buckle - Medium 121.92cm (48")			New		
SB02-L	☐	Auto buckle - Long 152.40cm (60")					
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")					
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")					
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")					
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")					
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")					
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")					
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")					
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")					
SB03	☐	Hardware attachment for 2 extra points (1 pair)					
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")					
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")					
SB99	☐	Omit positioning belt					

1- Not recommended for seat width of 45.72cm (18") and more

HEADREST ¹							
AC01-S	☐	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")					
AC01-M	☐	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")					
AC01-L	☐	Profiled Headrest - Large 14x2.3cm(5.5"x8")					
AC28	☐	Headrest mounting hardware with tension bar					

1- Only available with push handles

ACCESSORIES							
AC03	☐	Anti-tippers					
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐			
AC27	☐	Tip assist ²					
AC06	☐	Cane and crutch holder	Left ☐	Right ☐			
AC08	☐	Spoke guards ³	Available: 22", 24", 26"				
AC10	☐	Seat pouch					
AC30	☐	Newton multi-purpose pouch					
AC26	☐	BodyPoint Mobility Bag					
AC23	☐	Newton Backpack					
AC12	☐	Velcro Closing Strap					
AC13	☐	Tool kit					
AC16	☐	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐		
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐
						M/XXXL ☐	

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerity

Colors																				
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐	\$55	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
																				
Highlight ²	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color	Specify PPG Paintit color:																	additional 10 days lead-time

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_intl/products/folding-wheelchairs/veloce

Frame length chart

Frame length	Seat Depth						
	35.56(14")	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")
STD	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")
1	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")	N/A
2	40.64(16")	43.18(17")	45.72(18")	48.26(19")	50.80(20")	N/A	N/A
-1	N/A	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")	48.26(19")
-2	N/A	N/A	N/A	38.10(15")	40.64(16")	43.18(17")	45.72(18")

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4

Notes :