

# Vigour Back Support

---

Lightweight yet tough to meet the daily lifestyle demands of active wheelchair users. Back support angle and seat depth offers exceptional adjustability and can be facilitated with the user in the wheelchair. Ventilated spacer-fabric and shell system increases air-flow through the back offering daily comfort.

## FUNCTIONAL BENEFITS

Solid full height back support available in three heights for maximum trunk and shoulder support.

Designed for the individual who requires:

- A solid back support
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability



# Vigour Mid Back Support

---

The Spex Vigour Mid Back Support offers the perfect balance of stability and freedom for active wheelchair users. With its lightweight yet durable design, it provides essential lower trunk support while allowing for efficient self-propulsion. Its compatibility with lateral supports and additional positioning accessories makes it a versatile choice for those needing tailored postural support.

## FUNCTIONAL BENEFITS

Designed for movement, built for support.

- Compatible with lateral supports – The Vigour Mid is compatible with accessories to allow for additional trunk stability
- Lower trunk coverage – Offers support at the base to assist with pelvic positioning, especially when paired with the Spex Tessellated Shape System
- Highly adjustable – Seat depth and back angle can be fine-tuned while the user remains in the wheelchair
- Available in scripted and custom sizes – Ensuring the perfect fit for every user
- Planar support – A simple, flat back support allowing for easier customisation with accessories






Dual-layer cover reduces **SHEAR FORCES** with air-mesh side panels to counter moisture and temperature

Can be upgraded with **DIFFERENT FOAM PADDING OPTIONS** for excellent results

**SLOTS & HOLE POSITIONS** on shell allow easy installation of head supports, lateral supports and harnesses

Add laterals as required from our **LARGE MODULAR RANGE**

## SIZE CHART



- Squares are standard sizes.
- Non-standard or custom sizes are available on request.

## PRODUCT DETAILS

| Width   | Depth | Capacity       | Product Weight                          |
|---------|-------|----------------|-----------------------------------------|
| 10"     | 2"    | 200lbs (90kg)  | Vigour Back (16x16)<br><b>1.3kg</b>     |
| 11"-15" | 2.5"  |                | Vigour Mid Back (16x16)<br><b>1.0kg</b> |
| 16"-17" | 2.75" | 300lbs (136kg) |                                         |
| 18"-20" | 3"    |                |                                         |

- Note: This is the weight of the back support only and does not include the hardware. For hardware weight, refer to page 29.

## PRODUCT PERFORMANCE

|                      |  |  |   |   |  |
|----------------------|--|--|---|---|--|
| Postural Control     |  |  |   | 4 |  |
| Lateral Stability    |  |  | 3 |   |  |
| Postural Asymmetries |  |  | 3 |   |  |
| Changing Needs       |  |  | 3 |   |  |
| Pressure Care        |  |  | 3 |   |  |

- Note: This is a guide only and ratings are dependent on correct use of product.

- Vigour Back Support
- Vigour Mid Back Support