

Gravity Chair



INDOOR SEATING

SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		



Size Options

Measurements	Small	Medium	Large	X-Large	XL-Wide	User's Measurements
Order Code	3770-0021-XXX	3770-0022-XXX	3770-0023-XXX	3770-0024-XXX	3770-0025-XXX	
A Seat Depth	24cm	30cm	38cm	43cm	43cm	
B Seat Width at Back	21cm	23cm	27cm	30cm	41cm	
C Seat Width at Front	28cm	30cm	37cm	39cm	50cm	
D Back Height	60cm	69cm	74cm	81cm	81cm	
E Seat Front Height	20cm	23cm	26cm	29cm	29cm	
F Arm Support Height	29cm	35cm	39cm	43cm	43cm	
G Total Width of Chair	46cm	49cm	57cm	60cm	71cm	
H Total Height of Chair	55cm	61cm	67cm	71cm	72cm	
I Total Depth of Chair	60cm	71cm	82cm	91cm	91cm	
J Tilt Wedge Length	47cm	57cm	67cm	74cm	74cm	
K Tilt Wedge Height	10cm	12cm	14cm	17cm	17cm	
Weight Limit	100kg					

Gravity Chair Package includes:


Fully wipe-down, including internal contact surfaces where user is seated.

- Foam chair with cover
- Coloured liner
- Standard head support
- Tilt wedge
- Velcro side pocket



Our product codes are meaningful!

Here's how to generate a product code to place your order.








3 7 7 0 - 0 0 X X - 0 X X






EXAMPLE

A medium Gravity Chair with a green wipe-down liner, would have the code 3770-0022-055.


Cover/Liner Options

Classic Version		Wipe Down Liner	
Air breathable fabric on the user contact surfaces and rear back panel, with wipe-down covers in all other areas 3770-001X-XXX		Fully wipe-down, including internal contact surfaces where user is seated 3770-002X-XXX	
Red Liner	Green Liner	Blue Liner	
-022 	-055 	-066 	

Head Supports


Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm internal width)		Small (17cm internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
Wipe Down Cover 2010-4110-209		Large (21cm internal width)		Large (21cm internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
Wipe Down Cover 2010-4110-209		Large (21cm internal width)		Large (21cm internal width)	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Mobile Base

For ease of maneuverability in the home or classroom					
Size 1 3780-8701-000	Size 2 3780-8702-000	Size 3 3780-8703-000	Size 4 3780-8704-000	Size 5 3780-8705-000	

Tray Table

Encourages hand-to-eye coordination

<p>To fit Size Small and Medium 3770-7201-000</p>	<p>To fit size Large and X-Large 3770-7202-000</p>	<p>To fit size XX-Large 3770-7203-000</p>	
---	--	---	---

Foot Support


<p>Small (Suitable for Small or Medium Gravity chairs) 3770-1031-000</p>		<p>Large (Suitable for Large, XL and XLW Gravity Chairs) 3770-1032-000</p>	
---	---	---	---

Abductor Pommel

<p>Small (Suitable for Small or Medium Gravity chairs) 3770-1021-000</p>		<p>Large (Suitable for Large, XL and XLW Gravity Chairs) 3770-1022-000</p>	
---	---	---	---

Side Filler Pads

For reducing width of Gravity Chair by up to 6cm, adjustable in 1cm increments

<p>Small 3770-1011-000</p>	<p>Medium 3770-1012-000</p>	<p>Large 3770-1013-000</p>	<p>XL/XLW 3770-1014-000</p>	
---------------------------------------	--	---------------------------------------	--	---

Depth Adjustment Liner

For reducing the depth of the Gravity Chair

<p>Small 3770-0321-000</p>	<p>Medium 3770-0322-000</p>	<p>Large 3770-0323-000</p>	<p>XL 3770-0324-000</p>	<p>XLW 3770-0325-000</p>	
---------------------------------------	--	---------------------------------------	------------------------------------	-------------------------------------	---

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

26 Pardoe Street
Devonport TAS 7310
Australia

1300 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Gravity_20230711