

QUICKIE[®]

iXpress

All the convenience & transportability of a manual wheelchair

All the power & speed of a power chair

Quickie iXpress:
The power you need...



The power you need...

The iXpress offers a high performing, reliable Power Add On System, whether you are going shopping, going to the movies or just hanging around with your family & friends, the iXpress will truly extend lifes little pleasures.

Designed for Performance & Reliability



iXpress fitted on Quickie 2



iXpress fitted on Quickie Iris

With the Quickie iXpress you can Xpress the distance travelled, Xpress the speed of travel and Xpress all your transport needs.



Electronics Options:

- Japanese designed, highly reliable electronics
- Selectable speeds for easy transitions between environments
- Easy to use ergonomic joystick controller
- Full function attendant control with innovative glide to steer function

Quickie iXpress Specifications Adult

SPECIFICATIONS	QUICKIE 2	QUICKIE IRIS
Drive system	Direct rear wheel drive	Direct rear wheel drive
Controller	JWX-1	JWX-1, JUX-1
Motor	24 V, 120 Wx2	24 V, 120 Wx2
Battery	Nickel Metal Hydride	Nickel Metal Hydride
Range	16 km per charge	16 km per charge
Charging time	3 hrs	3 hrs
Attendant Control	Yes	Yes
Drive wheel size	24"	16", 24"
Speed	up to 6kph	up to 6kph
Weight capacity	120 kg	120 kg
Battery weight	2.9 kg	2.9 kg
System weight	17.4 kg	17.4 kg
iXpress Warranty	2yrs	2yrs
Frame Warranty	Limited Lifetime	Limited Lifetime