

# CentroGlide Lift-Anterior-Reach (LAR)



## Freedom of movement.

Based on our leading powered wheelchair, the CentroGlide, this variant adds a power lift and anterior reach function. The controls allow the user to vertically raise or lower the seat to the right height to be comfortable at a desk, high bench or low coffee table.

CentroGlide LAR has the additional anterior tilt function and a reach forward option to allow for a standing transfer smoothly.

## New levels of independence.

CentroGlide LAR enables remarkable manoeuvrability in tight spaces, and stability over kerbs. The precisely controlled tilt, lift and anterior reach movements offer unprecedented anatomical function and comfort.

## Designed for comfort and ease of use.



### Customised seating.

For ultimate comfort, many seating options are available, including customised cushioning and backrests.



### Active suspension.

Our patented 6-wheel APT suspension improves traction and aids ascending and descending kerbs safely without affecting the stability and angle of the seat.



### Wheels.

CentroGlide has 14" drive wheels, 8" and 7" castor wheels as standard. Other sizes available.



### Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



### Electronics.

Fully programmable R-net joystick control system. Compatible with many other 3rd party speciality electronics.



### Power lift.

Power lift can lift 50mm to 150mm before going into anterior tilt.



### Power anterior tilt.

The seat can be tilted forward from -4° to 45°



### Battery range.

The powerful battery gives you a range of up to 40km on a single charge.



# CentroGlide LAR sets the benchmark.

## Features.

- Power base with removable seating module
- High quality German 2 or optional 4 pole high torque motors
- Vehicle tie down points
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) kerb climbing ability
- 80mm ground clearance
- Max ramp incline/decline 10° (as per AS3695 testing)
- 4.1/3.50-6 (14") drive wheels
- Electronic brakes with manual park
- R-net LCD Std controls

## Safety.

- Full compliance AS 3695 150kg
- Crash tested to AS 3696:19 to max user weight 150kg

## Arrange a test drive.

Chat to a Glide Wheelchair specialist and experience the extraordinary features of the CentroGlide in person. Find your nearest distributor at [glide.com.au/nearest](http://glide.com.au/nearest)

## Specifications.



**Recommended User Weight**  
Up to 150kg AS 3695 full compliance



**Overall Length**  
900mm (foot plate removed)



**Seat Widths**  
31 to 61cm



**Overall Width**  
620mm (wheel base)



**Seat Depths**  
31 to 58cm



**Seat Height**  
46cm (top of rail using 0-46 set up)



**Back Height**  
1030mm Top of push handle (custom height available)



**Arm Rests**  
Flip up (adjustable)



**Battery**  
Up to 40km range



**Leg Rests**  
Power elevating centre mount"



**Frame Colour**  
Wide choice



**Brakes**  
Electronic and manual park