

Dandy Stander

A popular standing aid for children with mild to moderate postural needs. The compact design makes it user friendly in any home and school environment.

Complete Solution for Instant Standing!



No Tools Adjustment!

Knobs make adjusting the frame and supports quick and easy.



Robust frame

Keeps the child's centre of gravity close to the floor for safety. The Dandy Stander has separate chest and pelvic supports, which are adjustable in circumference, height, and direction. This is beneficial as once the child has progressed, the chest support can then be eliminated.

COMPLETE PACKAGE

The Dandy stander comes as a complete package fitted with:

- Knee Supports
- Body Supports
- Tray table
- Foot positioners



FEEDBACK

“ Emma is a clever four year old who has Spina Bifida. She has always been nervous of standing frames... until we tried the Dandy Standing frame! She loves this frame because it's fun-looking, colourful, and feels safe!

Because she has good upper body strength, the chest support was removed for unrestricted activities and learning.

Her mother and Physiotherapist love how adjustable the Dandy is, and how it doesn't look 'medicalised'. ”

Size	Code	Age	Max User Weight	Approx User Height	Tilt Range	Footrest - Table Height	Footrest - Pelvic Support	Footrest - Chest Support	Kneepads - Footrest
1	5061-0000-000	1-4yrs	35kg	85 cm	0° - 19°	51-71 cm	39-57 cm	39-57 cm	18-31 cm
2	5062-0000-000	4-8yrs	45kg	125 cm	0° - 14°	62-82 cm	44-70 cm	60-86 cm	23-44 cm
3	5063-0000-000	6-15yrs	145 cm 55kg	145 cm	0° - 11°	71-93 cm	56-91 cm	76-110 cm	23-54 cm