

Endeavour Stroller

The Endeavour stroller provides families and caregivers the opportunity to enjoy life with a new perspective of adventure.

A robust outdoor stroller solution for larger users up to 113kg



The Endeavour folds up compactly with either the wheels removed or kept on. This makes it easy to transport and for stowing in the boot of your car!



EXPLORE THE OUTDOORS!

The Endeavour is perfect for neighbourhood walks, long runs, mountain trails, beaches, family barbecues, theme parks, and more!



ACCESSORIES

Standard package includes:
Red stroller with sun canopy, 5-point harness, adjustable foot boxes, carry bag, and pneumatic bicycle wheels.

Front Castor Wheels

Castor wheel attachment for greater maneuverability when required.
2662-8410-000

Rain Cover

Plastic covering for protection from the rain.
2&3: 2662-7600-210
4: 2662-7604-214



Travel Bag

For easy storage
2&3: 2662-8180-000
4: 2662-8184-000



Size	Code	Approx Age of Child	Max Load	Seat Depth	Seat Width	Shoulder Height	Footplate Height	Backrest Angle	Approx Stroller Weight	Approx Stroller Dimensions
2	2662-0002-020	11 yrs	68 kg	38 cm	38 cm	30-58 cm	38-56 cm	30-90°	21 kg	155 x 61 x 109cm
3	2662-0003-020	14 yrs	91 kg	38 cm	46 cm	38-66 cm	38-61 cm	30-90°	21 kg	155 x 69 x 109cm
4	2662-0004-020	teenager - large adult	113 kg	46 cm	53 cm	43-74 cm	43-61 cm	30-90°	25 kg	175 x 74 x 135cm