


Hibbot Walker



**INNOVATIVE
WALKING DEVICE**


Client Details

Date:	
Client/User Name:	
NDIS/Funding:	
Further Information	
<hr/> <hr/>	

Measurements Required

Date of Measurement	
Weight	kg

1	Circumference around the hips, widest part	cm
2	Length from spina iliaca anterior superior (SIAS) to floor	cm
3	Total height of the user	cm



Growth Progress Guide

Wheel Size (Inch)	Pelvic Height (cm)	Axle Shift (cm)	Date
16"	25	7	
	26	9	
	27	10	
	28	12	
	29	13	
	30	15	
	31	16	
	32	18	
	33	19	
	34	21	
20"	35	22	
	36	7	
	37	9	
	38	11	
	39	12	
	40	14	
	41	16	
	42	17	
	43	19	
	44	21	
24"	45	22	
	46	7	
	47	9	
	48	11	
	49	12	
	50	14	
	51	16	
	52	17	
	53	19	
	54	21	
29"	55	22	
	56	7	
	57	9	
	58	10	
	59	11	
	60	12	
	61	14	
	62	15	
	63	16	
	64	17	
65	19		
66	20		
67	21		
68	22		

Brace Size (cm)	XS		S		M		L	
Pelvic-circum	38	43	43	48	48	54	54	61

Example Patient:
 A 2 year old girl, 82 cm in height with a pelvic height of 32 cm will have 16" wheels shifted at 18 cm back from the base line.
 The circumference around the pelvic is 40 cm, which means she requires a XS corset.

Sizing

Specifications	Braces	Wheels
User Height: 60 - 125 cm User Weight: 30 kg maximum Maximum Total Width: 70 cm Maximum Total Length: 110 cm Maximum Total Height: 106 cm	X - Small: 38 - 43 cm Small: 43 - 48 cm Medium: 48 - 54 cm Large: 54 - 61 cm	16" x 1.75" 20" x 1.75" 24" x 1.75" 29" x 1.75"

Select Package

All Options Include (must be configured)

Hibbot Walking Device
 Pair of Wheels
 Brace

5810-0010-000

Select Wheels (Required)

Choose Wheel Size (Pair)

16" (M4: 25-35 cm) 5810-0116-000	<input type="radio"/>	24" (M4: 46-55 cm) 5810-0124-000	<input type="radio"/>
20" (M4: 36-45 cm) 5810-0120-000	<input type="radio"/>	29" (M4: 56-68 cm) 5810-0129-000	<input type="radio"/>

Select Hibbot Brace (Required)

Choose Brace Size

XS (M3: 38-43 cm) 5810-0190-000	<input type="radio"/>	M (M3: 48-54 cm) 5810-0192-000	<input type="radio"/>
S (M3: 43-48 cm) 5810-0191-000	<input type="radio"/>	L (M3: 54-61 cm) 5810-0193-000	<input type="radio"/>

Select Weights (Recommended)

Weights (Set of Two)

Set of weights (2 x 1.25kg & 2 x 2.5kg)
5810-8431-000

Select Additional Wheels

Choose Wheel Size (Pair)			
16" (M4: 25-35 cm) 5810-9516-000	<input type="text" value="INPUT QUANTITY"/>	24" (M4: 46-55 cm) 5810-9524-000	<input type="text" value="INPUT QUANTITY"/>
20" (M4: 36-45 cm) 5810-9520-000	<input type="text" value="INPUT QUANTITY"/>	29" (M4: 56-68 cm) 5810-9529-000	<input type="text" value="INPUT QUANTITY"/>

Select Additional Hibbot Brace

Choose Brace Size			
XS (M3: 38-43 cm) 5810-1900-000	<input type="text" value="INPUT QUANTITY"/>	M (M3: 48-54 cm) 5810-1902-000	<input type="text" value="INPUT QUANTITY"/>
S (M3: 43-48 cm) 5810-1901-000	<input type="text" value="INPUT QUANTITY"/>	L (M3: 54-61 cm) 5810-1903-000	<input type="text" value="INPUT QUANTITY"/>

Additional Client Notes

Notes

medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: HibbotScriptForm20200428