

# Jay Ion™ Cushion

The JAY® Ion™ is a comfortable skin protection and positioning cushion composed of a multi-layered foam base, moisture resistant inner cover, and an anti-microbial outer cover.

### COMMON CLIENT CHARACTERISTICS

- May be at moderate risk of skin breakdown
- May be able to perform independent weight shifts
- May have moderate positioning needs
- May be in a wheelchair for extended periods of time



Ion Rating Scale	
Side-to-Side Stability	3
Forward Stability	3
Pressure Management	3

Component	Feature	Design Intent	Clinical Value
BASE	Multilayered foam base with optimized densities by width	Optimizes cushion performance by user weight ranges	Encourages improved orthopedic alignment and accommodates individual user needs
	Soft Visco™ memory foam top layer	Enhances comfort	Increases sitting tolerance
INNER COVER	Inner cover with Aquaquard™ zipper and anti-wicking thread	Protects foam from moisture build up and is easy to clean	Helps with incontinence management
COVER	Foam spacer layer	Dissipates heat and moisture	Reduces risk of skin breakdown
	X-static™ silver thread stretch outer cover	Reduce surface tension and resist bacterial growth	Helps prevent bacterial growth and reduces risk of skin breakdown

Technical Specifications	
Waterlow	<15
Width Range	14" to 24"
Depth Range	14" to 24"
Weight	159kgs (18"x18")
Height	4"
Weight Capacity	136kgs 227kgs (22-24" widths)
JAY Your Way Modifications Avail.	Yes
Positioning	No

		WIDTH												
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
DEPTH	12"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	13"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	14"	N/A	N/A	•	•	•	•	•	N/A	N/A	N/A	N/A	N/A	N/A
	15"	N/A	N/A	•	•	•	•	•	•	N/A	N/A	N/A	N/A	N/A
	16"	N/A	N/A	•	•	•	•	•	•	•	N/A	N/A	N/A	N/A
	17"	N/A	N/A	•	•	•	•	•	•	•	•	N/A	N/A	N/A
	18"	N/A	N/A	•	•	•	•	•	•	•	•	•	N/A	N/A
	19"	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•	•
	20"	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•	•
	21"	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•	•
	22"	N/A	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•
	23"	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	•	•	•
24"	N/A	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	•	•	

• Standard size available    JYW Size available through JAY Your Way    N/A Size not available