


JAY 2DEEP CONTOUR™

The JAY® J2[™] Deep Contour is a pre-contoured foam cushion that features a 3" fluid immersion in a JAY Flow[™] fluid tripad for superior skin protection, and an easy to modify base for optimal stability.



JAY 2 DEEP CONTOUR™ CUSHION

The JAY® J2™ Deep Contour is a pre-contoured foam cushion that features a 3" fluid immersion in a JAY Flow™ fluid tripad for superior skin protection, and an easy to modify base for optimal stability.

Contoured Foam Base

The Jay 2 Deep Contour cushion features a contoured foam base made from closed cell foam, this makes it easy to modify in the field or from our factory. Additionally extensive array of positioning components is available.

Jay Flow Fluid Bladder

The Jay 2 Deep Contour is available with a Jay Flow Fluid bladder for clients at high risk of skin breakdown. It features 3" of fluid immersion for maximum skin protection.

Stretch, Air Exchange or Incontinent cover options

Moisture-Resistant Cover with No-Slip Bottom

SPECIFICATIONS



JAY 2 Deep Contour Cushion	
Product Width:	14 - 24"
Weight:	2.8 kgs
Product Depth	14 - 20"
Product Height	4.5"
Weight Capacity	113 kgs
Cushion Type	Polyuerethene base with fluid bladder
Outer Cover	Ballistic stretch, Air Exchange, incontinence
Warranty	24 months

