
# 

### THE MOST VERSATILE JAY BACK FOR ALL SHAPES, SIZES, AND LIFESTYLES

The JAY<sup>®</sup> J3<sup>™</sup> wheelchair backrest product line offers a multitude of width, height, and contour depth back shells to fit the back to the client. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.



## **J3<sup>®</sup>Back**



Five Contour Shapes provide posterior and lateral pelvic and lateral thoracic stability for clients with varying postures and functional needs:

• **PA** - 2.2"

- PL 3"
- PD, PDL, PDC 6"
- Four Support Height Levels provide support for clients with wide ranges of torso heights and levels of injury:
- Lower Thoracic 6.5" to 9.5"
- Mid Thoracic 12" to 15"
- Upper Thoracic 16.5" to 19.5"
- Shoulder 21" to 24"



- Four Hardware Options to accommodate a span of client priorities and wheelchair configurations:
- SH easy to install, 2-point, quick-release hardware
- **RH** 2-point, quick-release hardware that is 25% lighter than SH hardware and optimized for ultra lightweight wheelchairs
- TH 2-point, quick-release hardware optimized for tilt-in-space wheelchairs, with lengthened bracket and attachment pin that accommodate up to 500 lbs.
- 4-point quick-release hardware that distributes force across four mounting locations



#### Carbon Fiber Option for

the ultimate in lightweight performance. Weighs only 2.2 lbs. with RH hardware.

#### **J3 BACK SPECS**

HCPCS Code E2613, E2614, E2615, E2620

Width Range 12" to 26"

Height

Lower Thoracic: 6.5", 8", 9.5" Mid Thoracic: 12", 13.5", 15" Upper Thoracic: 16.5", 18", 19.5" Shoulder Height: 21", 22.5", 24"

Weight Capacity

12" to 20" widths: 300 lbs. 22" to 26" widths: 500 lbs.

Hardware Adjustability

Up to +/- 20° angle adjustability Up to +2" width adjustability Up to +2.25" depth adjustability

Lateral Contour Support PA: 2.2" PL: 3" PD. PDL. PDC: 6"

Optional Accessories Lateral Thoracic Pads Chest Straps/Harnesses Whitmyer® Head Supports Spine Align Posture Conforming Inserts





800.333.4000 www.SunriseMedical.com