

Jay Union™ Cushion

The JAY® Union™ is a versatile, comfortable, skin protection and positioning cushion composed of a dynamic fluid and foam layering system, moisture resistant inner cover and an anti-microbial outer cover.

COMMON CLIENT CHARACTERISTICS

- May be at moderate to high risk of skin breakdown
- May be able to perform independent weight shifts
- May have moderate positioning needs
- May be in a wheelchair for extended periods of time
- Requires No Maintenance Cushion
- Requires Soft Comfortable Cushion periods of time



Reduced Profile Option

Union Rating Scale		
Side-to-Side Stability		4
Forward Stability		4
Pressure Management		3.5

Component	Feature	Design Intent	Clinical Value
BASE	Pre-contoured, high-resiliency foam with additional lateral pelvic support	Provides positioning of the pelvis and thighs	Encourages improved orthopedic alignment
	Soft Visco™ memory foam top layer with extra layer over pelvic loading area	Enhances comfort	Increases sitting tolerance
INSERT	JAY® Flow™ fluid tripod	Consistently ensures proper fluid placement by reducing fluid migration	Reduces risk of skin breakdown
INNER COVER	Inner cover with Aquaquard™ zipper and anti-wicking thread	Protects foam base from moisture build up and is easy to clean	Helps with incontinence management
COVER	Foam spacer layer	Dissipates heat and moisture	Reduces risk of skin breakdown
	X-static™ silver thread stretch outer cover	Reduces surface tension and resists bacterial growth	Helps prevent bacterial growth and reduces risk of skin breakdown

Technical Specifications	
Waterlow	<18
Width Range	14" to 24"
Depth Range	14" to 24"
Weight	2.4kgs (16"x18")
Height	4.5" (3.5" Reduced Profile)
Weight Capacity	136kgs 227kgs (22-24" widths)
JAY Your Way Modifications Avail.	Yes
Positioning	No

		WIDTH												
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
DEPTH	12"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	13"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	14"	N/A	N/A	•	•	•	•	•	N/A	N/A	N/A	N/A	N/A	N/A
	15"	N/A	N/A	•	•	•	•	•	•	N/A	N/A	N/A	N/A	N/A
	16"	N/A	N/A	•	•	•	•	•	•	•	N/A	N/A	N/A	N/A
	17"	N/A	N/A	•	•	•	•	•	•	•	•	N/A	N/A	N/A
	18"	N/A	N/A	•	•	•	•	•	•	•	•	•	N/A	•
	19"	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•	•
	20"	N/A	N/A	N/A	•	•	•	•	•	•	•	•	•	•
	21"	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	•	•	•
	22"	N/A	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	•	•
	23"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	•
24"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	•	•	•	•	•	

• Standard size available JYW Size available through JAY Your Way N/A Size not available