


Jenx Monkey



PRONE STANDER

Client Details

Client/User Name:	
Date:	
DOB:	
User Measurements	
A: Hip Width	
B: Axilla Height	
C: Chest Width	
D: Overall User Height	
User Weight	



Product Specifications			
Average Age	9 Months - 4 Years	Chest Height	51cm - 75cm
Weight Capacity	25 kg	Max Height	-
Hip Width	18.5cm - 28.5cm	Tilt Range (90° = Upright)	20° - 90°
Hip Height	36cm - 60cm	Frame Width	58cm
Chest Width	18.5cm - 28.5cm	Frame Length	88.5cm

Frame

Includes

Jenx Prone Stander, Chest Supports, Hip Supports, Leg Supports, Foot Straps, Sandals (size 1), Tray, Tray table bowl, Tray table infill

5050-0001-000



Accessories

Activity Frame	Grab Bar	Size 3 Sandal Raisers
For attaching toys or other items	Perfect for children who need to fix one hand to increase functionality	For leg length discrepancy
3615-7475-000	5050-7201-000	5356-6091-000
350mm High		4 x 15mm Blocks
		

medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Jenx_Monkey_20211108