



Quickie 7^R
2018
EIR11

medifab

Mark For:

Date: _____
Therapist: _____
Customer Name: _____
Contact: _____
Address: _____
City: _____ ST: _____ PC: _____
Phone: () _____ Fax: () _____
Confirmation Email: _____

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
ST: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME

Base Model

<input type="checkbox"/>	132WM88	Quickie 7 ^R - Adjustable (120kgs)
--------------------------	---------	--

Performance Package (No substitutions, discount given on separate line)

<input type="checkbox"/>	132QD01	Single Sided Fork
		Ultra Lightweight Wheel
		Schwalbe Marathon Plus Evolution
		Compact Wheel Lock
		Retail Value
		Package Price

Ultra Lite Package (No substitutions, discount given on separate line)

<input type="checkbox"/>	132QD02	Ultra Lite Seat Sling
		Ultra Lightweight Wheel
		Schwalbe Speed Run Tire
		Compact Wheel Lock
		Retail Value
		Package Price

Frame Angle

See **FIGURE 1** Measure from front frame to floor.

<input type="checkbox"/>	132FA1	75°	NC
<input type="checkbox"/>	132FA8	80°	NC
<input type="checkbox"/>	132FA2	85°	NC

Frame Inset Per Side

See **FIGURE 2** Frame Inset is the distance from the outside of the seat tube to the outside of the footrest.

<input type="checkbox"/>	132FF1	0"	NC
<input type="checkbox"/>	132FF2	1" (NA w/ 12" wide)	NC
<input type="checkbox"/>	132FF3	2" (NA with 12" & 13" wide)	NC

Frame Style

See **FIGURE 3** Measurement for set up.

<input type="checkbox"/>	132FS9	Standard	NC
<input type="checkbox"/>	132FS34	Classic	NC

* **Standard Frame** - Seat Sling to front frame bend angle is 1.5" shorter than the classic position. **Classic Frame** - Traditional seat sling to front frame bend angle dimensions. Front of back cane post to front edge of seat sling will remain the same for both Standard and Classic.

Frame - Direct Mount Wheel Lock (Optional)

Only available with Compact Wheel Lock, 0° & 3° Camber and 24-25 rear wheels. No mounting clamp hardware, NA with Aluminium or carbon fiber sideguards. NA with mountain tires or rear suspension. NA in 12"-15" seat depths. NA with height adjustable armrests.

<input type="checkbox"/>	132WH13	Direct Mount Wheel Lock	NC
--------------------------	---------	-------------------------	----

Transit

<input type="checkbox"/>	132TO45	Transit System
--------------------------	---------	----------------

Attitude Bike

<input type="checkbox"/>	132IN04	Attitude Hand Bundle	See Order Form
--------------------------	---------	----------------------	----------------

Figure 1




Figure 2




Figure 3




2. SIZE / SEATING

Depth Adjustable Backrest (Optional)

See **Figure 4** Depth Adjust Backrest extends the frame rear of the seat sling, select a factory setting.

<input type="checkbox"/>	132DP02	0 Setting Backrest at standard position, field adjustable to -1" or -2"	
<input type="checkbox"/>	132DP03	-1 Setting Backrest -1" offset, field adjustable to flush or -2"	
<input type="checkbox"/>	132DP04	-2 Setting Backrest -2" offset, field adjustable to flush or -1"	

Figure 4



Seat Width

See **FIGURE 5** Measure from outside of seat tube to the outside of opposite seat tube.

<input type="checkbox"/>	132S3	12"	NC	<input type="checkbox"/>	132S8	17"	NC
<input type="checkbox"/>	132S4	13"	NC	<input type="checkbox"/>	132S9	18"	NC
<input type="checkbox"/>	132S5	14"	NC	<input type="checkbox"/>	132S10	19"	NC
<input type="checkbox"/>	132S6	15"	NC	<input type="checkbox"/>	132S48	20"	
<input type="checkbox"/>	132S7	16"	NC				

Seat Depth

See **FIGURE 6** Measure from front of back post to front edge of seat sling. Seat sling is 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

<input type="checkbox"/>	132SD5	12"	NC	<input type="checkbox"/>	132SD10	17"	NC
<input type="checkbox"/>	132SD6	13"	NC	<input type="checkbox"/>	132SD11	18"	NC
<input type="checkbox"/>	132SD7	14"	NC	<input type="checkbox"/>	132SD12	19"	NC
<input type="checkbox"/>	132SD8	15"	NC	<input type="checkbox"/>	132SD13	20"	
<input type="checkbox"/>	132SD9	16"	NC				

Figure 5





Figure 6



Seating Systems

<input type="checkbox"/>	132OP01	Standard Seat Sling	NC
<input type="checkbox"/>	132OP02	Ultra Seat Sling Available 14" to 18" Widths and Depths. Ultra Seat Sling will be ~1/4" less than requested seat height value.	
<input type="checkbox"/>	132OP04	Omit Seat Sling	NC

Cushion

<input type="checkbox"/>	132CU1	2" Cushion	
--------------------------	--------	------------	--

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	132JC20	JAY Bundled Cushion Cushion Name or Part # Width: _____ Depth: _____	Price List
--------------------------	---------	--	------------


3. CASTERS / SEAT HEIGHT

Front Seat Height

See **FIGURE 7** Measure from the front of the seat rail to the floor.

<input type="checkbox"/>	132SH14	16.0"	NC	<input type="checkbox"/>	132SH20	18.5"	NC
<input type="checkbox"/>	132SH15	16.5"	NC	<input type="checkbox"/>	132SH21	19.0"	NC
<input type="checkbox"/>	132SH16	17.0"	NC	<input type="checkbox"/>	132SH23	19.5"	NC
<input type="checkbox"/>	132SH17	17.5"	NC	<input type="checkbox"/>	132SH24	20.0"	NC
<input type="checkbox"/>	132SH18	18.0"	NC	<input type="checkbox"/>	132SH26	20.5"	NC
				<input type="checkbox"/>	132SH27	21.0"	NC

Figure 7



Caster Wheels

<input type="checkbox"/>	132CW1	3" Micro	NC
<input type="checkbox"/>	132CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	132CW31	4" Micro Lighted	
<input type="checkbox"/>	132CW72	4" x 1" Aluminium Poly	
<input type="checkbox"/>	132CW44	4" x 1.25 Semi Pneumatic	
<input type="checkbox"/>	132CW64	4" x 1.5" Aluminium Soft Roll	
<input type="checkbox"/>	132CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	132CW51	5" Micro Lighted	
<input type="checkbox"/>	132CW6	5" Soft Roll (Grey) -NA w/Frog Legs	
<input type="checkbox"/>	132CW65	5" x 1" Aluminium Poly	
<input type="checkbox"/>	132CW7	5" x 1.5" Semi Pneumatic	
<input type="checkbox"/>	132CW66	5" x 1.5" Aluminium Soft Roll	
<input type="checkbox"/>	132CW34	6" Poly	NC
<input type="checkbox"/>	132CW33	6" Pneumatic	
<input type="checkbox"/>	132CW35	6" x 1.5" Semi Pneumatic	
<input type="checkbox"/>	132CW68	6" x 1.5" Aluminium Soft Roll	

Caster Fork and Suspension

<input type="checkbox"/>	132CF65	Multiple Position	NC
<input type="checkbox"/>	132CF69	Multiple Position Single Sided Fork	
<input type="checkbox"/>	132CF53	Frog Legs 0 weight to 45kgs	
<input type="checkbox"/>	132CF54	Frog Legs 46kgs weight to 86kgs	
<input type="checkbox"/>	132CF55	Frog Legs 87kgs weight to 120kgs	

3. CASTERS / SEAT HEIGHT (CONT)

Caster Housing Position

See **FIGURE 8**

<input type="checkbox"/>	132LH01	Setting 1 Suggested for 3" & 4" casters. For 2" Inset NA w/ setting 1	NC
<input type="checkbox"/>	132LH02	Setting 2 Suggested for 5" & 6" casters.	NC
<input type="checkbox"/>	132LH03	Setting 3 NA w/ 6" casters	NC




Figure 8

4. REAR WHEELS / SEAT HEIGHT

Rear Seat Height

RSH 13"-15.5" only w/ 20" Rear Wheel.

RSH 15"-15.5" Only w/ 20, 22 or 24" w/HP Clincher Rear Wheel.

RSH 16"-16.5" NA w/ 24, 25 & 26" Rear Wheels and Rear Suspension

<input type="checkbox"/>	132RH11	13.0" (NA w/ Rear Suspension)	NC
<input type="checkbox"/>	132RH12	13.5" (NA w/ Rear Suspension)	NC
<input type="checkbox"/>	132RH13	14.0"	NC
<input type="checkbox"/>	132RH14	14.5"	NC
<input type="checkbox"/>	132RH15	15.0"	NC
<input type="checkbox"/>	132RH16	15.5"	NC
<input type="checkbox"/>	132RH17	16.0"	NC
<input type="checkbox"/>	132RH18	16.5"	NC
<input type="checkbox"/>	132RH19	17.0"	NC
<input type="checkbox"/>	132RH26	17.5"	NC
<input type="checkbox"/>	132RH20	18.0"	NC
<input type="checkbox"/>	132RH27	18.5"	NC
<input type="checkbox"/>	132RH21	19.0"	NC
<input type="checkbox"/>	132RH37	19.5"	NC
<input type="checkbox"/>	132RH22	20.0" (NA w/ Rear Suspension)	NC

Camber

<input type="checkbox"/>	132CI1	0°	NC
<input type="checkbox"/>	132CI2	3°	NC
<input type="checkbox"/>	132CI4	6° (NA w/ 5 spk Mag wheel)	NC

Rear Wheel

Ultra Lightweight only available w/ Aluminium Handrim. Contact Built 4 Me for Spinergy colour options

<input type="checkbox"/>	132RW103	20" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	132RW104	22" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	132RW102	24" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	132RW90	20" Lite Spoke	
<input type="checkbox"/>	132RW91	22" Lite Spoke	
<input type="checkbox"/>	132RW92	24" Lite Spoke	NC
<input type="checkbox"/>	132RW93	25" Lite Spoke	NC
<input type="checkbox"/>	132RW94	26" Lite Spoke	
<input type="checkbox"/>	132RW06	22" Ultra Lightweight Spoke Wheel	
<input type="checkbox"/>	132RW10	24" Ultra Lightweight Spoke Wheel	
<input type="checkbox"/>	132RW16	25" Ultra Lightweight Spoke Wheel	
<input type="checkbox"/>	132RW28	22" Spinergy SPOX Black Spokes	
<input type="checkbox"/>	132RW95	24" Spinergy SPOX Black Spokes	
<input type="checkbox"/>	132RW96	25" Spinergy SPOX Black Spokes	
<input type="checkbox"/>	132RW97	26" Spinergy SPOX Black Spokes	
<input type="checkbox"/>	132RW32	22" Spinergy LX Black Spokes	
<input type="checkbox"/>	132RW98	24" Spinergy LX Black Spokes	
<input type="checkbox"/>	132RW99	25" Spinergy LX Black Spokes	
<input type="checkbox"/>	132RW40	26" Spinergy LX Black Spokes	
<input type="checkbox"/>	132RW50	24" Mountain	

Power Add On

<input type="checkbox"/>		Quickie Xtender	See Order Form
<input type="checkbox"/>		Quickie iXpress	See Order Form

Please see corresponding order form for pricing, options and compatibility

Rear Wheel Tire

<input type="checkbox"/>	132RT1	Pneumatic	NC
<input type="checkbox"/>	132RT2	Pneumatic Airless Insert	
<input type="checkbox"/>	132RT3	Full Poly Tire	NC
<input type="checkbox"/>	132RT38	SHOX Performance Poly	
<input type="checkbox"/>	132RT39	High-Pressure Clincher Black	
<input type="checkbox"/>	132RT21	Schwalbe Marathon Plus Evolution	
<input type="checkbox"/>	132RT27	Knobby V-TRAK by Kenda	
<input type="checkbox"/>	132RT12	Mountain Only available in 24"	
<input type="checkbox"/>	132RT28	Schwalbe Speed Run Tire	

Front to Rear Seat Height Matrix

12" Seat Depth	2.5" Max Front to Rear Seat Hgt.
13", 14" Seat Depth	3.0" Max Front to Rear Seat Hgt.
15", 16" Seat Depth	3.5" Max Front to Rear Seat Hgt.
17", 18" Seat Depth	4.0" Max Front to Rear Seat Hgt.
19", 20" Seat Depth	4.5" Max Front to Rear Seat Hgt.

Tire	Wheel			
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy
Pneumatic	20,22,24	20,22,24,26		22,24,26
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26		22,24,26
Full Poly	20,22,24	20,22,24		
SHOX Performance Poly	20,22,24	20,22,24,25		24,25
High-Pressure Clinch Blk		22,24,26	22,24	22,24,26
Schwalbe Marathon		24,25,26	22,24,25	24,25,26
Knobby V Trak by Kenda		24		24
Schwalbe Speed Run		24,25	24,25	24,25

20" = 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590

Handrim

Natural Fit Handrims only available in 22", 24", & 25".

Natural Fit w/ 5 spoke MAG only available in 24". See APPENDIX 11 for matrix

<input type="checkbox"/>	132HR1	Aluminium Anodized	NC
<input type="checkbox"/>	132HR4	Plastic Coated NA w/ 25" & 26" Lite Spoke Wheels	
<input type="checkbox"/>	132HR22	Natural Fit - Standard Grip	
<input type="checkbox"/>	132HR23	Natural Fit - Super Grip	
<input type="checkbox"/>	132HR25	Natural Fit - No Thumb Grip	
<input type="checkbox"/>	132HR30	Natural Fit LT - Standard Grip	
<input type="checkbox"/>	132HR31	Natural Fit LT - Super Grip	
<input type="checkbox"/>	132HR32	Natural Fit LT - No Thumb Grip	
<input type="checkbox"/>	132HR26	The Surge - Oval w/Gription Strip	
<input type="checkbox"/>	132HR27	The Surge LT - Oval w/Gription Strip	
<input type="checkbox"/>	132HR28	Q-Grip Neoprene Coated	
<input type="checkbox"/>	132HR24	Spinergy Flex Rim Only available on Spinergy LX wheels 24" & 25"	
<input type="checkbox"/>	132HR7	Omit Handrim NA w/ Ultra Lightweight Wheels	NC

4. REAR WHEELS / SEAT HEIGHT (CONT)

Axle			
<input type="checkbox"/>	132AX1	Quick Release Steel	NC
<input type="checkbox"/>	132AX13	Quad Quick Release Steel Only available with Mountain Whls	

Axle Plate			
<input type="checkbox"/>		Adjustable COG (1"-4")	NC
<input type="checkbox"/>		Adjustable Extended COG (0"-3")	NC

Rear Suspension (Optional)			
<input type="checkbox"/>	4-Link Rear Suspension COG (0"-3") Single Post only available up to 2" COG. Plastic Side Guards only available up to 2.5" COG.		
<input type="checkbox"/>	132SV08	Weight from 0 kgs to 79kgs	
<input type="checkbox"/>	132SV09	Weight from 80kgs to 120kgs	

Center of Gravity Preset			
<i>Measure from front of back post to the center of rear axle.</i>			
<input type="checkbox"/>	132CC12	0"	NC
<input type="checkbox"/>	132CC1	.5"	NC
<input type="checkbox"/>	132CC2	1.0"	NC
<input type="checkbox"/>	132CC3	1.5"	NC
<input type="checkbox"/>	132CC4	2.0"	NC
<input type="checkbox"/>	132CC5	2.5"	NC
<input type="checkbox"/>	132CC6	3.0"	NC
<input type="checkbox"/>	132CC7	3.5"	NC
<input type="checkbox"/>	132CC8	4.0"	NC

5. WHEEL LOCKS

Wheel Lock			
<i>Not all wheel locks available in all seat depths, wheel size, or COG</i>			
<input type="checkbox"/>	132WL1	High Mount Push	NC
<input type="checkbox"/>	132WL2	High Mount Pull	NC
<input type="checkbox"/>	132WL45	Ergo Scissor	NC
<input type="checkbox"/>	132WL30	Compact	
<input type="checkbox"/>	132WL10	Omit Wheel Lock	NC

Wheel Lock Options			
<input type="checkbox"/>	132WL12	Do Not Mount	NC

Extension Handle			
<i>Ext. handles only available with high mount push or pull wheel locks</i>			
<input type="checkbox"/>	132EH7	Extension Handle (Pair)	

6. FOOTREST OPTIONS

Footrest			
<input type="checkbox"/>	132F229	Aluminium Tubular	NC
<input type="checkbox"/>	132F274	Aluminium Tubular w/ Snap On Cover	
<input type="checkbox"/>	132F276	Angle Adjustable Flip Back	
<input type="checkbox"/>	132F325	Angle Adjustable FB Carbon Fiber	
<input type="checkbox"/>	132F51	Platform Footrest	
<input type="checkbox"/>	132F111	Ultra Lite Platform Footrest NA w/ 14" wide and 2" inset, NA w/ 19-20 wide w/ 0" inset	
<input type="checkbox"/>	132F323	Platform Carbon Fiber footrest Not available on 19"-20" W	

Adaptive Platform Spacer			
<i>NA with Aluminium Tubular or Ultra Lite Platform.</i>			
<input type="checkbox"/>	132AD1	Adaptive Platform 2" Spacer	
<input type="checkbox"/>	132AD2	Adaptive Platform 4" Spacer	
<input type="checkbox"/>	132AD3	Adaptive Platform 6" Spacer	

Footrest Options			
<input type="checkbox"/>	132DI08	High Mount Ext Tubes (NA w/ Tubular footrest)	
<input type="checkbox"/>	132DI03	Leg Strap	
<input type="checkbox"/>	132DI04	Extra Leg Strap	
<input type="checkbox"/>	132DI05	Gel Hanger Cover	
<input type="checkbox"/>	132DI07	Neoprene Impact Guards	
<input type="checkbox"/>	132DI02	Wear Pads	
<input type="checkbox"/>	132DI06	Position Plates NA for Tubular Footrest	

7. BACKREST

Back Type

Freestyle Backrest not available on Active Rigid chairs. See option 3 for details

<input type="checkbox"/>	132BT3	Non-Folding Angle Adjustable	NC
<input type="checkbox"/>	132BT96	Folding w/ Double Lock Down Angle Adj	NC
<input type="checkbox"/>	132BT121	Freestyle Backrest Non- folding	
<input type="checkbox"/>	132BT122	Freestyle Backrest Folding	

Height Adjustable Integral Push Handles are only available with folding back. NA with non-folding back.

Choose Option 1 or 2

1. Aluminium Height Adjustable Backrest

<input type="checkbox"/>	132B110	Setting	10 " - 14" Back Height	NC
<input type="checkbox"/>	132B127	Setting	12 " - 16" Back Height	NC
<input type="checkbox"/>	132B111	Setting	14 " - 18" Back Height	NC
<input type="checkbox"/>	132B137	Setting	16 " - 20" Back Height	NC

OR

2. Ultra Lightweight Aluminium Fixed Backrest

Select back height. 8", 9", or 10" only available with omit push handles.
17", 18", or 19" only available with fold down push handles

<input type="checkbox"/>	132BB21	8"	<input type="checkbox"/>	132BB17	14"
<input type="checkbox"/>	132BB20	9"	<input type="checkbox"/>	132BB18	15"
<input type="checkbox"/>	132BB13	10"	<input type="checkbox"/>	132BB19	16"
<input type="checkbox"/>	132BB14	11"	<input type="checkbox"/>	132BB10	17"
<input type="checkbox"/>	132BB15	12"	<input type="checkbox"/>	132BB11	18"
<input type="checkbox"/>	132BB16	13"	<input type="checkbox"/>	132BB12	19"

Height Adjustable Integral Push Handles are only available with folding back. NA with non-folding back.

Towel Bar

<input type="checkbox"/>	132QT01	Standard Style - 4" Offset	NC
<input type="checkbox"/>	132QT02	Low Profile - 2" Offset Limited w/Solid Backs, Not suggested with 19" & 20" wide	

Push Handles

<input type="checkbox"/>	132PH6	Integral Push Handles	NC
<input type="checkbox"/>	132PH5	Fold Down Push Handles NA w/ 10"-14" adj height & 8"-10" fixed height NA w/ standard upholstery	
<input type="checkbox"/>	132PH2	Omit Push Handles	NC

Backrest Upholstery

<input type="checkbox"/>	132BU1	Standard Upholstery	NC
<input type="checkbox"/>	132BU70	Backrest Tension Adj Ballistic Nylon *Reference 186BU70 for billing	
<input type="checkbox"/>	132BU71	Backrest Tension Adj 3DX Vented *Reference 186BU71 for billing	
<input type="checkbox"/>	132BU88	Standard Upholstery EXO	
<input type="checkbox"/>	132BU89	Backrest Tension Adj. EXO	
<input type="checkbox"/>	132BU14	Omit Back Upholstery	

JAY J3 Back

<input type="checkbox"/>	132JB01	J3 Back To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #s. NA w/Freestyleback	Price List
--------------------------	---------	---	------------

Military Service Patches

Only available on tension adjustable upholstery and when ordered separate through JAY J3 Back OF

<input type="checkbox"/>	132VU1	US Marines Patch
<input type="checkbox"/>	132VU2	US Army Patch
<input type="checkbox"/>	132VU3	US Navy Patch
<input type="checkbox"/>	132VU4	US Air Force Patch
<input type="checkbox"/>	132VU12	Coast Guard Patch
<input type="checkbox"/>	132VU5	Purple Heart Patch
<input type="checkbox"/>	132VU8	Purple Heart with Marines Patch
<input type="checkbox"/>	132VU9	Purple Heart with Army Patch
<input type="checkbox"/>	132VU10	Purple Heart with Navy Patch
<input type="checkbox"/>	132VU11	Purple Heart with Air Force Patch
<input type="checkbox"/>	132VU13	Purple Heart with Coast Guard Patch

Backrest Options

Telescoping Push Handles have Double "D" tubing to eliminate rotation.

<input type="checkbox"/>	132OD5	Narrow Backrest Inset Backrests 1" narrower than selected seat width. NA with 12" Wide. NA w/ Aluminium or Carbon Side Guards
<input type="checkbox"/>	132A192	Telescoping Push Handles Range of 13-19" from sling to handle
<input type="checkbox"/>	132A195	Telescoping Push Handles - Tall (25" from seat to push handle)

Freestyle Backrest Option 3

Accessories such as head supports & lateral supports are not available with Freestyle backrest. The arch back comes standard with a small lumbar positioning component

NA w/ any Side Guard, chair Seat depth will be reduced by 1" with backrest. Available with Plastic Side Guard + Limited Carbon Side Guard

Backrest Type

<input type="checkbox"/>	132JB06	JAY J3 Back LT Aluminium model	
<input type="checkbox"/>	132JB07	JAY J3 Back LT Carbon Fiber model	

Backrest Height

Medium size backrest not available in Carbon Fiber.

<input type="checkbox"/>	132QF04	Small 6.5" shell height - 9" Back Height	NC
<input type="checkbox"/>	132QF05	Medium 8" shell height-10.5" Back Height	NC
<input type="checkbox"/>	132QF06	Tall 9.5" shell height - 12" Back Height	NC

Backrest Colour


Not applicable with carbon fiber back

<input type="checkbox"/>	132GZ01	Black	NC
<input type="checkbox"/>	132GZ02	Matching Frame- Standard Colourss	NC
<input type="checkbox"/>	132GZ03	Matching Frame- Custom Colour	

Backrest Width

<input type="checkbox"/>	132JQ01	14" Backrest Width	NC
<input type="checkbox"/>	132JQ02	16" Backrest Width	NC
<input type="checkbox"/>	132JQ03	18" Backrest Width	NC
<input type="checkbox"/>	132JQ04	20" Backrest Width	NC

Height Ranges	Dim "A"	Dim "B"
Small	12.0" - 14.5"	3.0" - 5.5"
Medium	12.0" - 16.0"	1.5" - 5.5"
Tall	12.0" - 17.5"	0.5" - 5.5"



8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away

Some limitations when used with the Angle Adjustable Back in the rear most setting

Please select a setting, will default to 11" if no selection is made

<input type="checkbox"/>	132AA18	E0973	Padded Swing-Away STD (9" to 14")	
	Setting	<input type="checkbox"/> 9"	<input type="checkbox"/> 10"	<input type="checkbox"/> 11"
		<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"

Armrest - Height Adjustable

<input type="checkbox"/>	132AA11		Single Post Height Adjustable - Std Range 8" - 13". Adjust in 1" increments	
<input type="checkbox"/>	132AA29		Single Post Height Adjustable - Low Range 6.5" - 10.5". Adjust in 1" increments	

Armypad - For Height Adjustable

<input type="checkbox"/>	132AM23		Desk Length (10") Classic	NC
<input type="checkbox"/>	132AM21		Desk Length (10") Contour	NC
<input type="checkbox"/>	132AM24		Full Length (14") Classic	NC
<input type="checkbox"/>	132AM22		Full Length (14") Contour	NC

Side Guards

Increases overall width of chair. Some interference may exist w/ padded SA arms at lowest settings Not available w/Height Adjustable Arms

<input type="checkbox"/>	132SG33		Plastic Standard	
<input type="checkbox"/>	132SG29		Plastic Low	
<input type="checkbox"/>	132SG35		Aluminium NA w/High Push or Pull Wheel Locks	
<input type="checkbox"/>	132SG36		Aluminium w/ Fender NA w/ High Push or Pull Wheel Locks	
<input type="checkbox"/>	132SG12		Carbon Fiber NA w/ High Push or Pull Wheel Locks	
<input type="checkbox"/>	132SG41		Carbon Fiber Long w/LG Fender NA w/High Push or Pull Wheel Locks	

Side Guard - Colour Option

Available with Aluminium or Aluminium with Fender. **Not applicable to carbon**

fiber options

<input type="checkbox"/>	132GE01		Black	NC
<input type="checkbox"/>	132GE02		Matching Frame- Standard Colours	NC
<input type="checkbox"/>	132GE03		Matching Frame- Custom Colour	

9. KEY OPTIONS AND ACCESSORIES

Anti Tips

<input type="checkbox"/>	132AT37		Anti Tips (Pair)	
<input type="checkbox"/>	132AT14		Anti Tips (Single) - Flip Up	
		<input type="checkbox"/> Left	<input type="checkbox"/> Right (facing from rear)	

Additional Pack Accessories

<input type="checkbox"/>	132A1		Backpack - Adult	
<input type="checkbox"/>	132A2		Backpack - Kids	
<input type="checkbox"/>	1321A01		Seat Pouch Black	

Positioning Belt

<input type="checkbox"/>	132A824		Belt 2" Aircraft Buckle	
<input type="checkbox"/>	132A825		Belt 2" Aircraft Buckle Padded	
<input type="checkbox"/>	132A826		Belt 1.5" Auto Buckle Universal	
<input type="checkbox"/>	132A804		Positioning Belt Clamps Only	

Spoke Guards

Only available with Lite spoke

<input type="checkbox"/>	1324A01		Quickie Black Spoke Guard	
<input type="checkbox"/>	1324A02		Quickie Clear Spoke Guard	

Miscellaneous

<input type="checkbox"/>	132A2100		Crutch Holder	
<input type="checkbox"/>	1328A21		Tool Kit	
<input type="checkbox"/>	132NP01		Gel Swing Away Pad	

10. Colour

Standard Frame Colour

<input type="checkbox"/>	132C140		Aztec Gold	NC
<input type="checkbox"/>	132C118		Bead Blast	NC
<input type="checkbox"/>	132C2		Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	132C141		Black Cherry	NC
<input type="checkbox"/>	132C10		Black Opal	NC
<input type="checkbox"/>	132C38		Blue Opal	NC
<input type="checkbox"/>	132C36		Candy Blue	NC
<input type="checkbox"/>	132C14		Candy Purple	NC
<input type="checkbox"/>	132C35		Candy Red	NC
<input type="checkbox"/>	132C142		Electric Blue	NC
<input type="checkbox"/>	132C143		Evergreen	NC
<input type="checkbox"/>	132C44		Glossy White - Limited Edition	NC
<input type="checkbox"/>	132C126		Glow	NC
<input type="checkbox"/>	132C144		Green Apple	NC
<input type="checkbox"/>	132C147		Hot Sparkle Pink	NC
<input type="checkbox"/>	132C148		Mauve Pink	NC
<input type="checkbox"/>	132C103		Rootbeer	NC
<input type="checkbox"/>	132C114		Sparkle Silver	NC
<input type="checkbox"/>	132C145		Sunrise Orange	NC
<input type="checkbox"/>	132C113		Titanium Colour Paint	NC
<input type="checkbox"/>	132C8		Yellow	NC

Matte Frame Colours

<input type="checkbox"/>	132C116		Matte Black	NC
<input type="checkbox"/>	132C151		Matte Black Cherry	NC
<input type="checkbox"/>	132C150		Matte Electric Blue	NC
<input type="checkbox"/>	132C149		Matte Evergreen	NC
<input type="checkbox"/>	132C152		Matte Purple	NC

Custom Frame Colour

<input type="checkbox"/>	132C112		Desert Camouflage	
<input type="checkbox"/>	132C43		Mossy Oak Camouflage	
<input type="checkbox"/>	132C42		Stars & Stripes	
<input type="checkbox"/>	132C108		Zebra (Choose above base colour) Std with White unless above base paint is chosen NA w/ Blk, Matte Black, Blk Opal, Titanium or Bead Blast)	

10. COLOUR (CONT)

Anodized Colour Package

Forks, Backcane Bracket, Rear Axle Clamp, and Freestyle mounting bracket (if selected). NA w/ Frog Legs & Carbotecture forks.

<input type="checkbox"/>	132EB07	Black	NC
<input type="checkbox"/>	132EB08	Blue	
<input type="checkbox"/>	132EB13	Gold	
<input type="checkbox"/>	132EB09	Lime Green	
<input type="checkbox"/>	132EB10	Purple	
<input type="checkbox"/>	132EB11	Red	
<input type="checkbox"/>	132EB12	Ti Grey	


Frame Graphics

Graphic designs on front frame tubes. NA with transit & Matte Black Cherry.

NA w short height frames

<input type="checkbox"/>	1328A24	Bloom	
<input type="checkbox"/>	1328A25	Dragon	
<input type="checkbox"/>	1328A26	Ghost Flames	
<input type="checkbox"/>	1328A27	Retro Squares	
<input type="checkbox"/>	1328A28	Tribal	

Bloom design




Dragon design



Ghost Flames design



Retro Squares design



Tribal design



Built For Me

To inquire about any unique customization for Quickie products, contact:
Phone: 0800 543 343



Commonly Requested Options Include:

Swing Away Knee Adductors for Oval Tubing
One Arm Drive
Unilateral Wheel Locks
Seat Pan
Caster Pin Lock
Grade Aids
Active Rigid w/Folding Back
Spinergy Colour Spokes
Custom Colour casters (Alum soffroll only)
24", 25" Glance wheels
Lifestyle custom paint package (see website)



11. APPENDIX

Handrim & Wheel Matrix

Handrims	Wheels				
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain
Aluminium Anodized	20,22,24	20,22,24,26	22,24,25	22,24,25,26	24
Plastic Coat	20,22,24	22,24	N/A	22,24,25,26	24
Natural Fit - Standard Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Natural Fit - No Thumb Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Natural Fit - Super Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Natural Fit LT - Std Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Natural Fit LT - Super Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Natural Fit LT - No Thumb Grip	24	22,24,25,26	N/A	22,24,25,26	N/A
Surge - Oval w/Gription Strip	N/A	22,24,25,26	N/A	22,24,25,26	N/A
Surge LT - Oval w/Gription Strip	N/A	22,24,25,26	N/A	22,24,25,26	N/A
Q Grip Neoprene Coated	N/A	24,25	N/A	24,25	N/A
Spinergy Flex Rim	N/A	N/A	N/A	24,25	N/A
Omit Handrim	20,22,24	N/A	N/A	22,24,25,26	24