

ORDERING INFORMATION: Phone: 0800 543 343 • Email: solutions@medifab.com • www.medifab.co.nz

Account No: _____
 Date: _____
 P.O. #: _____
 Therapist: _____
 Customer Name: _____
 Rep Name & Phone: _____

Ship To: _____
 Name: _____
 Address: _____
 City / State: _____
 Postcode: _____
 Phone: _____

QUICKIE iXPRESS POWER ADD ON

- 1568-9341-000 24V, 120W x 2 Motor

QUICKIE iXPRESS POWER ADD ON RETRO KIT

Includes Refurbishment Package of:-
 - iXpress Power Add on Kit & Refurbishment of existing chair for wear & tear items
 - Clean & service of chair: needs determined by SRM assessment to restore to
 - Includes if required: forks, casters, seat upholstery, back upholstery, arm pads &
 - Wheelchair submitted for refurbishment must meet a minimum level of functional
 - Chair suitability and assessment is to be performed by & agreed usage of the
 - Not to be used to alter size or functionality of the chair unless a growth kit is
- Please obtain a RMA# from Sunrise Medical prior to delivering the chair to Sunrise
 0 24V, 120W x 2 Motor

BASE FRAME

Select the wheelchair frame from the list below.
 The chair is not included in the kit.

Please also attach the completed Order Form for the model selected.

- Quickie 2
 Quickie 2HP
 Quickie GPV
 Zippie IRIS
 Quickie IRIS

JOYSTICK HAND CONTROL

Please see Figure 1
 Only available in 16" & 24" wheel sizes

- 1568-9310-024 JWX-1 Joystick (24" Wheels)
 1568-9310-016 JWX-1 Joystick (16" Wheels)

Fig 1: Joystick



JOYSTICK HAND CONTROL MOUNT TYPE

- Standard STD
 Swing Away (Armrests are required on the chair for this mount)

JOYSTICK HAND CONTROL MOUNT LOCATION

- Left
 Right

ATTENDANT CONTROL

Please see Figure 2

- 1568-9315-000 Dual Control (i.e. Attendant & Joystick)
 1568-9313-000 Attendant Control Only (i.e. no Joystick) (RHS only)
 Attendant Control (Joystick Only)
 Right Hand Side
 Left Hand Side

Fig 2: Attendant Control



DISCLAIMER: Not all features & options offered are compatible with all configurations of the wheelchair.

BATTERY

- 1568-9381-000 NiMH Battery + Charger
 1568-9385-000 Li Ion Battery + Charger

Fig 3: Battery



SPARE BATTERY

- 1568-9381-900 Spare NiMH Power Pack
 1568-9383-900 Spare Li Ion Power Pack

SPARE MOUNTING BRACKET - HAND CONTROL TO BE SWITCHED TO BACK CANE MOUNT WITH NORMAL MOUNT-ON ARMREST

- 1568-9318-100 Hand Control used as an Attendant Control Unit Dual Mount (can be swapped to user control)

- Maximum weight capacity is 120kgs on Q2, GPV Quickie IRIS
- Maximum weight capacity is 75 kgs Zippie IRIS
- Transit option is not tested with the iXpress
- 3 & 4" casters are N/A with the iXpress
- iXpress N/A with Q2 depth-adjustable backrest
- The Quickie iXpress kit will come with a new axle plate that will allow use of std manual wheels. Brake assembly may require adjustment.
- Hemi height frames are not compatible with the Quickie iXpress
- iXpress Kit comes standard with anti-tips.
- STFH's will be determined by base chair configuration

Specifications	
Drive System	Direct rear wheel drive
Controller	JWX-1, JUX-1
Motor	24 V, 120 W x 2
Battery	Nickel Metal Hydride
Range	16kms per charge
Charging Time	3 hours
Attendant Control	Yes
Drive Wheel Size	16" or 24"
Speed	up to 6kph
Weight Capacity	75, 120kgs
Battery Weight	2.9kgs
System Weight	17.4kgs