

Poz'In'Form®



CLINICAL LYING & POSITIONING SOLUTIONS

Distributed by **medifab®**

SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Therapist Name:	
Clients DOB:		Address:	
Gender:			
Contact Telephone:			
Parent/Caregiver:			
Application Reference:		Therapist Telephone:	
Date:		Therapist Mobile:	
NDIS/Funding:		Therapist Email:	
Notes:			

- Poz'In'Form® elements and modules are entirely removable and washable. They are delivered with a Pharmatex or Lenzing FR® cover.
- The Pharmatex bi-elastic cover is made of polyurethane and is recommended for use in hospitals. Pharmatex is soft to the touch and very comfortable next to the skin, while still easy to clean.
- The Lenzing FR® cover is made of viscose, coated with a waterproof polyurethane film. Lenzing FR® is a major innovation in positioning of the long-term bedridden, preventing maceration. Hygroscopic properties of viscose allow moisture absorption. Moreover elasticity and flexibility of Lenzing FR® are particularly suitable for problems regularly observed in classical positioning.
- Lenzing spare covers can be purchased to cover Pharmatex pillows to a warmer and more absorbent feel.



- Pharmatex Cover (Wipe-down)



- Lenzing FR® Cover (Breathable)

Poz'In'Form Kit

Includes:

- Transport bag
- Universal Cushion 55cm x 40cm
- Half moon cushion
- Semi-lateral decubitus cushion
- Cylindrical cushion
- Ring cushion
- Heel cushion with strap
- Occiput head cushion
- Pharmatex covers or Lenzing covers



With Pharmatex Cover	With Lenzing FR® Cover
7009-0001-269	7009-0001-279
Available for trial	Not available for trial

Poz'In'Form Delta

With Pharmatex Cover	With Lenzing FR® Cover
32cm wide x 56cm long x 28cm high	32cm wide x 56cm long x 28cm high
7009-0551-269	7009-0551-279
Available for trial	Not available for trial
+ Spare Cover 7009-0551-269-C	+ Spare Cover 7009-0551-279-C



NOTE: Pharmatex is stocked option, Lenzing indent order only (not available for trials).

Poz'In'Form Universal

With Pharmatex Cover

25cm wide x 35cm long	40cm wide x 55cm long	15cm wide x 40cm long
7009-1251-269 <input type="radio"/>	7009-1252-269 <input type="radio"/>	7009-1253-269 <input type="radio"/>
Available for trial	Available for trial	Available for trial
+ Spare Cover 7009-1251-269-C <input type="radio"/>	+ Spare Cover 7009-1252-269-C <input type="radio"/>	+ Spare Cover 7009-1253-269-C <input type="radio"/>

With Lenzing FR® Cover

25cm wide x 35cm long	40cm wide x 55cm long	15cm wide x 40cm long
7009-1251-279 <input type="radio"/>	7009-1252-279 <input type="radio"/>	7009-1253-279 <input type="radio"/>
Not available for trial	Not available for trial	Not available for trial
+ Spare Cover 7009-1251-279-C <input type="radio"/>	+ Spare Cover 7009-1252-279-C <input type="radio"/>	+ Spare Cover 7009-1253-279-C <input type="radio"/>



Poz'In'Form Cylindrical

With Pharmatex Cover ☐

With Lenzing FR® Cover ☐

75cm long x 18cm diameter	75cm long x 18cm diameter
7009-1501-269	7009-1501-279
Available for trial	Not available for trial
+ Spare Cover 7009-1501-269-C <input type="radio"/>	+ Spare Cover 7009-1501-279-C <input type="radio"/>



Poz'In'Form Abduction

With Pharmatex Cover ☐

With Lenzing FR® Cover ☐

35cm long x 15cm high	35cm long x 15cm high
Small Width: 10cm	Small Width: 10cm
Large Width: 30cm	Large Width: 30cm
7009-1551-269	7009-1551-279
Available for trial	Not available for trial
+ Spare Cover 7009-1551-269-C <input type="radio"/>	+ Spare Cover 7009-1551-279-C <input type="radio"/>



Poz'In'Form Half-Moon

With Pharmatex Cover <input type="radio"/>	With Lenzing FR® Cover <input type="radio"/>
200cm long x 30cm wide	200cm long x 30cm wide
7009-1851-269	7009-1851-279
Available for trial	Not available for trial
+ Spare Cover 7009-1851-269-C <input type="radio"/>	+ Spare Cover 7009-1851-279-C <input type="radio"/>




Poz'In'Form Semi-Lateral Decubitus

With Pharmatex Cover <input type="radio"/>	With Lenzing FR® Cover <input type="radio"/>
185cm long x 35cm wide	185cm long x 35cm wide
7009-3001-269	7009-3001-279
Available for trial	Not available for trial
+ Spare Cover 7009-3001-269-C <input type="radio"/>	+ Spare Cover 7009-3001-279-C <input type="radio"/>



Poz'In'Form Ring & Half Ring

With Pharmatex Cover		With Lenzing FR® Cover	
Ring <input type="radio"/>	Half Ring <input type="radio"/>	Ring <input type="radio"/>	Half Ring <input type="radio"/>
135cm long x 22cm wide	105cm long x 22cm wide	135cm long x 22cm wide	105cm long x 22cm wide
7009-4081-269	7009-4051-269	7009-4081-279	7009-4051-279
Available for trial	Available for trial	Not available for trial	Not available for trial
+ Spare Cover 7009-4081-269-C <input type="radio"/>	+ Spare Cover 7009-4051-269-C <input type="radio"/>	+ Spare Cover 7009-4081-279-C <input type="radio"/>	+ Spare Cover 7009-4051-279-C <input type="radio"/>



Poz'In'Form Occiput

With Pharmatex Cover <input type="radio"/>	With Lenzing FR® Cover <input type="radio"/>
50cm long x 9cm high x 40cm wide	50cm long x 9cm high x 40cm wide
7009-4751-269	7009-4751-279
Available for trial	Not available for trial
+ Spare Cover 7009-4751-269-C <input type="radio"/>	+ Spare Cover 7009-4751-269-C <input type="radio"/>



Poz'In'Form Heel Cushion

With Pharmatex Cover <input checked="" type="radio"/>	With Lenzing FR® Cover <input checked="" type="radio"/>
85cm long x 25cm wide	85cm long x 25cm wide
7009-5201-269	7009-5201-279
Available for trial	Not available for trial
+ Spare Cover 7009-5201-269-C <input type="radio"/>	+ Spare Cover 7009-5201-279-C <input type="radio"/>



Poz'In'Form Heel Wedge

With Pharmatex Cover <input checked="" type="radio"/>	With Lenzing FR® Cover <input checked="" type="radio"/>
48cm long x 12cm high x 30cm wide	48cm long x 12cm high x 30cm wide
7009-5401-269	7009-5401-279
Available for trial	Not available for trial
+ Spare Cover 7009-5401-269-C <input type="radio"/>	+ Spare Cover 7009-5401-279-C <input type="radio"/>



Poz'In'Form Hand Wedge

With Pharmatex Cover <input checked="" type="radio"/>	With Lenzing FR® Cover <input checked="" type="radio"/>
39cm long x 9cm high x 22cm wide	39cm long x 9cm high x 22cm wide
7009-5411-269	7009-5411-279
Available for trial	Not available for trial
+ Spare Cover 7009-5411-269-C <input type="radio"/>	+ Spare Cover 7009-5411-279-C <input type="radio"/>





Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Poziform_20181019_NZ