

Quickie 2 Folding Wheelchair

VERSATILE & MODULAR.
THE INDUSTRY'S FAVOURITE
FOR OVER 25 YEARS!

Q2

Standard Modular frame

Q2 Lite

Ultra lightweight, modular frame

Q2 HP

High Performance



MODULAR FRAME DESIGN



The versatility of the Quickie 2 is in its modular frame that adapts with you as your function, body and environment change. Independent, interchangeable frame components accommodate changing conditions without having to replace the entire wheelchair.

QUICKIE 2 LITE



Love the Quickie 2, but looking for something with an edge? The Quickie 2's minimalist counterpart, the Quickie 2 Lite, features an extremely lightweight, modular frame with a rounded rear for unique style and either fixed or swing-away front hangers.

Manual Wheelchair guide

KRYPTON R



7^R



5^R



KRYPTON F



Seat Width	32 - 46 cm (12-18")	12-20"	12-20"	32-46 cm (12.8-18.4")
Seat Depth	34 - 48 cm (14-20")	12-20"	12-22"	34-48 cm (13.6 - 19.2")
Front Seat Height	43- 55 cm (17-22")	16-21"	16-21"	43-55 cm (17.2-22")
Rear Seat Height	38 - 50 cm (15-20")	13-20"	13-20"	43-55 cm (17.2-22")
Frame Style	Carbon Fibre Open Frame	Open Series 7000 Aluminium	Open Series 5000 Aluminium	Carbon Fibre Open frame
Back Height	25 - 45 cm (10-18")	8-20"	10-20"	25 -47.5 cm (10-19")
Back Support Angle	-31° - 15°	-27° - 8°	-17° - 8°	75° to 103°
Frame / Hanger Angle	92° / 100°	75° / 80° / 85°	75° / 80° / 85°	92° / 100°
Inset	0" / 1"	0" / 1 1/2"	0" / 1 1/2"	0-1 cm
Camber	0 ° / 3 °	0 ° / 3 ° / 6 °	0°, 3°, 6°	0° / 2° / 4°
Max. User Weight	125 kg	120 kg	120 kg	110 kg
Product Weight	3.65 kg transit / 6.2kg	8.5 kg	9.8 kg	From 5.8kgs transit / 8.3kgs
Transit Approved	No	Yes*	Yes*	No
Frame Colour	5 Finishes	29 colours	28 colours	5 Finishes
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime
ARTG	100934	100167	100167	100934
Folding	No	No	No	Yes

*Conforms to ANDI/RESNA WC-19 Standards. Please refer to the Quickie & Zippie WC-19 Transit Securement Supplement

QUICKIE 2 LITE



QUICKIE 2



QXi



M6



IRIS



Seat Width	14-20"	12-22"	14-20"	22-30"	14-22"
Seat Depth	14-20"	12-20"	14-20"	18-22"	14-22"
Front Seat Height	14.5-21.5"	13.5 to 22"	12.5-19.5"	17-20"	12.5-19"
Rear Seat Height	13-19.5"	13 to 22"	12-19"	17-20"	12.5-19"
Frame Style	Traditional	- Fixed or SA front	Traditional	Heavy Duty	Rotation In Space
Back Height	9-21"	9-21"	13-20"	17-20"	15-24"
Back Support Angle	-14° - +5° (Optional)	-15° - +3° (Optional)	-15° - 3°	N/A	-5° to +30°
Frame / Hanger Angle	70° / 80°	60° / 70° / 80° / 90°	60° / 70° / 80°	N/A	60° / 70° / 80° / 90° / Elevating
Inset	0 / 1" each side	N/A	N/A	N/A	N/A
Camber	(Adjustable) Set at 3°	(Adjustable) Set at 3°	(Adjustable) Set at 3°	0°	0°
Max. User Weight	120 kg	136 kg Std / 158 kg HD	120kg	295 kg	136 kg / 181kgs HD
Product Weight	10.6 kg	11.5kg	13.25kg	24 kg	17.6 kg
Transit Approved	Yes	Yes	Yes	No	Yes
Frame Colour	24 colours	24 colours	24 colours	24 colours	25 colours
Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime
ARTG	100167	100167	100167	100167	100167
Folding	Yes	Yes	Yes	Yes	Yes