

# Savant Head Support

Savant is a unique head support which is cool, light-weight and can easily be shaped by hand. It provides support at three levels on the head to control lateral flexion and rotation. A headband can be used for anterior support.

## TOTAL HEAD CONTROL



Savant is constructed with a unique steel core that allows selective shaping without tools.



The mounting disk accommodates a 1" (25 mm) ball. The disk can be removed to attach custom hardware using M5 screws. Compatible with Spex head support mechanisms.



Savant comes with a headband that can be used for anterior support. The Savant can be worn with or without the headband.

A side pad for extra lateral support can be added if required: 1280-4210-100.

## BENEFITS FOR THE USER AND CAREGIVER



- Cool and comfortable head support, the modular pad design allows for air circulation
- Low heat retention
- Easy to put on and remove
- Low profile; unobtrusive appearance
- Easy to keep clean (sponge clean with soap and water, then air dry)
- Total control for the head

## OTHER BENEFITS



- Easy to fit (head support can be shaped by hand)
- Attaches to standard three point hardware
- Optional use of headband (included) for additional anterior support
- Pads are easily reshaped by hand when the occupant is repositioned

Size	Code	On the client		On the head support		
		Circumference at Forehead	1/2 Circumference behind ears	Inside distance along top support	Inside distance along middle support	Overall Height
<b>PAEDIATRIC</b>	1280-4101-000	14-18" / 36-46 cm	5-6" / 13-15 cm	9" / 23 cm	5" / 13 cm	4.5" / 12 cm
<b>SMALL</b>	1280-4102-000	18-22" / 46-56 cm	6-7" / 15-18 cm	12" / 30 cm	7" / 18 cm	5" / 13 cm
<b>ADULT</b>	1280-4104-000	22-24" / 56-61 cm	8-9" / 20-30 cm	14" / 36 cm	8" / 20 cm	6" / 15 cm