



## HD PRODUCT LINE

QUICKIE® heavy duty (HD) wheelchairs are designed to meet your needs. These high-capacity manual and power wheelchairs have reinforced frames and durable components ensuring strength and patient safety. Select QUICKIE chairs can accommodate 350 to 650 lbs. weight capacity. JAY® cushions and backs can accommodate 500 to 650 lbs.



LIVE  
WITH  
OUT  
LIMITS



## TABLE OF CONTENTS:

# HEAVY DUTY WHEELCHAIR OPTIONS

QUICKIE wheelchairs are designed for you and by you. In development, we consider how you interact with your wheelchair to maintain maximum efficiency, performance, fit, and style.



## MANUAL FOLDING

QUICKIE® 2 HD: **350 lbs. weight capacity** PG. 3

QUICKIE M6®: **650 lbs. weight capacity** PG. 4

## MANUAL RIGID

QUICKIE GP® Series: **350 lbs. weight capacity** PG. 5  
(Available through Built-4-Me)

## MANUAL TILT-IN-SPACE

QUICKIE IRIS® HD: **400 lbs. weight capacity** PG. 6

## POWER WHEELCHAIRS

QUICKIE P-222 SE: **350 lbs. weight capacity** PG. 7

QUICKIE QM-715®: **400 lbs. weight capacity** PG. 7  
(with tilt)  
**350 lbs. weight capacity**  
(with tilt and recline)



# HEAVY DUTY SEATING OPTIONS



## CUSHIONS

JAY® wheelchair cushions are designed to stabilize the pelvis, protect skin from the risk of breakdown, and position the body. Premium materials and positioning components create low maintenance and highly customizable seating solutions. Weight capacity dependent on cushion size.

JAY BasicPRO®: **500 lbs. weight capacity**

JAY GO: **500 lbs. weight capacity**

JAY ION®: **500 lbs. weight capacity**

JAY UNION®: **500 lbs. weight capacity**

JAY 2® Series **500 lbs. weight capacity**

(JAY 2 & JAY 2 Deep)

JAY 3®: **500 lbs. weight capacity**

JAY FUSION®: **500 lbs. weight capacity**

JAY 2® Plus: **650 lbs. weight capacity**

## BACKS

JAY wheelchair backs are designed to stabilize the pelvis and torso. Their contoured back shells, soft foam, and adjustable components offer stability and comfort to reduce fatigue and maximize function. Ride in confidence knowing we've got your back! Weight capacity dependent on hardware.

JAY J3®: **500 lbs. weight capacity**

JAY J2® PLUS: **650 lbs. weight capacity**

JAY GO: **500 lbs. weight capacity**

Back page will contain all JAY items listed here.



## FOLDING



VERSATILE AND  
MODULAR  
DESIGN ADAPTS  
TO YOUR  
NEEDS

### Key Features and Specs:

- ◆ Ultra-lightweight folding chair in swing-away or fixed front frame
- ◆ Modular frame option. Can upgrade from STD to HD if weight increases
- ◆ One free growth kit within first 5 years
- ◆ Transit option available
- ◆ Seat width 18" – 22"
- ◆ Seat depth up to 20"

## QUICKIE 2 HD

### Modular Frame Design

The versatility of the QUICKIE 2 HD is in its modular wheelchair frame that adapts with you as your function, body, and environment change. Independent, interchangeable frame components accommodate changing conditions without having to replace the entire wheelchair. You're always getting the best from your QUICKIE 2 HD.

### 7000 Series Aluminum

7000 Series Aerospace Aluminum makes the redesigned crossbrace more rigid, stronger, and lighter. The stronger frame results in an increased weight capacity of 300 lbs. (350 lbs. for HD), while the overall chair weight is reduced.



**HD weight capacity**  
**350 lbs.**

**Standard weight capacity**  
**300 lbs.**

## FOLDING



CUSTOM  
HIGH  
STRENGTH  
MODULAR  
FRAME

### Key Features and Specs:

- ◆ A variety of interchangeable frame components allow the M6 to adjust as your needs change
- ◆ Seat width 22" – 30"
- ◆ Seat depth 18" – 22"
- ◆ Seat height 17" – 20"
- ◆ Weight 53 lbs. with legrests

## QUICKIE M6



### Modular, High-Strength Frame

The M6's modular frame, modeled after the industry favorite QUICKIE 2, features high strength steel for increased wheelchair durability.



### Locking Crossbrace and Stabilizer Bar

The M6 utilizes a double locking crossbrace and stabilizer bar to minimize frame flex and create a more rigid ride quality. This ensures that more energy transfers directly from your arms or feet to the axles for the most efficient propulsion.

**HD weight capacity**  
**650 lbs.**

# BARIATRIC MID WHEEL XCEL 2

The mid-wheel Heavy Duty **Quickie® XCEL 2™** features the Advance Geometric Design (A.G.D) which helps enhance the stability, mobility and freedom for the end user. Engineered without compromise, the XCEL 2™ offers a weight capacity of 249kgs, a bariatric CG power tilt bariatric power recline option, and heavy duty elevating legrests or center mount for the best comfort.



Heavy Duty Tilt



## TILT-IN-SPACE



## SUPERIOR STANDARD FOR ROTATION-IN- SPACE TECHNOLOGY

### Key Features and Specs:

- ◆ One of the lightest weight tilt-in-space wheelchairs today (39 lbs.)
- ◆ Fixed or width adjustable strut tubes
- ◆ One free growth kit within first 5 years (only on non-SE models)
- ◆ Ultra low seat-to-floor height as low as 12.5"
- ◆ Seat width up to 24" and seat depths up to 22"
- ◆ Transit option available on HD model
- ◆ Reverse configuration available

## QUICKIE IRIS HD

### Intelligent Rotation In Space

The superior standard for rotational tilt-in space wheelchair technology, the QUICKIE IRIS features up to 55° of intelligent rotation in space technology, the broadest range of positioning and caregiver options, numerous adjustments, and it's one of the lightest tilt-in-space wheelchairs available today.



**HD weight capacity**  
400 lbs.

**Standard weight capacity**  
300 lbs.

## POWER

HIGH  
PERFORMING  
AND RELIABLE**Quickie P-222 SE****Key Features and Specs:**

- ◆ 8.5 mph, 4-pole standard motor package
- ◆ PosiDrive™ technology increases the overall strength of the drive train and features easy-to-service hubs
- ◆ 25" turning radius for excellent maneuverability

**Quickie Q700 HD****Key Features and Specs:**

- ◆ SpiderTrac® suspension – one of the most aggressive suspensions on the market today
- ◆ Multiple seating system options, including SEDEO® ERGO™
- ◆ 5 mph motors, with SureTrac® option
- ◆ Seat widths up to 24"  
Seat depths up to 22"  
(up to 24" through Built-4-Me)

## QUICKIE P-222 SE &amp; Q700 HD

**Quickie P-222 SE**

This power wheelchair offers the ultimate indoor and outdoor maneuverability and performance. Unique manual conversion kit for convenient transport available.

**Q700 HD**

QM-7 Series power wheelchairs combine function and luxury for a unique driving experience featuring SpiderTrac® 2.0 Suspension and SEDEO ERGO Seating

**P-222 SE weight capacity**

**350 lbs.**

**QM-715 HD weight capacity**

**400 lbs.** (With tilt)

**QM-715 HD weight capacity**

**350 lbs.** (with tilt and recline)

# JAY® Heavy Duty Seating Selection

## JAY CUSHIONS

These JAY cushions with up to 21" seat width have a standard weight capacity of 300 lbs. The weight capacity increases to 500 lbs. with 22" – 24" seat widths. J2 PLUS has a weight capacity of 650 lbs. on all of its sizes, available in 20" – 26" seat widths and 18" – 22" seat depths.



**BasicPRO®**  
E2602  
Up to 500 lbs.



**GO**  
E2604  
Up to 500 lbs.



**Ion®**  
E2608  
Up to 500 lbs.



**Union®**  
E2608  
Up to 500 lbs.



**Fusion®**  
E2623  
Up to 500 lbs.



**J3®**  
E2623, E2625  
Up to 500 lbs.



**J2® Series**  
E2623, E2625  
Up to 500 lbs.




**J2® PLUS**  
K0669  
Up to 650 lbs.

## JAY BACKS



**J3® PA**  
E2613, E2614



4-point hardware      TH hardware

- ◆ Available up to 500 lbs. capacity with 4-point hardware or TH hardware
- ◆ 2.2" of contour
- ◆ Widths available 22", 24" and 26" (dependent on height)
- ◆ Heights available 12", 15", 18", and 21"



**J2® PLUS**  
E2613, E2614

- ◆ Available up to 650 lbs. capacity
- ◆ 15 degrees of angle adjustability
- ◆ Widths available 20" – 26"



Vented Backshell



**GO**  
E2611, E2612

- ◆ Available up to 500 lbs. capacity
- ◆ Widths available 12" – 26"
- ◆ Heights available  
Short 10" – 14"  
Tall 15" – 19"

To learn about customizing JAY seating products, visit the JAY Your Way product page on [www.SunriseMedical.com](http://www.SunriseMedical.com).

Available from Medifab in New Zealand  
[www.medifab.co.nz](http://www.medifab.co.nz)    [solutions@medifab.com](mailto:solutions@medifab.com)    0800 543 343

medifab