

WHITMYER® HEADS UP™

 \equiv The premium solution for anterior head positioning



WHITMYER® HEADS UP™

The Whitmyer Heads Up is an innovative head support designed to maximize clients' interaction with their environment by providing anterior, posterior and lateral support of the head and neck. It features swing-out, adjustable arms to easily accommodate client transfers and varying client presentations.

Innovative Anterior Stabilization Arms

- Provide anterior and lateral stabilization of the shoulders and trunk for optimal head control
- Swing-out function allows for easy transfers while maintaining proper head support
- · Lower lateral adjustments accommodate varying body shapes

Available with Two Industry-Proven Posterior Support Options

- Contoured Cradle
- Sub-Occipital (shown to right)


Compatible with Full Whitmyer Hardware Offering

- Removable
- Swing-Away
- Fixed











SPECIFICATIONS

JAY

Client Neck Size (maximum)	Infant: 2.8" Small: 3.5" Medium: 4.1" Large: 5.3"	Attaching Hardware	AXYS ONYX LINX PRO
Pad Height	Short, Tall		COBRA XTRA COBRA XTRA Flip-Back
Arm Length	Short Infant/Small Pads: 5.6" Medium/Large Pads: 7.3" Long Infant/Small Pads: 7" Medium/Large Pads: 9"	Avail. Cover Materials	Lycra® Reverse Dartex®
		Avail. Foam Materials	Standard Comfort
		Occipital Pad Sizes	5", 6", 8", 10"

