

BasiX<sup>2</sup><sup>™</sup>



 Breezy<sup>®</sup>

**NEW**

## Attendant Controlled Drum Brake Option



### Optional Extras



Half Folding Backrest



24" Spoke Wheel Drum brakes



Height adjustable armrest with sliding pad for long and short position



Step Tube




Elevating legrest



Stump pad



Cast Support



Tray Table



Headrests



Crutch Holder



12" Transit Wheel Drum Brake



Attendant Control Drum Brake

# BasiX<sup>2</sup><sup>TM</sup>

*Ultra lightweight frame designed for personalised comfort*



*Available in both Self Propelled and Transit models, the Breezy BasiX<sup>2</sup> is a lightweight wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look.*

- Adapted seat heights accommodating for foot propulsion
- Armrests conveniently flip back, for easier transfers
- Easy accessible wheel locks prevent the wheelchair from rolling forward or backward
- Angle Adjustable Footplates
- Folding Frame for convenient transportation and storage



Frame Design includes seat depth adjustability (41 to 46 cm).



Tension adjustable backrest from 41cm to 46 cm.



Swing Away and Adjustable Anti-Tips



Height adjustable push handles





## Technical Specifications:



Max. User Weight:	125 kgs
Chair weight	14.9 kgs (12" wheel) / 17.5kg (24" wheel)
Overall width	Seat Width + 19cm / Seat Width + 20cm (drumbrake version)
Seat width	14, 16, 18 & 20"
Seat Depth	16 - 18" - adjustable
Seat Height	19" (adj. 18" - 20")
Back Height	16" - 18"
Armrest Height	8" - 9"
Armrest height adjustable	8" - 11"
Height adjustable push handles	5"
Hangers:	70° (lower leg length adjustable from 14" - 20") Angle adjustable footplate (+/- 15°)
Colour:	Robust hammerite - makes scratches less visible Red - on selected SKU's
Crash Test Approved	Yes
TGA#	145834
Warranty	Limited Lifetime