

# BINGO Evo



## PAEDIATRIC ADAPTIVE STROLLER

Distributed by  
**medifab**

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

April 2024

# Client Details

<b>Client/User Name:</b>		<b>Notes</b>
<b>Date:</b>		
<b>NDIS/Funding:</b>		

Measurements	Size 1	Size 2	Size 3
BINGO Evo Order Code	2881-0015-000	2882-0015-000	2883-0008-000
<b>Mobility Base with Swivelling Wheels</b>			
Push handle height (cm)	75-110	75-110	
Overall width (cm)	63	68	
Wheel diameter, front/swivel (cm)	19	19	
Wheel diameter, rear (cm)	29	29	
Tilt in space (Forward and rear facing)	45°	45° (Size 3 - Forward facing only)	
Folded size, mobility base (cm) (with rear wheels on)	80 x 63 x 38	80 x 68 x 38	
Folded size, mobility base (cm) (with rear wheels removed)	80 x 63 x 31	80 x 68 x 31	
Weight, mobility base (kg) (swivelling front wheels)	8.6	8.7	
Load capacity, mobility base (kg)	60	60	
<b>Seat Unit</b>			
Seat width (cm)	18-30	23-35	23-35
Seat depth (cm)	16-34	22-40	30-45
Lower leg length (cm)	19-40	19-40	19-40
Back support height (cm)	45-63	54-75	62-83
Height of slots for shoulder straps (cm)	28-46	34-55	43-64
Seat to floor height (cm) (dependant on tilt position of seat)	46	46	46
Back support recline	90-160°	90-160°	90-135°
Weight of seat (kg)	6.4	7.4	8.4
Load capacity of seat unit (kg)	50	50	50
Weight of seat and mobility base combined (kg)	15	16.1	17.1




## Stroller package includes:

- Sun canopy and rain cover
- Matte black foldable stroller frame with black wheels
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)
- Fabric bag (max. 3kg)
- Adjustable suspension set
- Receiver for grab rail or therapy tray
- Elevating foot support hanger (angle adjustable 90-160°)

**Note:** Upholstery not included in package



Items outlined in blue are required as standard package

# Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Upholstery (required)					
<b>Size 1</b> 2881-8801-347		<b>Size 1</b> 2881-8801-349		<b>Size 1</b> 2881-8801-346	
<b>Size 2</b> 2882-8802-347		<b>Size 2</b> 2882-8802-349		<b>Size 2</b> 2882-8802-346	
<b>Size 3</b> 2883-8803-347		<b>Size 3</b> 2883-8803-349		<b>Size 3</b> 2883-8803-346	

Hip Support Upholstery (required)					
<b>Size 1</b> 2881-1075-347		<b>Size 1</b> 2881-1075-349		<b>Size 1</b> 2881-1075-346	
<b>Size 2</b> 2882-1076-347		<b>Size 2</b> 2882-1076-349		<b>Size 2</b> 2882-1076-346	
<b>Size 3</b> 2883-1077-347		<b>Size 3</b> 2883-1077-349		<b>Size 3</b> 2883-1077-346	

# Head Supports

Spacer Fabric Head Support Pads		Head Support Speakers	
Required for Head Support Speakers to be used.		Requires Spacer Fabric Head Support Pads.	
<b>Size 1</b> 2880-4471-346 11cm (D) x 16cm (H)		2880-8330-000	
<b>Size 2</b> 2880-4472-346 13cm (D) x 18cm (H)			

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
<b>Merino/Spacer Fabric</b> 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		<b>Merino/Spacer Fabric</b> 2010-4210-099		<b>Merino/Spacer Fabric</b> 2010-4310-099	
		<b>Wipe Down Cover</b> 2010-4210-209		<b>Wipe Down Cover</b> 2010-4310-209	
		Large (21cm Internal width)		Large (21cm Internal width)	
<b>Wipe Down Cover</b> 2010-4110-209		<b>Merino/Spacer Fabric</b> 2010-4211-099		<b>Merino/Spacer Fabric</b> 2010-4311-099	
		<b>Wipe Down Cover</b> 2010-4211-209		<b>Wipe Down Cover</b> 2010-4311-209	

# Lateral Supports

## Standard

Height and angle adjustable. All sizes are compatible with all BINGO seat units

**Size 1**  
8cm(H) x 10cm(L)  
2880-3010-346

**Size 2**  
10cm(H) x 12cm(L)  
2880-3020-346

**Size 3**  
12cm(H) x 13cm(L)  
2880-3030-346



## Spex Narrow Straight

**Size 0**  
85mm (H) x  
115mm (D)  
2880-3008-100



## Spex Broad

**Size 1**  
105mm (H) x  
130mm (D)  
2880-3008-310



## Spex Teardrop


**Size 0**  
150mm (H) x  
105mm (D)  
2880-3008-500



# Sacral Support

## Sacral foam insert including spacer


22cm (W) x 9cm H x 1cm thick  
2880-2510-060



# Seat Options


## Thigh Wedge

Velcro attached under upholstery.  
22.5cm (W) x 10.5cm (D).  
Height: 2cm tapering down to 1cm  
(Recommended for  
size 1 BINGO Evo only).  
2006-0505-000



## Trochanter Wedges

Option to contour standard  
flat cushion. Velcro attached  
under upholstery.  
(Recommended for  
size 1 BINGO Evo only).  
2006-0660-000



# Abductor Pommel

## Swing down abductor pommel for increased abduction of the legs

All sizes are compatible with all BINGO seat units.



**Size 1**  
9cm(H) x 6cm(D)  
2880-1651-346


**Size 2**  
1cm(H) x 8cm(D)  
2880-1652-346

**Size 3**  
12.5cm(H) x 9cm(D)  
2880-1653-346



# Belt Options

Spex 2-Point Hip Belt		Wonderseat 4-Point Hip Belt	
<p><b>X-small</b> 2 point slim line, dual pull with neoprene sleeve padding. 1405-6220-027</p>		<p>Stabilises and controls the pelvis. Can be quickly removed for washing. 1406-5420-027</p>	


Groin Strap		
<p><b>To fit Size 1</b> 2880-6361-346</p>	<p><b>To fit Size 2 &amp; 3</b> 2880-6362-346</p>	

# Harness Options



Wonderseat Padded Centrepoint Shoulder Harness	
<p><b>Size 0</b> Pad height 22cm 1409-5660-017</p>	
<p><b>Size 1</b> Pad height 28cm 1409-5661-017</p>	

5-Point Harness		Spex 5-Point Harness	
<p>Includes front pull adjusters and neck padding. 2880-6511-346</p>		<p>Includes tee pad for babies and very small children. (Recommended for size 1 BINGO Evo only). 2006-6551-000</p>	

# Foot Positioning



Footplate Spacer		Ankle Footfasts	
<p>Spacer block to reduce lower leg length. 4cm thickness.</p>			
<p><b>Size 1</b> 2880-5321-000</p>	<p><b>Size 2 &amp; 3</b> 2880-5322-000</p>		
		<p>Requires technician to install. Purchase only. Inner circumference listed below:</p>	
		<p><b>Size 1</b> 16-25cm 2880-6101-000</p>	<p><b>Size 2</b> 21-29cm 2880-6102-000</p>

# Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			
<b>To fit Size 1</b> 2880-7111-000	<b>To fit Size 2 &amp; 3</b> 2880-7112-000		
		For smaller trunk cut out. Fits both tray sizes. 2880-7180-000	

Grab Rail	
Cannot be used with therapy tray	
<b>To fit Size 1</b> 2880-7251-346	<b>To fit Size 2 &amp; 3</b> 2880-7252-346
	

# Critical Care Options

IV Pole	Extra Bracket	Oxygen Bottle Holder
With mobility bracket for attachment to mobility base. Fits all sizes. 2880-8275-000		
	For QR mounting of IV pole to second base (Cobra or other indoor base) 2880-8275-050	With mounting bracket for attachment to mobility base. 120mm diameter for 'C Type' bottle. Fits all sizes. 2880-8253-000



Chassis Mount Tray	
<b>Size 1</b> 2880-8058-000	<b>Size 2 &amp; 3</b> 2880-8059-000
	

# Tubie Fun Covers

(Spacer Fabric Headsupport Pads required, must be used in conjunction with standard covers selected on page 2).

Personalise and protect your Bingo Evo Stroller with our brand-new washable covers from Tubie Fun!

Cover set includes seat liner with waterproof layer along with covers for hip support, shoulder wing and spacer fabric head support pad. Size 1 cover set includes size 1 head support pad covers, size 2 cover set includes size 2 head support pad covers.

Cover Set - Llamas		Cover Set - Planets	
<b>Size 1</b> 2880-9561-121		<b>Size 1</b> 2880-9561-036	
<b>Size 2</b> 2880-9562-121		<b>Size 2</b> 2880-9562-036	
<b>5 Point Harness Cover - Llamas</b> 2880-6511-121		<b>5 Point Harness Cover - Planets</b> 2880-6511-036	



# Hand Strap

Wrist safety tether strap

2005-6930-000



# Swivel Lock

To hold castors straight when brisk walking

For all sizes  
2880-8440-000



# Cobra Indoor Base

Cobra Hi-Low Indoor Frame with 75mm Twin Castors

The seat unit of the BINGO Evo easily transfers out of the stroller base and into the Cobra indoor base for feeding and activity times. Lower down the floor for playing and activities, and raise up to a table for meal time! Tilt settings can still be adjusted in the same way as the stroller base.

**Size 1**  
2880-8711-000

**Size 2 & 3**  
2880-8721-000



# Separate Seat Unit

Use with Cobra Base

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

**Size 1**  
2881-0011-000

**Size 2**  
2882-0012-000

**Size 3**  
2883-0013-000



# Twin Bingo Stroller

Connection mount that allows two strollers to be connected one in front of the other for an additional sibling. The strollers are able to be easily detached and used separately where required, providing flexibility when travelling with or without siblings.

This is not available for use with a size 3 Bingo Evo. This can be used to connect either: bases of the same size (and corresponding sized Bingo Spex strollers), or the Bingo Mini and Evo Size 1. It is not suitable for use between any other size combinations. Please fill out an additional script form to order the 2nd stroller unit with a unique configuration and accessories.

Please note that this can only be fitted at the time of ordering both stroller units. (Medifab cannot retrofit this attachment).

2880-8630-000



### Attendant Brake

Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected.

2880-8452-000

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: BINGO\_20240416\_NZ