

Carrot 3000



CAR SEAT

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE


Distributed by
medifab

Client Details

| | |
|-------------------|--|
| Client/User Name: | |
| Date: | |
| NDIS/Funding: | |

| Notes |
|-------|
| |


| User Measurements | |
|---------------------------|--|
| A Hip width | |
| B Seat depth | |
| C Seat to top of shoulder | |
| D Seat to top of head | |
| E Chest width | |
| F Seat to axilla | |
| G Leg length | |





Seat Specifications

| Age Range | User Weight Range | Seat Depth | Shoulder Height | Seat Weight |
|------------|-------------------|------------|---|-------------|
| 3-15 years | 13.6-49kg | 25-40cm | 32-47cm 45-60cm (when shoulder protector wings are added) | 8.7kg |

Packages

| Basic Package: | | |
|---|---------------|---|
| <ul style="list-style-type: none"> Carrot seat unit 5-point harness Headrest wings Tee-pad Free recline backrest | 4116-0010-259 |  |

| Behavioural Packages: | | |
|---|---|---|
| <ul style="list-style-type: none"> Carrot seat unit 5-point harness Headrest wings Tee-pad Free recline backrest | Small Package with Anti-Escape Buckle 4116-0010-026 |  |
| | Small Package with Pen-Key Buckle 4116-0010-025 |  |

Swivel Base Attachment

ISOfix essential.
 (Standard adjustment of mandatory
 footprop is 34-40cm).
 4116-0180-104



Swivel Base Protection Plate

Foam plate to protect car seat
 cover, mounts under swivel base.
 4116-0190-000




Seat Depth Growth

For growth, please select one of the following 3 rodded options, and associated extensions below, if required:
 (Rodded components cannot be used in combination with each other).

Option 1

Rodded 50mm Seat Depth Extension

4116-0281-259



50mm Growth Pad


4116-0285-259



Option 2

Rodded Waterfall Seat Depth Extension

4116-0283-259




50mm Growth Pad*

4116-0285-259



100mm Growth Pad*

4116-0289-259



Option 3

Rodded Foot Support

Small

Lower leg (D): 25-30 cm
 4116-5241-250

Medium

Lower leg (D): 30-40 cm
 4116-5242-250

Large

Lower leg (D): 40-50 cm
 4116-5243-250




50mm Growth Pad*

4116-0285-259



100mm Growth Pad*

4116-0289-259



Adaptor

For use of foot
 support with
 swivel base.
 4116-5120-000



*50mm and 100mm growth pads cannot be used in combination with each other.

Back Support Height Growth






| 50mm Height Extension | 100mm Height Extension | Shoulder Protector Wings |
|---|---|--|
| <p>50mm back support height extension infill pad. 4116-2481-259</p>  | <p>100mm back support height extension infill pad. 4116-2482-259</p>  | <p>(For children over 45cm shoulder height). 4116-2753-259</p>  |

*50mm and 100mm infill pads cannot be used in combination with each other












Installation Accessories

| Tether Strap (Optional) | ISOfix Lower Straps (Optional) | Airplane Fixing Strap |
|--|---|--|
| <p>No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6953-017</p>  | <p>No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6910-000</p>  | <p>Required for securing to aircraft seat. 4116-6960-017</p>  |

Cover Options

| Wipe-down | | |
|--|---|--|
| <p>Wipe-down seat base 4116-0476-209</p> |  | <p>Tee-pad, Standard 4116-6553-209</p> |
| <p>Wipe-down backrest (Main unit only). 4116-2476-209</p> |  | <p>Tee-pad, Large 4116-6558-209</p> |
|  | | |
| Removable Incontinence Cover | Air Mesh Cushion Pad | |
| <p>Attaches using Velcro. Unable to be attached to wipe down covers. 4116-0475-209</p> |  | <p>Attaches using Velcro. Unable to be attached to wipe down covers. 4116-8880-099</p>  |

Accessories (continued)

| Tray / Forearm Support | | Tilt Wedge | | Upgraded Lateral Supports | |
|---|---|---|---|---|---|
| 4116-5534-259 |  | Provides 20° tilt. 4116-0291-259 |  | Firm padding for added trunk control. 4116-3250-090 |  |
| Seat Belt Holder | | | Shoulder Strap Retainer | | |
| Designed to keep seat belt away from the neck. 4116-6790-017 |  | Small Chest Circumference 55-70cm 4116-6721-000 |  | Medium Chest Circumference 70-90cm 4116-6722-000 | |
| | | Large Chest Circumference 90-100cm 4116-6723-000 | | | |
| | | | | | |
| Medial Thigh Support | | Hinge Locking Belt | | | |
| Detachable & Washable. Attaches to chest support and crotch strap. 4116-1510-000 |  | For easy carrying. 4116-8170-000 |  | | |
| Hip Belt | | Accessory Cross Strap | | | |
| Small 4116-6531-027 |  | Anti-Escape 4604-3004-026 |  | | |
| Medium 4116-6532-027 | | Pen Release 4604-3004-025 | | | |
| Large 4116-6533-027 | | Side Release 4604-3112-019 | | | |
| Deeper Head Support Wings | | Larger Tee Pad | | | |
| 4116-4653-090 |  | 4116-6542-049 |  | | |
| | | Cover (Optional) 4116-6558-259 | | | |
| | | Black Spextex Cover (Optional) 4116-6558-209 | | | |

medifab

Distributed by Medifab Ltd
solutions@medifab.com
medifab.com



DISTRIBUTED BY: