

Duro Stroller



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

August 2024

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Product Details

Specifications	Size 1	Size 2
Turquoise	2780-8801-349	2780-8802-349
Black	2780-8801-346	2780-8802-346

Seat Specifications		
Seat width (cm)	23-34	28-39
Seat depth (cm)	26-40	31-45
Seat angle (cm)	16 or 20°	16 or 20°
Back height (cm)	45-63	61-83
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length (cm)	20-44	20-44
Folded size (L x W x H) (cm)	84 x 63 x 32	84 x 68 x 38

Base Specifications		
Max weight (kg)	75	75
Front wheel diameter (cm)	20	20
Rear wheel diameter (cm)	30	30



Stroller package includes:

- Hoggi Duro Rehab-Buggy
- Hip Support, Aluminium
- Seat and Backrest Upholsteries
- Small Spex Hip Belt
- Canopy
- PU Wheels with Black Rim
- Swivelling Front Wheels
- Supporting Brackets
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support with Heel Stop
- Accessories Bag (Max weight: 3kg)
- Tie Down Kit (According to ISO 7176-19)
- Retainer for Accessories
- User Manual




Items outlined in blue are required as standard package

Upholstery Colour Options

Turquoise		Black	
Seat Upholstery (required)			
Size 1 2780-8801-349		Size 1 2780-8801-346	
Size 2 2780-8802-349		Size 2 2780-8802-346	

Head Supports

Head Support Pads	
Size 1 11cm (D) x 16cm (H) 2880-4471-346	
Size 2 13cm (D) x 18cm (H) 2880-4472-346	



Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
Wipe Down Cover 2010-4210-209		Small (17cm Internal width)		Small (17cm Internal width)	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
Wipe Down Cover 2010-4211-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Supports


Standard Lateral Supports		
Height and angle adjustable. All sizes are compatible with both Duro sizes.		
Size 1 8cm (H) x 10cm (L) 2880-3010-346	Size 2 10cm (H) x 12cm (L) 2880-3020-346	Size 3 12cm (H) x 13cm (L) 2880-3030-346




Lateral Supports (continued)

Teardrop Lateral Supports		Contoured Lateral Supports	
<p>Size 1 13.5cm (H) x 7.5cm (L) 2006-3110-099</p>	<p>Size 2 15.5cm (H) x 10.5cm (L) 2006-3120-099</p>		<p>Spex Discovery Axial design 9cm (H) x 13cm (L) 2006-3030-099</p> 

Harness Options

Spex Centrepoint Harness	5-Point Harness
<p>X-Small Recommended for size 1 Duro Total pad dimensions: 26cm (H) x 22cm (W) 2006-6660-017</p>	<p>With front pull adjusters and neck padding 2880-6511-346</p> 
<p>Small Recommended for size 1 & 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6661-017</p>	
<p>Medium Recommended for size 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6662-017</p>	






Hip Belt Options

Spex Push Button	Pen Popper
<p>X-Small Pad width 30cm 1405-6210-028</p>	<p>X-Small Pad width 30cm 1405-6210-025</p>
<p>Small Pad width 36cm 1405-6211-028</p>	<p>Small Pad width 36cm 1405-6211-025</p>
<p>Medium Pad width 46cm 1405-6212-028</p>	<p>Medium Pad width 46cm 1405-6212-025</p>




Tray & Grab Rail Options


Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			
To fit Size 1 2880-7111-000	To fit Size 2 2880-7112-000		

Grab Rail	
Cannot be used with therapy tray	
To fit Size 1 2880-7251-346	To fit Size 2 2880-7252-346
	

Hand Strap

Wrist safety tether strap	
2005-6930-000	

Swivel Lock

To hold castors straight when brisk walking	
2880-8440-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Duro_20240821_NZ