

# Ella Stander



**STANDARD & COMPLEX  
STANDING AID**


Distributed by  
**medifab**

SCRIPTING AND ORDERING GUIDE

October 2024

# Client details

Client/User name:	
DOB:	
Therapist:	
Date:	
<b>User Measurements</b>	
A Hip width	
B Axilla height	
C Chest width	
D Overall user height	
E Heel to knee	
F Hip height	
User weight	




# Specifications

	Size 1	Size 2	Size 3	Size 4
Approx Age	0 - 4 years	3 - 7 years	6 - 12 years	11 - 18 years
User Height	60 - 95 cm	90 - 130 cm	120 - 155 cm	130 - 180 cm
Max User Weight	30 kg	40 kg	55 kg	80 kg
Trunk Pad Size	20 x 15 cm	25 x 20 cm	28 x 24 cm	34 x 28 cm
Trunk Width - Inset Laterals	12 - 20 cm	17 - 25 cm	19 - 30 cm	25 - 41 cm
Trunk Width - Straight Laterals	16 - 24 cm	21 - 29 cm	25 - 36 cm	31 - 47 cm
Hip Abduction Angle (Minimum)	0 - 18° per side	0 - 18° per side	0 - 18° per side	0 - 18° per side
Supine Angle Adjustment	0 - 90°	0 - 90°	0 - 90°	0 - 90°
Prone Angle Adjustment	0 - 90°	0 - 90°	0 - 90°	0 - 90°
Knee Cuff Width Inside Padding (Can Be Modified With More Padding)	8 cm	10 cm	10 cm	15 cm
Tray Size Width x Length	50 x 28 cm	55 x 33 cm	61 x 35 cm	61 x 35 cm
Tray Angle Adjustment Supine	40°	40°	40°	40°
Tray Angle Adjustment Prone	45°	45°	45°	45°

# Package Options

Standard Configuration	Complex Configuration
------------------------	-----------------------



**Package Includes:**

1. Mobile base with powered tilt function, abduction set, four lockable castors, upright column and adjustment brackets for the feet and knees.
2. Sandals
3. Knee Cups with Buckle Strap
4. Pelvic Support
5. Trunk Support
6. Belt
7. Straight Laterals
8. Tray
9. Headrest with Hardware
10. Measuring Tape on Upright Extrusion
11. Tilt Angle Indicator

**Note:** Standard package includes wider brackets at the knees and feet (not pictured).

**Package Includes:**

1. Mobile base with powered tilt function, abduction set, four lockable castors, upright column and adjustment brackets for the feet and knees.
2. Sandals
3. Knee Cups with Buckle Strap
4. De-rotational Pelvic Support
5. De-rotational Trunk Support
6. Belt
7. Straight Swing-Away Laterals
8. Straight Flexi-Pad Swing-Away
9. 10" Head Rest with Hardware
10. Shoulder Protraction Pads
11. Tray (Not Shown)
12. Harness Fitting Kit (Not Shown)
13. Measuring Tape on Upright Extrusion
14. Tilt Angle Indicator

Product Codes	Size 1	Size 2	Size 3	Size 4
Standard Package	5221-0020-000	5222-0020-000	5223-0020-000	5224-0020-000
Complex Package	N/A	5222-0030-000	5223-0030-000	5224-0030-000

# Accessories

**Note:** This section is for additional accessories over and above what is included in base packages.

Head Support Options	
<b>Adjust-A-Plush Standard</b> 8" - 18" (W) x 5" (H) 1480-4330-281	
<b>Adjust-A-Plush Medium</b> 7" - 16" (W) x 4" (H) 1480-4320-281	
<b>Adjust-A-Plush Narrow</b> 6" - 14" (W) x 3" (H) 1480-4310-281	

Foot Supports (Pair)		Knee Supports (Pair)	
Includes sandals with straps.		Includes knee cups, pads and straps.	
<b>Size 1</b> 134x76cm 5220-6071-000		<b>Size 1</b> Diameter: 80mm (knee cups) 5220-1761-000	
<b>Size 2</b> 184x83cm 5220-6072-000		<b>Size 2</b> Diameter: 100mm (knee cups) 5220-1762-000	
<b>Size 3</b> 235x114cm 5220-6073-000		<b>Size 3</b> Diameter: 120mm (knee cups) 5220-1763-000	


Lateral supports - includes hardware and pads - for pelvic and trunk support (pair)		
<b>Swing-away</b>	<b>Size 2</b> swing-away straight 5220-3321-000	
	<b>Size 3 &amp; 4</b> swing-away inset 27mm per side 5220-3330-000	
	<b>Size 3 &amp; 4</b> swing-away straight 5220-3331-000	
<b>Flexi-pad (swing-away)</b>	<b>Size 2</b> Flexi-pad swing-away straight 5220-3421-000	
	<b>Size 3 &amp; 4</b> Flexi-pad swing-away inset 27mm per side 5220-3430-000	
	<b>Size 3 &amp; 4</b> Flexi-pad swing-away straight 5220-3431-000	

# Accessories Continued...



## Tray with bowl

<b>Size 1</b> 5220-7181-000	
<b>Size 2</b> 5220-7182-000	
<b>Size 3</b> 5220-7183-000	

Shoulder protraction pads (pair)		Harness fitting kit	
----------------------------------	--	---------------------	--

Includes pads and harness fitting kit.		Includes brackets and cam locks. 5220-6150-000	
<b>Size 1</b> 5220-2751-000			
<b>Size 2</b> 5220-2752-000			

Rotational Lumbar Support		Knee Support Rotation - Module	
---------------------------	--	--------------------------------	--

<b>Small</b> 5220-2601-000		Enables the rotation of the Knee Support according to the rotation of the footplate. 5220-7200-000	
<b>Medium</b> 5220-2602-000			
<b>Large</b> 5220-2603-000			

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Ella Stander\_20241030